


# Thoracolumbar

## Pedicular Screws



12NM TORQ



ANTI TORQUE




BONE AWL

# Thoracolumbar

Pedicular Screws



**COMPRESSOR**



**INSUTU BENDER LEFT-RIGHT**



**DISTRACTOR**



**LYSTHESIS CUTTER**




**MALLET**



**TAB 4MM-5MM-6MM**

# Thoracolumbar

## Pedicular Screws



**NUT INSERTER**



**NUT INSERTER**



**PERSUADER**

# Thoracolumbar

Pedicular Screws



**PIN STOPPER**



**PROBE FEELER CURVED**



**POLYAXIAL SCREW DRIVER**



**PROBE FEELER STRAIGHT**



**TRANSVER NUT DRIVER**



**T HANDLE**

# Thoracolumbar

## Pedicular Screws



**PROBE  
(CURVED/STRAIGHT)**



**ROCKER**



**ROD BENDER**


# Thoracolumbar Pedicular Screws



**ROD GRIPPER**



**ROD ROTATER**



**STRONG ROD HOLDER**



**ROD HOLDER FORCEPS**



**ROD PUSHER**



**STRAIGHT HANDLE**

# Thoracolumbar Pedicular Screws

MCMSC-4.5x25	MC MULTIAXIAL SCREW CEMENT TYPE 4.5X25 MM	TI6AL4V-ELI	4.5x25mm
MCMSC-4.5x30	MC MULTIAXIAL SCREW CEMENT TYPE 4.5X30 MM	TI6AL4V-ELI	4.5x30mm
MCMSC-4.5x35	MC MULTIAXIAL SCREW CEMENT TYPE 4.5x35 MM	TI6AL4V-ELI	4.5x35mm
MCMSC-4.5x40	MC MULTIAXIAL SCREW CEMENT TYPE 4.5X40 MM	TI6AL4V-ELI	4.5x40mm
MCMSC-4.5x45	MC MULTIAXIAL SCREW CEMENT TYPE 4.5x45 MM	TI6AL4V-ELI	4.5x45mm
MCMSC-4.5x50	MC MULTIAXIAL SCREW CEMENT TYPE 4.5X50 MM	TI6AL4V-ELI	4.5x50mm
MCMSC-5x30	MC MULTIAXIAL SCREW CEMENT TYPE 5X30 MM	TI6AL4V-ELI	5x30mm
MCMSC-5x35	MC MULTIAXIAL SCREW CEMENT TYPE 5x35 MM	TI6AL4V-ELI	5x35mm
MCMSC-5x40	MC MULTIAXIAL SCREW CEMENT TYPE 5X40 MM	TI6AL4V-ELI	5x40mm
MCMSC-5x45	MC MULTIAXIAL SCREW CEMENT TYPE 5x45 MM	TI6AL4V-ELI	5x45mm
MCMSC-5x50	MC MULTIAXIAL SCREW CEMENT TYPE 5X50 MM	TI6AL4V-ELI	5x50mm
MCMSC-5x55	MC MULTIAXIAL SCREW CEMENT TYPE 5X55 MM	TI6AL4V-ELI	5x55mm
MCMSC-5x60	MC MULTIAXIAL SCREW CEMENT TYPE 5X60 MM	TI6AL4V-ELI	5x60mm
MCMSC-5x65	MC MULTIAXIAL SCREW CEMENT TYPE 5X65 MM	TI6AL4V-ELI	5x65mm
MCMSC-5.5x30	MC MULTIAXIAL SCREW CEMENT TYPE 5.5X30 MM	TI6AL4V-ELI	5.5x30mm
MCMSC-5.5x35	MC MULTIAXIAL SCREW CEMENT TYPE 5.5x35 MM	TI6AL4V-ELI	5.5x35mm
MCMSC-5.5x40	MC MULTIAXIAL SCREW CEMENT TYPE 5.5X40 MM	TI6AL4V-ELI	5.5x40mm
MCMSC-5.5x45	MC MULTIAXIAL SCREW CEMENT TYPE 5.5x45 MM	TI6AL4V-ELI	5.5x45mm
MCMSC-5.5x50	MC MULTIAXIAL SCREW CEMENT TYPE 5.5X50 MM	TI6AL4V-ELI	5.5x50mm
MCMSC-5.5x55	MC MULTIAXIAL SCREW CEMENT TYPE 5.5X55 MM	TI6AL4V-ELI	5.5x55mm
MCMSC-5.5x60	MC MULTIAXIAL SCREW CEMENT TYPE 5.5X60 MM	TI6AL4V-ELI	5.5x60mm
MCMSC-5.5x65	MC MULTIAXIAL SCREW CEMENT TYPE 5.5X65 MM	TI6AL4V-ELI	5.5x65mm
MCMSC-6x30	MC MULTIAXIAL SCREW CEMENT TYPE 6X30 MM	TI6AL4V-ELI	6x30mm
MCMSC-6x35	MC MULTIAXIAL SCREW CEMENT TYPE 6x35 MM	TI6AL4V-ELI	6x35mm
MCMSC-6x40	MC MULTIAXIAL SCREW CEMENT TYPE 6X40 MM	TI6AL4V-ELI	6x40mm
MCMSC-6x45	MC MULTIAXIAL SCREW CEMENT TYPE 6x45 MM	TI6AL4V-ELI	6x45mm
MCMSC-6x50	MC MULTIAXIAL SCREW CEMENT TYPE 6X50 MM	TI6AL4V-ELI	6x50mm
MCMSC-6x55	MC MULTIAXIAL SCREW CEMENT TYPE 6X55 MM	TI6AL4V-ELI	6x55mm
MCMSC-6x60	MC MULTIAXIAL SCREW CEMENT TYPE 6X60 MM	TI6AL4V-ELI	6x60mm
MCMSC-6x65	MC MULTIAXIAL SCREW CEMENT TYPE 6X65 MM	TI6AL4V-ELI	6x65mm

# Thoracolumbar

## Pedicular Screws

MCMSC-6.5x30	MC MULTIAXIAL SCREW CEMENT TYPE 6.5X30 MM	TI6AL4V-ELI	6.5x30mm
MCMSC-6.5x35	MC MULTIAXIAL SCREW CEMENT TYPE 6.5X35 MM	TI6AL4V-ELI	6.5x35mm
MCMSC-6.5x40	MC MULTIAXIAL SCREW CEMENT TYPE 6.5X40 MM	TI6AL4V-ELI	6.5x40mm
MCMSC-6.5x45	MC MULTIAXIAL SCREW CEMENT TYPE 6.5X45 MM	TI6AL4V-ELI	6.5x45mm
MCMSC-6.5x50	MC MULTIAXIAL SCREW CEMENT TYPE 6.5X50 MM	TI6AL4V-ELI	6.5x50mm
MCMSC-6.5x55	MC MULTIAXIAL SCREW CEMENT TYPE 6.5X55 MM	TI6AL4V-ELI	6.5x55mm
MCMSC-6.5x60	MC MULTIAXIAL SCREW CEMENT TYPE 6.5X60 MM	TI6AL4V-ELI	6.5x60mm
MCMSC-6.5x65	MC MULTIAXIAL SCREW CEMENT TYPE 6.5X65 MM	TI6AL4V-ELI	6.5x65mm
MCMSC-7x30	MC MULTIAXIAL SCREW CEMENT TYPE 7X30 MM	TI6AL4V-ELI	7x30mm
MCMSC-7x35	MC MULTIAXIAL SCREW CEMENT TYPE 7X35 MM	TI6AL4V-ELI	7x35mm
MCMSC-7x40	MC MULTIAXIAL SCREW CEMENT TYPE 7X40 MM	TI6AL4V-ELI	7x40mm
MCMSC-7x45	MC MULTIAXIAL SCREW CEMENT TYPE 7X45 MM	TI6AL4V-ELI	7x45mm
MCMSC-7x50	MC MULTIAXIAL SCREW CEMENT TYPE 7X50 MM	TI6AL4V-ELI	7x50mm
MCMSC-7x55	MC MULTIAXIAL SCREW CEMENT TYPE 7X55 MM	TI6AL4V-ELI	7x55mm
MCMSC-7x60	MC MULTIAXIAL SCREW CEMENT TYPE 7X60 MM	TI6AL4V-ELI	7.5x60mm
MCMSC-7x65	MC MULTIAXIAL SCREW CEMENT TYPE 7X65 MM	TI6AL4V-ELI	7.5x65mm
MCMSC-7.5x30	MC MULTIAXIAL SCREW CEMENT TYPE 7.5X30 MM	TI6AL4V-ELI	7.5x30mm
MCMSC-7.5x35	MC MULTIAXIAL SCREW CEMENT TYPE 7.5X35 MM	TI6AL4V-ELI	7.5x35mm
MCMSC-7.5x40	MC MULTIAXIAL SCREW CEMENT TYPE 7.5X40 MM	TI6AL4V-ELI	7.5x40mm
MCMSC-7.5x45	MC MULTIAXIAL SCREW CEMENT TYPE 7.5X45 MM	TI6AL4V-ELI	7.5x45mm
MCMSC-7.5x50	MC MULTIAXIAL SCREW CEMENT TYPE 7.5X50 MM	TI6AL4V-ELI	7.5x50mm
MCMSC-7.5x55	MC MULTIAXIAL SCREW CEMENT TYPE 7.5X55 MM	TI6AL4V-ELI	7.5x55mm
MCMSC-7.5x60	MC MULTIAXIAL SCREW CEMENT TYPE 7.5X60 MM	TI6AL4V-ELI	7.5x60mm
MCMSC-7.5x65	MC MULTIAXIAL SCREW CEMENT TYPE 7.5X65 MM	TI6AL4V-ELI	7.5x65mm
MCMSC-8x30	MC MULTIAXIAL SCREW CEMENT TYPE 8X30 MM	TI6AL4V-ELI	8x30mm
MCMSC-8x35	MC MULTIAXIAL SCREW CEMENT TYPE 8X35 MM	TI6AL4V-ELI	8x35mm
MCMSC-8x40	MC MULTIAXIAL SCREW CEMENT TYPE 8X40 MM	TI6AL4V-ELI	8x40mm
MCMSC-8x45	MC MULTIAXIAL SCREW CEMENT TYPE 8X45 MM	TI6AL4V-ELI	8x45mm
MCMSC-8x50	MC MULTIAXIAL SCREW CEMENT TYPE 8X50 MM	TI6AL4V-ELI	8x50mm
MCMSC-8x55	MC MULTIAXIAL SCREW CEMENT TYPE 8X55 MM	TI6AL4V-ELI	8x55mm
MCMSC-8x60	MC MULTIAXIAL SCREW CEMENT TYPE 8X60 MM	TI6AL4V-ELI	8x60mm
MCMSC-8x65	MC MULTIAXIAL SCREW CEMENT TYPE 8X65 MM	TI6AL4V-ELI	8x65mm

# Thoracolumbar Pedicular Screws



MCPSS-4.5X25	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 4.5X25 MM	TI6AL4V-ELI	4.5x25mm
MCPSS-4.5X30	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 4.5X30 MM	TI6AL4V-ELI	4.5x30mm
MCPSS-4.5X35	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 4.5X35 MM	TI6AL4V-ELI	4.5x35mm
MCPSS-4.5X40	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 4.5X40 MM	TI6AL4V-ELI	4.5x40mm
MCPSS-4.5X45	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 4.5X45 MM	TI6AL4V-ELI	4.5x45mm
MCPSS-4.5X50	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 4.5X50 MM	TI6AL4V-ELI	4.5x50mm
MCPSS-5X30	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 5X30 MM	TI6AL4V-ELI	5x30mm
MCPSS-5X35	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 5X35 MM	TI6AL4V-ELI	5x35mm
MCPSS-5X40	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 5X40 MM	TI6AL4V-ELI	5x40mm
MCPSS-5X45	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 5X45 MM	TI6AL4V-ELI	5x45mm
MCPSS-5X50	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 5X50 MM	TI6AL4V-ELI	5x50mm
MCPSS-5X55	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 5X55 MM	TI6AL4V-ELI	5x55mm
MCPSS-5X60	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 5X60 MM	TI6AL4V-ELI	5x60mm
MCPSS-5X65	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 5X65 MM	TI6AL4V-ELI	5x65mm
MCPSS-5.5X30	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 5.5X30 MM	TI6AL4V-ELI	5.5x30mm
MCPSS-5.5X35	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 5.5X35 MM	TI6AL4V-ELI	5.5x35mm
MCPSS-5.5X40	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 5.5X40 MM	TI6AL4V-ELI	5.5x40mm
MCPSS-5.5X45	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 5.5X45 MM	TI6AL4V-ELI	5.5x45mm
MCPSS-5.5X50	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 5.5X50 MM	TI6AL4V-ELI	5.5x50mm
MCPSS-5.5X55	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 5.5X55 MM	TI6AL4V-ELI	5.5x55mm
MCPSS-5.5X60	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 5.5X60 MM	TI6AL4V-ELI	5.5x60mm

# Thoracolumbar Pedicular Screws



MCPSS-5.5X65	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 5.5X65 MM	Ti6AL4V-ELI	5.5x65mm
MCPSS-6X30	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 6X30 MM	Ti6AL4V-ELI	6x30mm
MCPSS-6X35	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 6X35 MM	Ti6AL4V-ELI	6x35mm
MCPSS-6X40	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 6X40 MM	Ti6AL4V-ELI	6x40mm
MCPSS-6X45	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 6X45 MM	Ti6AL4V-ELI	6x45mm
MCPSS-6X50	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 6X50 MM	Ti6AL4V-ELI	6x50mm
MCPSS-6X55	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 6X55 MM	Ti6AL4V-ELI	6x55mm
MCPSS-6X60	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 6X60 MM	Ti6AL4V-ELI	6x60mm
MCPSS-6X65	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 6X65 MM	Ti6AL4V-ELI	6x65mm
MCPSS-6.5X30	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 6.5X30 MM	Ti6AL4V-ELI	6.5x30mm
MCPSS-6.5X35	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 6.5X35 MM	Ti6AL4V-ELI	6.5x35mm
MCPSS-6.5X40	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 6.5X40 MM	Ti6AL4V-ELI	6.5x40mm
MCPSS-6.5X45	MC POLYAXIAL SPONDYLOLISTHESIS SCREW -6.5X45 MM	Ti6AL4V-ELI	6.5x45mm
MCPSS-6.5X50	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 6.5X50 MM	Ti6AL4V-ELI	6.5x50mm
MCPSS-6.5X55	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 6.5X55 MM	Ti6AL4V-ELI	6.5x55mm
MCPSS-6.5X60	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 6.5X60 MM	Ti6AL4V-ELI	6.5x60mm
MCPSS-6.5X65	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 6.5X65 MM	Ti6AL4V-ELI	6.5x65mm
MCPSS-7X30	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 7X30 MM	Ti6AL4V-ELI	7x30mm
MCPSS-7X35	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 7X35 MM	Ti6AL4V-ELI	7x35mm
MCPSS-7X40	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 7X40 MM	Ti6AL4V-ELI	7x40mm
MCPSS-7X45	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 7X45 MM	Ti6AL4V-ELI	7x45mm
MCPSS-7X50	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 7X50 MM	Ti6AL4V-ELI	7x50mm

# Thoracolumbar

## Pedicular Screws



MCPSS-7X55	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 7X55 MM	TI6AL4V-ELI	7x55mm
MCPSS-7X60	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 7X60 MM	TI6AL4V-ELI	7x60mm
MCPSS-7X65	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 7X65 MM	TI6AL4V-ELI	7x65mm
MCPSS-7.5X30	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 7.5X30 MM	TI6AL4V-ELI	7.5x30mm
MCPSS-7.5X35	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 7.5X35 MM	TI6AL4V-ELI	7.5x35mm
MCPSS-7.5X40	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 7.5X40 MM	TI6AL4V-ELI	7.5x40mm
MCPSS-7.5X45	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 7.5X45 MM	TI6AL4V-ELI	7.5x45mm
MCPSS-7.5X50	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 7.5X50 MM	TI6AL4V-ELI	7.5x50mm
MCPSS-7.5X55	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 7.5X55 MM	TI6AL4V-ELI	7.5x55mm
MCPSS-7.5X60	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 7.5X60 MM	TI6AL4V-ELI	7.5x60mm
MCPSS-7.5X65	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 7.5X65 MM	TI6AL4V-ELI	7.5x65mm
MCPSS-8X30	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 8X30 MM	TI6AL4V-ELI	8x30mm
MCPSS-8X35	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 8X35 MM	TI6AL4V-ELI	8x35mm
MCPSS-8X40	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 8X40 MM	TI6AL4V-ELI	8x40mm
MCPSS-8X45	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 8X45 MM	TI6AL4V-ELI	8x45mm
MCPSS-8X50	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 8X50 MM	TI6AL4V-ELI	8x50mm
MCPSS-8X55	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 8X55 MM	TI6AL4V-ELI	8x55mm
MCPSS-8X60	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 8X60 MM	TI6AL4V-ELI	8x60mm
MCPSS-8X65	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 8X65 MM	TI6AL4V-ELI	8x65mm
MCPSS-8X70	MC POLYAXIAL SPONDYLOLISTHESIS SCREW 8X70 MM	TI6AL4V-ELI	8x70mm

# Thoracolumbar

## Pedicular Screws



MCPS-5X30	MC POLYAXIAL SCREW 5X30MM	TI6AL4V-ELI	5x30mm
MCPS-5X35	MC POLYAXIAL SCREW 5X35MM	TI6AL4V-ELI	5x35mm
MCPS-5X40	MC POLYAXIAL SCREW 5X40MM	TI6AL4V-ELI	5x40mm
MCPS-5X45	MC POLYAXIAL SCREW 5X45MM	TI6AL4V-ELI	5x45mm
MCPS-5X50	MC POLYAXIAL SCREW 5X50MM	TI6AL4V-ELI	5x50mm
MCPS-5X55	MC POLYAXIAL SCREW 5X55MM	TI6AL4V-ELI	5x55mm
MCPS-5X60	MC POLYAXIAL SCREW 5X60MM	TI6AL4V-ELI	5x60mm
MCPS-5X65	MC POLYAXIAL SCREW 5X65MM	TI6AL4V-ELI	5x65mm
MCPS-5.5X30	MC POLYAXIAL SCREW 5.5X30MM	TI6AL4V-ELI	5.5x30mm
MCPS-5.5X35	MC POLYAXIAL SCREW 5.5X35MM	TI6AL4V-ELI	5.5x35mm
MCPS-5.5X40	MC POLYAXIAL SCREW 5.5X40MM	TI6AL4V-ELI	5.5x40mm
MCPS-5.5X45	MC POLYAXIAL SCREW 5.5X45MM	TI6AL4V-ELI	5.5x45mm
MCPS-5.5X50	MC POLYAXIAL SCREW 5.5X50MM	TI6AL4V-ELI	5.5x50mm
MCPS-5.5X55	MC POLYAXIAL SCREW 5.5X55MM	TI6AL4V-ELI	5.5x55mm
MCPS-5.5X60	MC POLYAXIAL SCREW 5.5X60MM	TI6AL4V-ELI	5.5x60mm
MCPS-5.5X65	MC POLYAXIAL SCREW 5.5X65MM	TI6AL4V-ELI	5.5x65mm
MCPS-6X30	MC POLYAXIAL SCREW 6X30MM	TI6AL4V-ELI	6x30mm
MCPS-6X35	MC POLYAXIAL SCREW 6X35MM	TI6AL4V-ELI	6x35mm
MCPS-6X40	MC POLYAXIAL SCREW 6X40MM	TI6AL4V-ELI	6x40mm
MCPS-6X45	MC POLYAXIAL SCREW 6X45MM	TI6AL4V-ELI	6x45mm
MCPS-6X50	MC POLYAXIAL SCREW 6X50MM	TI6AL4V-ELI	6x50mm
MCPS-6X55	MC POLYAXIAL SCREW 6X55MM	TI6AL4V-ELI	6x55mm
MCPS-6X60	MC POLYAXIAL SCREW 6X60MM	TI6AL4V-ELI	6x60mm

# Thoracolumbar

## Pedicular Screws



MCPS-6X65	MC POLYAXIAL SCREW 6X65MM	Ti6AL4V-ELI	6x65mm
MCPS-6.5X30	MC POLYAXIAL SCREW 6.5X30MM	Ti6AL4V-ELI	6.5x30mm
MCPS-6.5X35	MC POLYAXIAL SCREW 6.5X35MM	Ti6AL4V-ELI	6.5x35mm
MCPS-6.5X40	MC POLYAXIAL SCREW 6.5X40MM	Ti6AL4V-ELI	6.5x40mm
MCPS-6.5X45	MC POLYAXIAL SCREW 6.5X45MM	Ti6AL4V-ELI	6.5x45mm
MCPS-6.5X50	MC POLYAXIAL SCREW 6.5X50MM	Ti6AL4V-ELI	6.5x50mm
MCPS-6.5X55	MC POLYAXIAL SCREW 6.5X55MM	Ti6AL4V-ELI	6.5x55mm
MCPS-6.5X60	MC POLYAXIAL SCREW 6.5X60MM	Ti6AL4V-ELI	6.5x60mm
MCPS-6.5X65	MC POLYAXIAL SCREW 6.5X65MM	Ti6AL4V-ELI	6.5x65mm
MCPS-7X30	MC POLYAXIAL SCREW 7X30MM	Ti6AL4V-ELI	7x30mm
MCPS-7X35	MC POLYAXIAL SCREW 7X35MM	Ti6AL4V-ELI	7x35mm
MCPS-7X40	MC POLYAXIAL SCREW 7X40MM	Ti6AL4V-ELI	7x40mm
MCPS-7X45	MC POLYAXIAL SCREW 7X45MM	Ti6AL4V-ELI	7x45mm
MCPS-7X50	MC POLYAXIAL SCREW 7X50MM	Ti6AL4V-ELI	7x50mm
MCPS-7X55	MC POLYAXIAL SCREW 7X55MM	Ti6AL4V-ELI	7x55mm
MCPS-7X60	MC POLYAXIAL SCREW 7X60MM	Ti6AL4V-ELI	7x60mm
MCPS-7X65	MC POLYAXIAL SCREW 7X65MM	Ti6AL4V-ELI	7x65mm
MCPS-7.5X30	MC POLYAXIAL SCREW 7.5X30MM	Ti6AL4V-ELI	7.5x30mm
MCPS-7.5X35	MC POLYAXIAL SCREW 7.5X35MM	Ti6AL4V-ELI	7.5x35mm
MCPS-7.5X40	MC POLYAXIAL SCREW 7.5X40MM	Ti6AL4V-ELI	7.5x40mm
MCPS-7.5X45	MC POLYAXIAL SCREW 7.5X45MM	Ti6AL4V-ELI	7.5x45mm
MCPS-7.5X50	MC POLYAXIAL SCREW 7.5X50MM	Ti6AL4V-ELI	7.5x50mm
MCPS-7.5X55	MC POLYAXIAL SCREW 7.5X55MM	Ti6AL4V-ELI	7.5x55mm

# Thoracolumbar

## Pedicular Screws



MCPS-7.5X60	MC POLYAXIAL SCREW 7.5X60MM	Ti6AL4V-ELI	7.5x60mm
MCPS-7.5X65	MC POLYAXIAL SCREW 7.5X65MM	Ti6AL4V-ELI	7.5x65mm
MCPS-8X30	MC POLYAXIAL SCREW 8X30MM	Ti6AL4V-ELI	8x30mm
MCPS-8X35	MC POLYAXIAL SCREW 8X35MM	Ti6AL4V-ELI	8x35mm
MCPS-8X40	MC POLYAXIAL SCREW 8X40MM	Ti6AL4V-ELI	8x40mm
MCPS-8X45	MC POLYAXIAL SCREW 8X45MM	Ti6AL4V-ELI	8x45mm
MCPS-8X50	MC POLYAXIAL SCREW 8X50MM	Ti6AL4V-ELI	8x50mm
MCPS-8X55	MC POLYAXIAL SCREW 8X55MM	Ti6AL4V-ELI	8x55mm
MCPS-8X60	MC POLYAXIAL SCREW 8X60MM	Ti6AL4V-ELI	8x60mm
MCPS-8X65	MC POLYAXIAL SCREW 8X65MM	Ti6AL4V-ELI	8x65mm

# Thoracolumbar

## Pedicular Screws



MCTTB-50	MC CROSSLINK 50 MM	Ti6AL4V-ELI	50mm
MCTTB-60	MC CROSSLINK 60 MM	Ti6AL4V-ELI	60mm
MCTTB-70	MC CROSSLINK 70 MM	Ti6AL4V-ELI	70mm
MCTTB-80	MC CROSSLINK 80 MM	Ti6AL4V-ELI	80mm
MCTTB-90	MC CROSSLINK 90 MM	Ti6AL4V-ELI	90mm

# Thoracolumbar

## Pedicular Screws



MCPMS-4.5X25	MC POSTERIOR MONOAXIAL SCREW 4.5X25	TI6AL4V-ELI	4.5x25mm
MCPMS-4.5X30	MC POSTERIOR MONOAXIAL SCREW 4.5X30	TI6AL4V-ELI	4.5x30mm
MCPMS-4.5X35	MC POSTERIOR MONOAXIAL SCREW 4.5X35	TI6AL4V-ELI	4.5x35mm
MCPMS-4.5X40	MC POSTERIOR MONOAXIAL SCREW 4.5X40	TI6AL4V-ELI	4.5x40mm
MCPMS-4.5X45	MC POSTERIOR MONOAXIAL SCREW 4.5X45	TI6AL4V-ELI	4.5x45mm
MCPMS-4.5X50	MC POSTERIOR MONOAXIAL SCREW 4.5X50	TI6AL4V-ELI	4.5x50mm
MCPMS-5X30	MC POSTERIOR MONOAXIAL SCREW 5X30	TI6AL4V-ELI	5x30mm
MCPMS-5X35	MC POSTERIOR MONOAXIAL SCREW 5X35	TI6AL4V-ELI	5x35mm
MCPMS-5X40	MC POSTERIOR MONOAXIAL SCREW 5X40	TI6AL4V-ELI	5x40mm
MCPMS-5X45	MC POSTERIOR MONOAXIAL SCREW 5X45	TI6AL4V-ELI	5x45mm
MCPMS-5X50	MC POSTERIOR MONOAXIAL SCREW 5X50	TI6AL4V-ELI	5x50mm
MCPMS-5X55	MC POSTERIOR MONOAXIAL SCREW 5X55	TI6AL4V-ELI	5x55mm
MCPMS-5X60	MC POSTERIOR MONOAXIAL SCREW 5X60	TI6AL4V-ELI	5x60mm
MCPMS-5X65	MC POSTERIOR MONOAXIAL SCREW 5X65	TI6AL4V-ELI	5x65mm
MCPMS-5.5X30	MC POSTERIOR MONOAXIAL SCREW 5.5X30	TI6AL4V-ELI	5.5x30mm
MCPMS-5.5X35	MC POSTERIOR MONOAXIAL SCREW 5.5X35	TI6AL4V-ELI	5.5x35mm
MCPMS-5.5X40	MC POSTERIOR MONOAXIAL SCREW 5.5X40	TI6AL4V-ELI	5.5x40mm
MCPMS-5.5X45	MC POSTERIOR MONOAXIAL SCREW 5.5X45	TI6AL4V-ELI	5.5x45mm
MCPMS-5.5X50	MC POSTERIOR MONOAXIAL SCREW 5.5X50	TI6AL4V-ELI	5.5x50mm
MCPMS-5.5X55	MC POSTERIOR MONOAXIAL SCREW 5.5X55	TI6AL4V-ELI	5.5x55mm
MCPMS-5.5X60	MC POSTERIOR MONOAXIAL SCREW 5.5X60	TI6AL4V-ELI	5.5x60mm
MCPMS-5.5X65	MC POSTERIOR MONOAXIAL SCREW 5.5X65	TI6AL4V-ELI	5.5x65mm
MCPMS-6X30	MC POSTERIOR MONOAXIAL SCREW 6X30	TI6AL4V-ELI	6x30mm
MCPMS-6X35	MC POSTERIOR MONOAXIAL SCREW 6X35	TI6AL4V-ELI	6x35mm
MCPMS-6X40	MC POSTERIOR MONOAXIAL SCREW 6X40	TI6AL4V-ELI	6x40mm
MCPMS-6X45	MC POSTERIOR MONOAXIAL SCREW 6X45	TI6AL4V-ELI	6x45mm

# Thoracolumbar

## Pedicular Screws



MCPMS-6X50	MC POSTERIOR MONOAXIAL SCREW 6X50	Ti6AL4V-ELI	6x50mm
MCPMS-6X55	MC POSTERIOR MONOAXIAL SCREW 6X55	Ti6AL4V-ELI	6x55mm
MCPMS-6X60	MC POSTERIOR MONOAXIAL SCREW 6X60	Ti6AL4V-ELI	6x60mm
MCPMS-6X65	MC POSTERIOR MONOAXIAL SCREW 6X65	Ti6AL4V-ELI	6x65mm
MCPMS-6.5X30	MC POSTERIOR MONOAXIAL SCREW 6.5X30	Ti6AL4V-ELI	6.5x30mm
MCPMS-6.5X35	MC POSTERIOR MONOAXIAL SCREW 6.5X35	Ti6AL4V-ELI	6.5x35mm
MCPMS-6.5X40	MC POSTERIOR MONOAXIAL SCREW 6.5X40	Ti6AL4V-ELI	6.5x40mm
MCPMS-6.5X45	MC POSTERIOR MONOAXIAL SCREW 6.5X45	Ti6AL4V-ELI	6.5x45mm
MCPMS-6.5X50	MC POSTERIOR MONOAXIAL SCREW 6.5X50	Ti6AL4V-ELI	6.5x50mm
MCPMS-6.5X55	MC POSTERIOR MONOAXIAL SCREW 6.5X55	Ti6AL4V-ELI	6.5x55mm
MCPMS-6.5X60	MC POSTERIOR MONOAXIAL SCREW 6.5X60	Ti6AL4V-ELI	6.5x60mm
MCPMS-6.5X65	MC POSTERIOR MONOAXIAL SCREW 6.5X65	Ti6AL4V-ELI	6.5x65mm
MCPMS-7X30	MC POSTERIOR MONOAXIAL SCREW 7X30	Ti6AL4V-ELI	7x30mm
MCPMS-7X35	MC POSTERIOR MONOAXIAL SCREW 7X35	Ti6AL4V-ELI	7x35mm
MCPMS-7X40	MC POSTERIOR MONOAXIAL SCREW 7X40	Ti6AL4V-ELI	7x40mm
MCPMS-7X45	MC POSTERIOR MONOAXIAL SCREW 7X45	Ti6AL4V-ELI	7x45mm
MCPMS-7X50	MC POSTERIOR MONOAXIAL SCREW 7X50	Ti6AL4V-ELI	7x50mm
MCPMS-7X55	MC POSTERIOR MONOAXIAL SCREW 7X55	Ti6AL4V-ELI	7x55mm
MCPMS-7X60	MC POSTERIOR MONOAXIAL SCREW 7X60	Ti6AL4V-ELI	7x60mm
MCPMS-7X65	MC POSTERIOR MONOAXIAL SCREW 7X65	Ti6AL4V-ELI	7x65mm
MCPMS-7.5X30	MC POSTERIOR MONOAXIAL SCREW 7.5X30	Ti6AL4V-ELI	7.5x30mm
MCPMS-7.5X35	MC POSTERIOR MONOAXIAL SCREW 7.5X35	Ti6AL4V-ELI	7.5x35mm
MCPMS-7.5X40	MC POSTERIOR MONOAXIAL SCREW 7.5X40	Ti6AL4V-ELI	7.5x40mm
MCPMS-7.5X45	MC POSTERIOR MONOAXIAL SCREW 7.5X45	Ti6AL4V-ELI	7.5x45mm
MCPMS-7.5X50	MC POSTERIOR MONOAXIAL SCREW 7.5X50	Ti6AL4V-ELI	7.5x50mm
MCPMS-7.5X55	MC POSTERIOR MONOAXIAL SCREW 7.5X55	Ti6AL4V-ELI	7.5x55mm

# Thoracolumbar

## Pedicular Screws



MCPMS-7.5X60	MC POSTERIOR MONOAXIAL SCREW 7.5X60	TI6AL4V-ELI	7.5x60mm
MCPMS-7.5X65	MC POSTERIOR MONOAXIAL SCREW 7.5X65	TI6AL4V-ELI	7.5x65mm
MCPMS-8X30	MC POSTERIOR MONOAXIAL SCREW 8X30	TI6AL4V-ELI	8x30mm
MCPMS-8X35	MC POSTERIOR MONOAXIAL SCREW 8X35	TI6AL4V-ELI	8x35mm
MCPMS-8X40	MC POSTERIOR MONOAXIAL SCREW 8X40	TI6AL4V-ELI	8x40mm
MCPMS-8X45	MC POSTERIOR MONOAXIAL SCREW 8X45	TI6AL4V-ELI	8x45mm
MCPMS-8X50	MC POSTERIOR MONOAXIAL SCREW 8X50	TI6AL4V-ELI	8x50mm
MCPMS-8X55	MC POSTERIOR MONOAXIAL SCREW 8X55	TI6AL4V-ELI	8x55mm
MCPMS-8X60	MC POSTERIOR MONOAXIAL SCREW 8X60	TI6AL4V-ELI	8x60mm
MCPMS-8X65	MC POSTERIOR MONOAXIAL SCREW 8X65	TI6AL4V-ELI	8x65mm





# Thoracolumbar

## Pedicular Screws

MCIC - 745	MC ILIAC SCREW 7X45 MM	Ti6AL4V-ELI	7x45mm
MCIC - 750	MC ILIAC SCREW 7X50 MM	Ti6AL4V-ELI	7x50mm
MCIC - 760	MC ILIAC SCREW 7X60 MM	Ti6AL4V-ELI	7x60mm
MCIC - 770	MC ILIAC SCREW 7X70 MM	Ti6AL4V-ELI	7x70mm
MCIC - 780	MC ILIAC SCREW 7X80 MM	Ti6AL4V-ELI	7x80mm
MCIC - 790	MC ILIAC SCREW 7X90 MM	Ti6AL4V-ELI	7x90mm
MCIC - 7100	MC ILIAC SCREW 7X100 MM	Ti6AL4V-ELI	7x100mm
MCIC - 845	MC ILIAC SCREW 8X45 MM	Ti6AL4V-ELI	8x45mm
MCIC - 850	MC ILIAC SCREW 8X50 MM	Ti6AL4V-ELI	8x50mm
MCIC - 860	MC ILIAC SCREW 8X60 MM	Ti6AL4V-ELI	8x60mm
MCIC - 870	MC ILIAC SCREW 8X70 MM	Ti6AL4V-ELI	8x70mm
MCIC - 880	MC ILIAC SCREW 8X80 MM	Ti6AL4V-ELI	8x80mm
MCIC - 890	MC ILIAC SCREW 8X90 MM	Ti6AL4V-ELI	8x90mm
MCIC - 8100	MC ILIAC SCREW 8X100 MM	Ti6AL4V-ELI	8x100mm
MCIC - 945	MC ILIAC SCREW 9X45 MM	Ti6AL4V-ELI	9x45mm
MCIC - 950	MC ILIAC SCREW 9X50 MM	Ti6AL4V-ELI	9x50mm
MCIC - 960	MC ILIAC SCREW 9X60 MM	Ti6AL4V-ELI	9x60mm
MCIC - 970	MC ILIAC SCREW 9X70 MM	Ti6AL4V-ELI	9x70mm
MCIC - 980	MC ILIAC SCREW 9X80 MM	Ti6AL4V-ELI	9x80mm
MCIC - 990	MC ILIAC SCREW 9X90 MM	Ti6AL4V-ELI	9x90mm
MCIC - 9100	MC ILIAC SCREW 9X100 MM	Ti6AL4V-ELI	9x100mm

# Thoracolumbar Pedicular Screws




# Thoracolumbar

## Pedicular Screws



MCOLC MC Open Lateral Connector



MCEDNFR-6X50MM	MC EXPANDABLE DOUBLE NECK FLEX ROD 6X50MM	TI6AL4V-ELI	6x50mm
MCEDNBFR-6X80MM	MC EXPANDABLE DOUBLE NECK FLEX ROD 6X80MM	TI6AL4V-ELI	6x80mm
MCEDNFR-6X100MM	MC EXPANDABLE DOUBLE NECK FLEX ROD 6X100MM	TI6AL4V-ELI	6x100mm
MCEDNFR-6X120MM	MC EXPANDABLE DOUBLE NECK FLEX ROD 6X120MM	TI6AL4V-ELI	6x120mm

# Thoracolumbar

## Pedicular Screws



MCFR-3.5X60	MC FLEX ROD 3.5X60 MM	TI6AL4V-ELI	3.5x60mm
MCFR-3.5X80	MC FLEX ROD 3.5X80 MM	TI6AL4V-ELI	3.5x80mm
MCFR-4X60	MC FLEX ROD 4X60 MM	TI6AL4V-ELI	4x60mm
MCFR-4X80	MC FLEX ROD 4X80 MM	TI6AL4V-ELI	4x80mm
MCFR-5.5X40	MC FLEX ROD 5.5X40 MM	TI6AL4V-ELI	5.5x40mm
MCFR-5.5X50	MC FLEX ROD 5.5X50 MM	TI6AL4V-ELI	5.5x50mm
MCFR-5.5X60	MC FLEX ROD 5.5X60 MM	TI6AL4V-ELI	5.5x60mm
MCFR-5.5X70	MC FLEX ROD 5.5X70 MM	TI6AL4V-ELI	5.5x70mm
MCFR-5.5X80	MC FLEX ROD 5.5X80 MM	TI6AL4V-ELI	5.5x80mm
MCFR-5.5X90	MC FLEX ROD 5.5X90 MM	TI6AL4V-ELI	5.5x90mm
MCFR-5.5X100	MC FLEX ROD 5.5X100 MM	TI6AL4V-ELI	5.5x100mm
MCFR-5.5X110	MC FLEX ROD 5.5X110 MM	TI6AL4V-ELI	5.5x110mm
MCFR-5.5X120	MC FLEX ROD 5.5X120 MM	TI6AL4V-ELI	5.5x120mm
MCFR-5.5X130	MC FLEX ROD 5.5X130 MM	TI6AL4V-ELI	5.5x130mm
MCFR-5.5X140	MC FLEX ROD 5.5X140 MM	TI6AL4V-ELI	5.5x140mm
MCFR-5.5X150	MC FLEX ROD 5.5X150 MM	TI6AL4V-ELI	5.5x150mm
MCFR-5.5X160	MC FLEX ROD 5.5X160 MM	TI6AL4V-ELI	5.5x160mm
MCFR-5.5X170	MC FLEX ROD 5.5X170 MM	TI6AL4V-ELI	5.5x170mm
MCFR-5.5X180	MC FLEX ROD 5.5X180MM	TI6AL4V-ELI	5.5x180mm
MCFR-5.5X190	MC FLEX ROD 5.5X190 MM	TI6AL4V-ELI	5.5x190mm
MCFR-5.5X200	MC FLEX ROD 5.5X200 MM	TI6AL4V-ELI	5.5x200mm
MCFR-5.5X210	MC FLEX ROD 5.5X210 MM	TI6AL4V-ELI	5.5x210mm
MCFR-5.5X220	MC FLEX ROD 5.5X220MM	TI6AL4V-ELI	5.5x220mm
MCFR-5.5X230	MC FLEX ROD 5.5X230 MM	TI6AL4V-ELI	5.5x230mm

# Thoracolumbar

## Pedicular Screws




MCFR-5.5X240	MC FLEX ROD 5.5X240 MM	TI6AL4V-ELI	5.5x240mm
MCFR-5.5X250	MC FLEX ROD 5.5X250 MM	TI6AL4V-ELI	5.5x250mm
MCFR-6X40	MC FLEX ROD 6X40 MM	TI6AL4V-ELI	6x40mm
MCFR-6X50	MC FLEX ROD 6X50 MM	TI6AL4V-ELI	6x50mm
MCFR-6X60	MC FLEX ROD 6X60 MM	TI6AL4V-ELI	6x60mm
MCFR-6X70	MC FLEX ROD 6X70 MM	TI6AL4V-ELI	6x70mm
MCFR-6X80	MC FLEX ROD 6X80 MM	TI6AL4V-ELI	6x80mm
MCFR-6X90	MC FLEX ROD 6X90 MM	TI6AL4V-ELI	6x90mm
MCFR-6X100	MC FLEX ROD 6X100 MM	TI6AL4V-ELI	6x100mm
MCFR-6X110	MC FLEX ROD 6X110 MM	TI6AL4V-ELI	6x110mm
MCFR-6X120	MC FLEX ROD 6X120 MM	TI6AL4V-ELI	6x120mm
MCFR-6X130	MC FLEX ROD 6X130 MM	TI6AL4V-ELI	6x130mm
MCFR-6X140	MC FLEX ROD 6X140 MM	TI6AL4V-ELI	6x140mm
MCFR-6X150	MC FLEX ROD 6X150 MM	TI6AL4V-ELI	6x150mm
MCFR-6X160	MC FLEX ROD 6X160 MM	TI6AL4V-ELI	6x160mm
MCFR-6X170	MC FLEX ROD 6X170 MM	TI6AL4V-ELI	6x170mm
MCFR-6X180	MC FLEX ROD 6X180MM	TI6AL4V-ELI	6x180mm
MCFR-6X190	MC FLEX ROD 6X190 MM	TI6AL4V-ELI	6x190mm
MCFR-6X200	MC FLEX ROD 6X200 MM	TI6AL4V-ELI	6x200mm
MCFR-6X210	MC FLEX ROD 6X210 MM	TI6AL4V-ELI	6x210mm
MCFR-6X220	MC FLEX ROD 6X220MM	TI6AL4V-ELI	6x220mm
MCFR-6X230	MC FLEX ROD 6X230 MM	TI6AL4V-ELI	6x230mm
MCFR-6X240	MC FLEX ROD 6X240 MM	TI6AL4V-ELI	6x240mm
MCFR-6X250	MC FLEX ROD 6X250 MM	TI6AL4V-ELI	6x250mm

# Thoracolumbar

Pedicular Screws



MCMD-01	MC MONOBLOCK DOMINO SINGLE
MCMD	MC MONOBLOCK DOMINO
MCMD-02	MC MONOBLOCK DOMINO DOUBLE
MCMA	MC MONOBLOCK AXIAL
MCMC	MC MONOBLOCK CONNECTOR
MCAMC	MC ANGLED MONOBLOCK CONNECTOR



# Thoracolumbar

## Pedicular Screws



MCPPSS-45X25	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 4.5X25 MM	Ti6AL4V-ELI	4.5x25mm
MCPPSS-45X30	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 4.5X30 MM	Ti6AL4V-ELI	4.5x30mm
MCPPSS-45X35	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 4.5X35 MM	Ti6AL4V-ELI	4.5x35mm
MCPPSS-45X40	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 4.5X40 MM	Ti6AL4V-ELI	4.5x40mm
MCPPSS-45X45	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 4.5X45 MM	Ti6AL4V-ELI	4.5x45mm
MCPPSS-45X50	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 4.5X50 MM	Ti6AL4V-ELI	4.5x50mm
MCPPSS-5X30	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 5X30 MM	Ti6AL4V-ELI	5x30mm
MCPPSS-5X35	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 5X35 MM	Ti6AL4V-ELI	5x35mm
MCPPSS-5X40	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 5X40 MM	Ti6AL4V-ELI	5x40mm
MCPPSS-5X45	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 5X45 MM	Ti6AL4V-ELI	5x45mm
MCPPSS-5X50	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 5X50 MM	Ti6AL4V-ELI	5x50mm
MCPPSS-5X55	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 5X55 MM	Ti6AL4V-ELI	5x55mm
MCPPSS-5X60	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 5X60 MM	Ti6AL4V-ELI	5x60mm
MCPPSS-5X65	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 5X65 MM	Ti6AL4V-ELI	5x65mm
MCPPSS-5.5X30	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 5.5X30 MM	Ti6AL4V-ELI	5.5x30mm
MCPPSS-5.5X35	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 5.5X35 MM	Ti6AL4V-ELI	5.5x35mm
MCPPSS-5.5X40	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 5.5X40 MM	Ti6AL4V-ELI	5.5x40mm
MCPPSS-5.5X45	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 5.5X45 MM	Ti6AL4V-ELI	5.5x45mm
MCPPSS-5.5X50	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 5.5X50 MM	Ti6AL4V-ELI	5.5x50mm
MCPPSS-5.5X55	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 5.5X55 MM	Ti6AL4V-ELI	5.5x55mm
MCPPSS-5.5X60	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 5.5X60 MM	Ti6AL4V-ELI	5.5x60mm
MCPPSS-5.5X65	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 5.5X65 MM	Ti6AL4V-ELI	5.5x65mm
MCPPSS-6X30	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 6X30 MM	Ti6AL4V-ELI	6x30mm
MCPPSS-6X35	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 6X35 MM	Ti6AL4V-ELI	6x35mm
MCPPSS-6X40	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 6X40 MM	Ti6AL4V-ELI	6x40mm

# Thoracolumbar Pedicular Screws




MCPPSS-6X45	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 6X45 MM	Ti6AL4V-ELI	6x45mm
MCPPSS-6X50	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 6X50 MM	Ti6AL4V-ELI	6x50mm
MCPPSS-6X55	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 6X55 MM	Ti6AL4V-ELI	6x55mm
MCPPSS-6X60	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 6X60 MM	Ti6AL4V-ELI	6x60mm
MCPPSS-6X65	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 6X65 MM	Ti6AL4V-ELI	6x65mm
MCPPSS-6.5X30	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 6.5X30 MM	Ti6AL4V-ELI	6.5x30mm
MCPPSS-6.5X35	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 6.5X35 MM	Ti6AL4V-ELI	6.5x35mm
MCPPSS-6.5X40	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 6.5X40 MM	Ti6AL4V-ELI	6.5x40mm
MCPPSS-6.5X45	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 6.5X45 MM	Ti6AL4V-ELI	6.5x45mm
MCPPSS-6.5X50	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 6.5X50 MM	Ti6AL4V-ELI	6.5x50mm
MCPPSS-6.5X55	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 6.5X55 MM	Ti6AL4V-ELI	6.5x55mm
MCPPSS-6.5X60	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 6.5X60 MM	Ti6AL4V-ELI	6.5x60mm
MCPPSS-6.5X65	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 6.5X65 MM	Ti6AL4V-ELI	6.5x65mm
MCPPSS-7X30	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 7X30 MM	Ti6AL4V-ELI	7x30mm
MCPPSS-7X35	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 7X35 MM	Ti6AL4V-ELI	7x35mm
MCPPSS-7X40	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 7X40 MM	Ti6AL4V-ELI	7x40mm
MCPPSS-7X45	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 7X45 MM	Ti6AL4V-ELI	7x45mm
MCPPSS-7X50	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 7X50 MM	Ti6AL4V-ELI	7x50mm
MCPPSS-7X55	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 7X55 MM	Ti6AL4V-ELI	7x55mm
MCPPSS-7X60	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 7X60 MM	Ti6AL4V-ELI	7x60mm
MCPPSS-7X65	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 7X65 MM	Ti6AL4V-ELI	7x65mm
MCPPSS-7.5X30	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 7.5X30 MM	Ti6AL4V-ELI	7.5x30mm
MCPPSS-7.5X35	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 7.5X35 MM	Ti6AL4V-ELI	7.5x35mm
MCPPSS-7.5X40	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 7.5X40 MM	Ti6AL4V-ELI	7.5x40mm
MCPPSS-7.5X45	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 7.5X45 MM	Ti6AL4V-ELI	7.5x45mm
MCPPSS-7.5X50	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 7.5X50 MM	Ti6AL4V-ELI	7.5x50mm

# Thoracolumbar

## Pedicular Screws



MCPSS-7.5X55	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 7.5X55 MM	Ti6AL4V-ELI	7.5x55mm
MCPSS-7.5X50	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 7.5X60 MM	Ti6AL4V-ELI	7.5x60mm
MCPSS-7.5X65	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 7.5X65 MM	Ti6AL4V-ELI	7.5x65mm
MCPSS-8X30	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 8X30 MM	Ti6AL4V-ELI	8x30mm
MCPSS-8X35	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 8X35 MM	Ti6AL4V-ELI	8x35mm
MCPSS-8X40	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 8X40 MM	Ti6AL4V-ELI	8x40mm
MCPSS-8X45	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 8X45 MM	Ti6AL4V-ELI	8x45mm
MCPSS-8X50	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 8X50 MM	Ti6AL4V-ELI	8x50mm
MCPSS-8X55	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 8X55 MM	Ti6AL4V-ELI	8x55mm
MCPSS-8X60	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 8X60 MM	Ti6AL4V-ELI	8x60mm
MCPSS-8X65	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 8X65 MM	Ti6AL4V-ELI	8x65mm
MCPSS-8X70	MC POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW 8X70 MM	Ti6AL4V-ELI	8x70mm



# Thoracolumbar Pedicular Screws



MCPPS-4.5X25	MC POSTERIOR POLYAXIAL SCREW 4.5X25 MM	TI6AL4V-ELI	4.5x25mm
MCPPS-4.5X30	MC POSTERIOR POLYAXIAL SCREW 4.5X30 MM	TI6AL4V-ELI	4.5x30mm
MCPPS-4.5X35	MC POSTERIOR POLYAXIAL SCREW 4.5X35 MM	TI6AL4V-ELI	4.5x35mm
MCPPS-4.5X40	MC POSTERIOR POLYAXIAL SCREW 4.5X40 MM	TI6AL4V-ELI	4.5x40mm
MCPPS-4.5X45	MC POSTERIOR POLYAXIAL SCREW 4.5X45 MM	TI6AL4V-ELI	4.5x45mm
MCPPS-4.5X50	MC POSTERIOR POLYAXIAL SCREW 4.5X50 MM	TI6AL4V-ELI	4.5x50mm
MCPPS-5X30	MC POSTERIOR POLYAXIAL SCREW 5X30 MM	TI6AL4V-ELI	5x30mm
MCPPS-5X35	MC POSTERIOR POLYAXIAL SCREW 5X35 MM	TI6AL4V-ELI	5x35mm
MCPPS-5X40	MC POSTERIOR POLYAXIAL SCREW 5X40 MM	TI6AL4V-ELI	5x40mm
MCPPS-5X45	MC POSTERIOR POLYAXIAL SCREW 5X45 MM	TI6AL4V-ELI	5x45mm
MCPPS-5X50	MC POSTERIOR POLYAXIAL SCREW 5X50 MM	TI6AL4V-ELI	5x50mm
MCPPS-5X55	MC POSTERIOR POLYAXIAL SCREW 5X55 MM	TI6AL4V-ELI	5x55mm
MCPPS-5X60	MC POSTERIOR POLYAXIAL SCREW 5X60 MM	TI6AL4V-ELI	5x60mm
MCPPS-5X65	MC POSTERIOR POLYAXIAL SCREW 5X65 MM	TI6AL4V-ELI	5x65mm
MCPPS-5.5X30	MC POSTERIOR POLYAXIAL SCREW 5.5X30 MM	TI6AL4V-ELI	5.5x30mm
MCPPS-5.5X35	MC POSTERIOR POLYAXIAL SCREW 5.5X35 MM	TI6AL4V-ELI	5.5x35mm
MCPPS-5.5X40	MC POSTERIOR POLYAXIAL SCREW 5.5X40 MM	TI6AL4V-ELI	5.5x40mm
MCPPS-5.5X45	MC POSTERIOR POLYAXIAL SCREW 5.5X45 MM	TI6AL4V-ELI	5.5x45mm
MCPPS-5.5X50	MC POSTERIOR POLYAXIAL SCREW 5.5X50 MM	TI6AL4V-ELI	5.5x50mm
MCPPS-5.5X55	MC POSTERIOR POLYAXIAL SCREW 5.5X55 MM	TI6AL4V-ELI	5.5x55mm
MCPPS-5.5X60	MC POSTERIOR POLYAXIAL SCREW 5.5X60 MM	TI6AL4V-ELI	5.5x60mm
MCPPS-5.5X65	MC POSTERIOR POLYAXIAL SCREW 5.5X65 MM	TI6AL4V-ELI	5.5x65mm
MCPPS-6X30	MC POSTERIOR POLYAXIAL SCREW 6X30 MM	TI6AL4V-ELI	6x30mm
MCPPS-6X35	MC POSTERIOR POLYAXIAL SCREW 6X35 MM	TI6AL4V-ELI	6x35mm
MCPPS-6X40	MC POSTERIOR POLYAXIAL SCREW 6X40 MM	TI6AL4V-ELI	6x40mm

# Thoracolumbar

## Pedicular Screws



MCPPS-6X45	MC POSTERIOR POLYAXIAL SCREW 6X45 MM	TI6AL4V-ELI	6x45mm
MCPPS-6X50	MC POSTERIOR POLYAXIAL SCREW 6X50 MM	TI6AL4V-ELI	6x50mm
MCPPS-6X55	MC POSTERIOR POLYAXIAL SCREW 6X55 MM	TI6AL4V-ELI	6x55mm
MCPPS-6X60	MC POSTERIOR POLYAXIAL SCREW 6X60 MM	TI6AL4V-ELI	6x60mm
MCPPS-6X65	MC POSTERIOR POLYAXIAL SCREW 6X65 MM	TI6AL4V-ELI	6x65mm
MCPPS-6.5X30	MC POSTERIOR POLYAXIAL SCREW 6.5X30 MM	TI6AL4V-ELI	6.5x30mm
MCPPS-6.5X35	MC POSTERIOR POLYAXIAL SCREW 6.5X35 MM	TI6AL4V-ELI	6.5x35mm
MCPPS-6.5X40	MC POSTERIOR POLYAXIAL SCREW 6.5X40 MM	TI6AL4V-ELI	6.5x40mm
MCPPS-6.5X45	MC POSTERIOR POLYAXIAL SCREW 6.5X45 MM	TI6AL4V-ELI	6.5x45mm
MCPPS-6.5X50	MC POSTERIOR POLYAXIAL SCREW 6.5X50 MM	TI6AL4V-ELI	6.5x50mm
MCPPS-6.5X55	MC POSTERIOR POLYAXIAL SCREW 6.5X55 MM	TI6AL4V-ELI	6.5x55mm
MCPPS-6.5X60	MC POSTERIOR POLYAXIAL SCREW 6.5X60 MM	TI6AL4V-ELI	6.5x60mm
MCPPS-6.5X65	MC POSTERIOR POLYAXIAL SCREW 6.5X65 MM	TI6AL4V-ELI	6.5x65mm
MCPPS-7X30	MC POSTERIOR POLYAXIAL SCREW 7X30 MM	TI6AL4V-ELI	7x30mm
MCPPS-7X35	MC POSTERIOR POLYAXIAL SCREW 7X35 MM	TI6AL4V-ELI	7x35mm
MCPPS-7X40	MC POSTERIOR POLYAXIAL SCREW 7X40 MM	TI6AL4V-ELI	7x40mm
MCPPS-7X45	MC POSTERIOR POLYAXIAL SCREW 7X45 MM	TI6AL4V-ELI	7x45mm
MCPPS-7X50	MC POSTERIOR POLYAXIAL SCREW 7X50 MM	TI6AL4V-ELI	7x50mm
MCPPS-7X55	MC POSTERIOR POLYAXIAL SCREW 7X55 MM	TI6AL4V-ELI	7x55mm
MCPPS-7X60	MC POSTERIOR POLYAXIAL SCREW 7X60 MM	TI6AL4V-ELI	7x60mm
MCPPS-7X65	MC POSTERIOR POLYAXIAL SCREW 7X65 MM	TI6AL4V-ELI	7x65mm
MCPPS-7.5X30	MC POSTERIOR POLYAXIAL SCREW 7.5X30 MM	TI6AL4V-ELI	7.5x30mm
MCPPS-7.5X35	MC POSTERIOR POLYAXIAL SCREW 7.5X35 MM	TI6AL4V-ELI	7.5x35mm
MCPPS-7.5X40	MC POSTERIOR POLYAXIAL SCREW 7.5X40 MM	TI6AL4V-ELI	7.5x40mm
MCPPS-7.5X45	MC POSTERIOR POLYAXIAL SCREW 7.5X45 MM	TI6AL4V-ELI	7.5x45mm

# Thoracolumbar Pedicular Screws



MCPPS-7.5X50	MC POSTERIOR POLYAXIAL SCREW 7.5X50 MM	Ti6AL4V-ELI	7.5x50mm
MCPPS-7.5X55	MC POSTERIOR POLYAXIAL SCREW 7.5X55 MM	Ti6AL4V-ELI	7.5x55mm
MCPPS-7.5X60	MC POSTERIOR POLYAXIAL SCREW 7.5X60 MM	Ti6AL4V-ELI	7.5x60mm
MCPPS-7.5X65	MC POSTERIOR POLYAXIAL SCREW 7.5X65 MM	Ti6AL4V-ELI	7.5x65mm
MCPPS-8X30	MC POSTERIOR POLYAXIAL SCREW 8X30 MM	Ti6AL4V-ELI	8x30mm
MCPPS-8X35	MC POSTERIOR POLYAXIAL SCREW 8X35 MM	Ti6AL4V-ELI	8x35mm
MCPPS-8X40	MC POSTERIOR POLYAXIAL SCREW 8X40 MM	Ti6AL4V-ELI	8x40mm
MCPPS-8X45	MC POSTERIOR POLYAXIAL SCREW 8X45 MM	Ti6AL4V-ELI	8x45mm
MCPPS-8X50	MC POSTERIOR POLYAXIAL SCREW 8X50 MM	Ti6AL4V-ELI	8x50mm
MCPPS-8X55	MC POSTERIOR POLYAXIAL SCREW 8X55 MM	Ti6AL4V-ELI	8x55mm
MCPPS-8X60	MC POSTERIOR POLYAXIAL SCREW 8X60 MM	Ti6AL4V-ELI	8x60mm
MCPPS-8X65	MC POSTERIOR POLYAXIAL SCREW 8X65 MM	Ti6AL4V-ELI	8x65mm



# Thoracolumbar

## Pedicular Screws



MCPFR-5.5X40	MC POSTERIOR FIXATION ROD 5.5X40MM	Ti6AL4V-ELI	5.5x40mm
MCPFR-5.5X50	MC POSTERIOR FIXATION ROD 5.5X50MM	Ti6AL4V-ELI	5.5x50mm
MCPFR-5.5X60	MC POSTERIOR FIXATION ROD 5.5X60MM	Ti6AL4V-ELI	5.5x60mm
MCPFR-5.5X70	MC POSTERIOR FIXATION ROD 5.5X70MM	Ti6AL4V-ELI	5.5x70mm
MCPFR-5.5X80	MC POSTERIOR FIXATION ROD 5.5X80MM	Ti6AL4V-ELI	5.5x80mm
MCPFR-5.5X90	MC POSTERIOR FIXATION ROD 5.5X90MM	Ti6AL4V-ELI	5.5x90mm
MCPFR-5.5X100	MC POSTERIOR FIXATION ROD 5.5X100MM	Ti6AL4V-ELI	5.5x100mm
MCPFR-5.5X110	MC POSTERIOR FIXATION ROD 5.5X110MM	Ti6AL4V-ELI	5.5x110mm
MCPFR-5.5X120	MC POSTERIOR FIXATION ROD 5.5X120MM	Ti6AL4V-ELI	5.5x120mm
MCPFR-5.5X130	MC POSTERIOR FIXATION ROD 5.5X130MM	Ti6AL4V-ELI	5.5x130mm
MCPFR-5.5X140	MC POSTERIOR FIXATION ROD 5.5X140MM	Ti6AL4V-ELI	5.5x140mm
MCPFR-5.5X150	MC POSTERIOR FIXATION ROD 5.5X150MM	Ti6AL4V-ELI	5.5x150mm
MCPFR-5.5X160	MC POSTERIOR FIXATION ROD 5.5X160MM	Ti6AL4V-ELI	5.5x160mm
MCPFR-5.5X170	MC POSTERIOR FIXATION ROD 5.5X170MM	Ti6AL4V-ELI	5.5x170mm
MCPFR-5.5X180	MC POSTERIOR FIXATION ROD 5.5X180MM	Ti6AL4V-ELI	5.5x180mm
MCPFR-5.5X190	MC POSTERIOR FIXATION ROD 5.5X190MM	Ti6AL4V-ELI	5.5x190mm
MCPFR-5.5X200	MC POSTERIOR FIXATION ROD 5.5X200MM	Ti6AL4V-ELI	5.5x200mm
MCPFR-5.5X210	MC POSTERIOR FIXATION ROD 5.5X210MM	Ti6AL4V-ELI	5.5x210mm
MCPFR-5.5X220	MC POSTERIOR FIXATION ROD 5.5X220MM	Ti6AL4V-ELI	5.5x220mm
MCPFR-5.5X230	MC POSTERIOR FIXATION ROD 5.5X230MM	Ti6AL4V-ELI	5.5x230mm
MCPFR-5.5X240	MC POSTERIOR FIXATION ROD 5.5X240MM	Ti6AL4V-ELI	5.5x240mm
MCPFR-5.5X250	MC POSTERIOR FIXATION ROD 5.5X250MM	Ti6AL4V-ELI	5.5x250mm
MCPFR-5.5X260	MC POSTERIOR FIXATION ROD 5.5X260MM	Ti6AL4V-ELI	5.5x260mm
MCPFR-5.5X270	MC POSTERIOR FIXATION ROD 5.5X270MM	Ti6AL4V-ELI	5.5x270mm
MCPFR-5.5X280	MC POSTERIOR FIXATION ROD 5.5X280MM	Ti6AL4V-ELI	5.5x280mm
MCPFR-5.5X290	MC POSTERIOR FIXATION ROD 5.5X290MM	Ti6AL4V-ELI	5.5x290mm

# Thoracolumbar

## Pedicular Screws



MCPFR-5.5X300	MC POSTERIOR FIXATION ROD 5.5X300MM	Ti6AL4V-ELI	5.5x300mm
MCPFR-5.5X310	MC POSTERIOR FIXATION ROD 5.5X310MM	Ti6AL4V-ELI	5.5x310mm
MCPFR-5.5X320	MC POSTERIOR FIXATION ROD 5.5X320MM	Ti6AL4V-ELI	5.5x320mm
MCPFR-5.5X330	MC POSTERIOR FIXATION ROD 5.5X330MM	Ti6AL4V-ELI	5.5x330mm
MCPFR-5.5X340	MC POSTERIOR FIXATION ROD 5.5X340MM	Ti6AL4V-ELI	5.5x340mm
MCPFR-5.5X350	MC POSTERIOR FIXATION ROD 5.5X350MM	Ti6AL4V-ELI	5.5x350mm
MCPFR-5.5X360	MC POSTERIOR FIXATION ROD 5.5X360MM	Ti6AL4V-ELI	5.5x360mm
MCPFR-5.5X370	MC POSTERIOR FIXATION ROD 5.5X370MM	Ti6AL4V-ELI	5.5x370mm
MCPFR-5.5X380	MC POSTERIOR FIXATION ROD 5.5X380MM	Ti6AL4V-ELI	5.5x380mm
MCPFR-5.5X390	MC POSTERIOR FIXATION ROD 5.5X390MM	Ti6AL4V-ELI	5.5x390mm
MCPFR-5.5X400	MC POSTERIOR FIXATION ROD 5.5X400MM	Ti6AL4V-ELI	5.5x400mm
MCPFR-5.5X410	MC POSTERIOR FIXATION ROD 5.5X410MM	Ti6AL4V-ELI	5.5x410mm
MCPFR-5.5X420	MC POSTERIOR FIXATION ROD 5.5X420MM	Ti6AL4V-ELI	5.5x420mm
MCPFR-5.5X430	MC POSTERIOR FIXATION ROD 5.5X430MM	Ti6AL4V-ELI	5.5x430mm
MCPFR-5.5X440	MC POSTERIOR FIXATION ROD 5.5X440MM	Ti6AL4V-ELI	5.5x440mm
MCPFR-5.5X450	MC POSTERIOR FIXATION ROD 5.5X450MM	Ti6AL4V-ELI	5.5x450mm
MCPFR-5.5X460	MC POSTERIOR FIXATION ROD 5.5X460MM	Ti6AL4V-ELI	5.5x460mm
MCPFR-5.5X470	MC POSTERIOR FIXATION ROD 5.5X470MM	Ti6AL4V-ELI	5.5x470mm
MCPFR-5.5X480	MC POSTERIOR FIXATION ROD 5.5X480MM	Ti6AL4V-ELI	5.5x480mm
MCPFR-5.5X490	MC POSTERIOR FIXATION ROD 5.5X490MM	Ti6AL4V-ELI	5.5x490mm
MCPFR-5.5X500	MC POSTERIOR FIXATION ROD 5.5X500MM	Ti6AL4V-ELI	5.5x500mm
MCPFR-5.5X600	MC POSTERIOR FIXATION ROD 5.5X600MM	Ti6AL4V-ELI	5.5x600mm
MCR-6X40	MC POSTERIOR FIXATION ROD 6X40 MM	Ti6AL4V-ELI	6x40mm
MCR-6X50	MC POSTERIOR FIXATION ROD 6X50 MM	Ti6AL4V-ELI	6x50mm
MCR-6X60	MC POSTERIOR FIXATION ROD 6X60 MM	Ti6AL4V-ELI	6x60mm
MCR-6X70	MC POSTERIOR FIXATION ROD 6X70 MM	Ti6AL4V-ELI	6x70mm

# Thoracolumbar

## Pedicular Screws



MCR-6X80	MC POSTERIOR FIXATION ROD 6X80 MM	Ti6AL4V-ELI	6x80mm
MCR-6X90	MC POSTERIOR FIXATION ROD 6X90 MM	Ti6AL4V-ELI	6x90mm
MCR-6X100	MC POSTERIOR FIXATION ROD 6X100 MM	Ti6AL4V-ELI	6x100mm
MCR-6X110	MC POSTERIOR FIXATION ROD 6X110 MM	Ti6AL4V-ELI	6x110mm
MCR-6X120	MC POSTERIOR FIXATION ROD 6X120 MM	Ti6AL4V-ELI	6x120mm
MCR-6X130	MC POSTERIOR FIXATION ROD 6X130 MM	Ti6AL4V-ELI	6x130mm
MCR-6X140	MC POSTERIOR FIXATION ROD 6X140 MM	Ti6AL4V-ELI	6x140mm
MCR-6X150	MC POSTERIOR FIXATION ROD 6X150 MM	Ti6AL4V-ELI	6x150mm
MCR-6X160	MC POSTERIOR FIXATION ROD 6X160 MM	Ti6AL4V-ELI	6x160mm
MCR-6X170	MC POSTERIOR FIXATION ROD 6X170 MM	Ti6AL4V-ELI	6x170mm
MCR-6X180	MC POSTERIOR FIXATION ROD 6X180 MM	Ti6AL4V-ELI	6x180mm
MCR-6X190	MC POSTERIOR FIXATION ROD 6X190 MM	Ti6AL4V-ELI	6x190mm
MCR-6X200	MC POSTERIOR FIXATION ROD 6X200 MM	Ti6AL4V-ELI	6x200mm
MCR-6X210	MC POSTERIOR FIXATION ROD 6X210 MM	Ti6AL4V-ELI	6x220mm
MCR-6X220	MC POSTERIOR FIXATION ROD 6X220 MM	Ti6AL4V-ELI	6x230mm
MCR-6X230	MC POSTERIOR FIXATION ROD 6X230 MM	Ti6AL4V-ELI	6x240mm
MCR-6X240	MC POSTERIOR FIXATION ROD 6X240 MM	Ti6AL4V-ELI	6x240mm
MCR-6X250	MC POSTERIOR FIXATION ROD 6X250 MM	Ti6AL4V-ELI	6x250mm
MCR-6X260	MC POSTERIOR FIXATION ROD 6X260 MM	Ti6AL4V-ELI	6x260mm
MCR-6X270	MC POSTERIOR FIXATION ROD 6X270 MM	Ti6AL4V-ELI	6x270mm
MCR-6X280	MC POSTERIOR FIXATION ROD 6X280 MM	Ti6AL4V-ELI	6x280mm
MCR-6X290	MC POSTERIOR FIXATION ROD 6X290 MM	Ti6AL4V-ELI	6x290mm
MCR-6X300	MC POSTERIOR FIXATION ROD 6X300 MM	Ti6AL4V-ELI	6x300mm
MCR-6X310	MC POSTERIOR FIXATION ROD 6X310 MM	Ti6AL4V-ELI	6x310mm
MCR-6X320	MC POSTERIOR FIXATION ROD 6X320 MM	Ti6AL4V-ELI	6x320mm
MCR-6X330	MC POSTERIOR FIXATION ROD 6X330 MM	Ti6AL4V-ELI	6x330mm

# Thoracolumbar

## Pedicular Screws



MCR-6X340	MC POSTERIOR FIXATION ROD 6X340 MM	Ti6AL4V-ELI	6x340mm
MCR-6X350	MC POSTERIOR FIXATION ROD 6X350 MM	Ti6AL4V-ELI	6x350mm
MCR-6X360	MC POSTERIOR FIXATION ROD 6X360 MM	Ti6AL4V-ELI	6x360mm
MCR-6X370	MC POSTERIOR FIXATION ROD 6X370 MM	Ti6AL4V-ELI	6x370mm
MCR-6X380	MC POSTERIOR FIXATION ROD 6X380 MM	Ti6AL4V-ELI	6x380mm
MCR-6X390	MC POSTERIOR FIXATION ROD 6X390 MM	Ti6AL4V-ELI	6x390mm
MCR-6X400	MC POSTERIOR FIXATION ROD 6X400 MM	Ti6AL4V-ELI	6x400mm
MCR-6X410	MC POSTERIOR FIXATION ROD 6X410 MM	Ti6AL4V-ELI	6x410mm
MCR-6X420	MC POSTERIOR FIXATION ROD 6X420 MM	Ti6AL4V-ELI	6x420mm
MCR-6X430	MC POSTERIOR FIXATION ROD 6X430 MM	Ti6AL4V-ELI	6x430mm
MCR-6X440	MC POSTERIOR FIXATION ROD 6X440 MM	Ti6AL4V-ELI	6x440mm
MCR-6X450	MC POSTERIOR FIXATION ROD 6X450 MM	Ti6AL4V-ELI	6x450mm
MCR-6X460	MC POSTERIOR FIXATION ROD 6X460 MM	Ti6AL4V-ELI	6x460mm
MCR-6X470	MC POSTERIOR FIXATION ROD 6X470 MM	Ti6AL4V-ELI	6x470mm
MCR-6X480	MC POSTERIOR FIXATION ROD 6X480 MM	Ti6AL4V-ELI	6x480mm
MCR-6X490	MC POSTERIOR FIXATION ROD 6X490 MM	Ti6AL4V-ELI	6x490mm
MCR-6X500	MC POSTERIOR FIXATION ROD 6X500 MM	Ti6AL4V-ELI	6x500mm
MCR-6X600	MC POSTERIOR FIXATION ROD 6X600 MM	Ti6AL4V-ELI	6x510mm
MCR-6.30X70	MC POSTERIOR FIXATION ROD 6.30X70 MM	Ti6AL4V-ELI	6.30x70mm
MCR-6.30X100	MC POSTERIOR FIXATION ROD 6.30X100 MM	Ti6AL4V-ELI	6.30x100mm
MCR-6.30x150	MC POSTERIOR FIXATION ROD 6.30X150 MM	Ti6AL4V-ELI	6.30x150mm
MCR-6.30X250	MC POSTERIOR FIXATION ROD 6.30X250 MM	Ti6AL4V-ELI	6.30x250mm
MCR-6.5X70	MC POSTERIOR FIXATION ROD 6.5X70 MM	Ti6AL4V-ELI	6.5x70mm
MCR-6.5X100	MC POSTERIOR FIXATION ROD 6.5X100 MM	Ti6AL4V-ELI	6.5x100mm
MCR-6.5X150	MC POSTERIOR FIXATION ROD 6.5X150 MM	Ti6AL4V-ELI	6.5x150mm
MCR-6.5X250	MC POSTERIOR FIXATION ROD 6.5X250 MM	Ti6AL4V-ELI	6.5x250mm

# Thoracolumbar

## Pedicular Screws



MCPPR-5.5X40	MC POSTERIOR PERKÜTAN ROD 5.5X40MM	T16AL4V-ELI	5.5x40mm
MCPPR-5.5X50	MC POSTERIOR PERKÜTAN ROD 5.5X50MM	T16AL4V-ELI	5.5x50mm
MCPPR-5.5X60	MC POSTERIOR PERKÜTAN ROD 5.5X60MM	T16AL4V-ELI	5.5x60mm
MCPPR-5.5X70	MC POSTERIOR PERKÜTAN ROD 5.5X70MM	T16AL4V-ELI	5.5x70mm
MCPPR-5.5X80	MC POSTERIOR PERKÜTAN ROD 5.5X80MM	T16AL4V-ELI	5.5x80mm
MCPPR-5.5X90	MC POSTERIOR PERKÜTAN ROD 5.5X90MM	T16AL4V-ELI	5.5x90mm
MCPPR-5.5X100	MC POSTERIOR PERKÜTAN ROD 5.5X100MM	T16AL4V-ELI	5.5x100mm
MCPPR-5.5X110	MC POSTERIOR PERKÜTAN ROD 5.5X110MM	T16AL4V-ELI	5.5x110mm
MCPPR-5.5X120	MC POSTERIOR PERKÜTAN ROD 5.5X120MM	T16AL4V-ELI	5.5x120mm
MCPPR-5.5X140	MC POSTERIOR PERKÜTAN ROD 5.5X140MM	T16AL4V-ELI	5.5x140mm

# Thoracolumbar

## Pedicular Screws



MCPEOT-26-32MM	MC POSTERIOR EXPANDABLE OFFSETT TRANSVERS 26-32MM	Ti6AL4V-ELI	26-32mm
MCPEOT-32-40MM	MC POSTERIOR EXPANDABLE OFFSETT TRANSVERS 32-40MM	Ti6AL4V-ELI	32-40mm
MCPEOT-40-46MM	MC POSTERIOR EXPANDABLE OFFSETT TRANSVERS 40-46MM	Ti6AL4V-ELI	40-46mm
MCPEOT-40-50MM	MC POSTERIOR EXPANDABLE OFFSETT TRANSVERS 40-50MM	Ti6AL4V-ELI	40-50mm
MCPEOT-60-70MM	MC POSTERIOR EXPANDABLE OFFSETT TRANSVERS 60-70MM	Ti6AL4V-ELI	60-70mm
MCPEOT-50-60MM	MC POSTERIOR EXPANDABLE OFFSETT TRANSVERS 50-60MM	Ti6AL4V-ELI	50-60mm
MCPEOT-80-90MM	MC POSTERIOR EXPANDABLE OFFSETT TRANSVERS 80-90MM	Ti6AL4V-ELI	80-90mm



MCPTC-32-40	MC POSTERIOR TRANSVERS CONNECTOR 32-40MM	Ti6AL4V-ELI	32-40 mm
MCPTC-50-60	MC POSTERIOR TRANSVERS CONNECTOR 50-60MM	Ti6AL4V-ELI	50-60mm
MCPTC-40-50	MC POSTERIOR TRANSVERS CONNECTOR 40-50MM	Ti6AL4V-ELI	40-50mm
MCPTC-60-70	MC POSTERIOR TRANSVERS CONNECTOR 60-70MM	Ti6AL4V-ELI	60-70mm
MCPTC-80-90	MC POSTERIOR TRANSVERS CONNECTOR 80-90MM	Ti6AL4V-ELI	80-90mm

# Thoracolumbar

## Pedicular Screws



MCITMSS-45x25	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 4.5X25	TI6AL4V-ELI	4.5x25mm
MCITMSS-45x30	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 4.5X30	TI6AL4V-ELI	4.5x30mm
MCITMSS-45x35	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 4.5X35	TI6AL4V-ELI	4.5x35mm
MCITMSS-45x40	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 4.5X40	TI6AL4V-ELI	4.5x40mm
MCITMSS-45x45	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 4.5X45	TI6AL4V-ELI	4.5x45mm
MCITMSS-45x50	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 4.5X50	TI6AL4V-ELI	4.5x50mm
MCITMSS-5x30	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 5X30	TI6AL4V-ELI	5x30mm
MCITMSS-5x40	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 5X35	TI6AL4V-ELI	5x35mm
MCITMSS-5x40	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 5X40	TI6AL4V-ELI	5x40mm
MCITMSS-5x45	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 5X45	TI6AL4V-ELI	5x45mm
MCITMSS-5x50	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 5X50	TI6AL4V-ELI	5x50mm
MCITMSS-5x55	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 5X55	TI6AL4V-ELI	5x55mm
MCITMSS-5x60	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 5X60	TI6AL4V-ELI	5x60mm
MCITMSS-5x65	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 5X65	TI6AL4V-ELI	5x65mm
MCITMSS-5.5x30	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 5.5X30	TI6AL4V-ELI	5.5x30mm
MCITMSS-5.5x35	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 5.5X35	TI6AL4V-ELI	5.5x35mm
MCITMSS-5.5x40	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 5.5X40	TI6AL4V-ELI	5.5x40mm
MCITMSS-5.5x45	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 5.5X45	TI6AL4V-ELI	5.5x45mm
MCITMSS-5.5x50	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 5.5X50	TI6AL4V-ELI	5.5x50mm
MCITMSS-5.5x55	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 5.5X55	TI6AL4V-ELI	5.5x55mm
MCITMSS-5.5x60	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 5.5X60	TI6AL4V-ELI	5.5x60mm
MCITMSS-5.5x65	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 5.5X65	TI6AL4V-ELI	5.5x65mm
MCITMSS-6x30	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 6X30	TI6AL4V-ELI	6x30mm
MCITMSS-6x35	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 6X35	TI6AL4V-ELI	6x35mm
MCITMSS-6x40	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 6X40	TI6AL4V-ELI	6x40mm

# Thoracolumbar Pedicular Screws



MCITMSS-6x45	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 6X45	TI6AL4V-ELI	6x45mm
MCITMSS-6x50	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 6X50	TI6AL4V-ELI	6x50mm
MCITMSS-6x55	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 6X55	TI6AL4V-ELI	6x55mm
MCITMSS-6x60	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 6X60	TI6AL4V-ELI	6x60mm
MCITMSS-6x65	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 6X65	TI6AL4V-ELI	6x65mm
MCITMSS-6.5x30	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 6.5X30	TI6AL4V-ELI	6.5x30mm
MCITMSS-6.5x35	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 6.5X35	TI6AL4V-ELI	6.5x35mm
MCITMSS-6.5x40	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 6.5X40	TI6AL4V-ELI	6.5x40mm
MCITMSS-6.5x45	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 6.5X45	TI6AL4V-ELI	6.5x45mm
MCITMSS-6.5x50	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 6.5X50	TI6AL4V-ELI	6.5x50mm
MCITMSS-6.5x55	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 6.5X55	TI6AL4V-ELI	6.5x55mm
MCITMSS-6.5x60	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 6.5X60	TI6AL4V-ELI	6.5x60mm
MCITMSS-6.5x65	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 6.5X65	TI6AL4V-ELI	6.5x65mm
MCITMSS-7x30	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 7X30	TI6AL4V-ELI	7x30mm
MCITMSS-7x35	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 7X35	TI6AL4V-ELI	7x35mm
MCITMSS-7x40	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 7X40	TI6AL4V-ELI	7x40mm
MCITMSS-7x45	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 7X45	TI6AL4V-ELI	7x45mm
MCITMSS-7x50	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 7X50	TI6AL4V-ELI	7x50mm
MCITMSS-7x55	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 7X55	TI6AL4V-ELI	7x55mm
MCITMSS-7x60	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 7X60	TI6AL4V-ELI	7x60mm
MCITMSS-7x65	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 7X65	TI6AL4V-ELI	7x65mm
MCITMSS-7.5x30	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 7.5X30	TI6AL4V-ELI	7.5x30mm
MCITMSS-7.5x35	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 7.5X35	TI6AL4V-ELI	7.5x35mm
MCITMSS-7.5x40	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 7.5X40	TI6AL4V-ELI	7.5x40mm
MCITMSS-7.5x45	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 7.5X45	TI6AL4V-ELI	7.5x45mm

# Thoracolumbar

## Pedicular Screws



MCITMSS-7.5x50	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 7.5X50	TI6AL4V-ELI	7.5x50mm
MCITMSS-7.5x55	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 7.5X55	TI6AL4V-ELI	7.5x55mm
MCITMSS-7.5x60	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 7.5X60	TI6AL4V-ELI	7.5x60mm
MCITMSS-7.5x65	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 7.5X65	TI6AL4V-ELI	7.5x65mm
MCITMSS-8x30	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 8X30	TI6AL4V-ELI	8x30mm
MCITMSS-8x35	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 8X35	TI6AL4V-ELI	8x35mm
MCITMSS-8x40	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 8X40	TI6AL4V-ELI	8x40mm
MCITMSS-8x45	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 8X45	TI6AL4V-ELI	8x45mm
MCITMSS-8x50	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 8X50	TI6AL4V-ELI	8x50mm
MCITMSS-8x55	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 8X55	TI6AL4V-ELI	8x55mm
MCITMSS-8x60	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 8X60	TI6AL4V-ELI	8x60mm
MCITMSS-8x65	MC I TYPE MONOAXIAL SPONDYLOLISTHESIS SCREW 8X65	TI6AL4V-ELI	8x65mm

# Thoracolumbar Pedicular Screws



MCITPSS-45x25	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 4.5X25	TI6AL4V-ELI	4.5x25mm
MCITPSS-45x30	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 4.5X30	TI6AL4V-ELI	4.5x30mm
MCITPSS-45x35	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 4.5X35	TI6AL4V-ELI	4.5x35mm
MCITPSS-45x40	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 4.5X40	TI6AL4V-ELI	4.5x40mm
MCITPSS-45x45	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 4.5X45	TI6AL4V-ELI	4.5x45mm
MCITPSS-45x50	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 4.5X50	TI6AL4V-ELI	4.5x50mm
MCITPSS-5x30	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 5X30	TI6AL4V-ELI	5x30mm
MCITPSS-5x35	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 5X35	TI6AL4V-ELI	5x35mm
MCITPSS-5x40	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 5X40	TI6AL4V-ELI	5x40mm
MCITPSS-5x45	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 5X45	TI6AL4V-ELI	5x45mm
MCITPSS-5x50	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 5X50	TI6AL4V-ELI	5x50mm
MCITPSS-5x55	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 5X55	TI6AL4V-ELI	5x55mm
MCITPSS-5x60	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 5X60	TI6AL4V-ELI	5x60mm
MCITPSS-5x65	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 5X65	TI6AL4V-ELI	5x65mm
MCITPSS-5.5x30	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 5.5X30	TI6AL4V-ELI	5.5x30mm
MCITPSS-5.5x35	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 5.5X35	TI6AL4V-ELI	5.5x35mm
MCITPSS-5.5x40	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 5.5X40	TI6AL4V-ELI	5.5x40mm
MCITPSS-5.5x45	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 5.5X45	TI6AL4V-ELI	5.5x45mm
MCITPSS-5.5x50	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 5.5X50	TI6AL4V-ELI	5.5x50mm
MCITPSS-5.5x55	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 5.5X55	TI6AL4V-ELI	5.5x55mm
MCITPSS-5.5x60	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 5.5X60	TI6AL4V-ELI	5.5x60mm
MCITPSS-5.5x65	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 5.5X65	TI6AL4V-ELI	5.5x65mm
MCITPSS-6x30	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 6X30	TI6AL4V-ELI	6x30mm
MCITPSS-6x35	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 6X35	TI6AL4V-ELI	6x35mm
MCITPSS-6x40	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 6X40	TI6AL4V-ELI	6x40mm

# Thoracolumbar

## Pedicular Screws



MCITPSS-6x45	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 6X45	TI6AL4V-ELI	6x45mm
MCITPSS-6x50	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 6X50	TI6AL4V-ELI	6x50mm
MCITPSS-6x55	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 6X55	TI6AL4V-ELI	6x55mm
MCITPSS-6x60	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 6X60	TI6AL4V-ELI	6x60mm
MCITPSS-6x65	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 6X65	TI6AL4V-ELI	6x65mm
MCITPSS-6.5x30	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 6.5X30	TI6AL4V-ELI	6.5x30mm
MCITPSS-6.5x35	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 6.5X35	TI6AL4V-ELI	6.5x35mm
MCITPSS-6.5x40	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 6.5X40	TI6AL4V-ELI	6.5x40mm
MCITPSS-6.5x45	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 6.5X45	TI6AL4V-ELI	6.5x45mm
MCITPSS-6.5x50	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 6.5X50	TI6AL4V-ELI	6.5x50mm
MCITPSS-6.5x55	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 6.5X55	TI6AL4V-ELI	6.5x55mm
MCITPSS-6.5x60	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 6.5X60	TI6AL4V-ELI	6.5x60mm
MCITPSS-6.5x65	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 6.5X65	TI6AL4V-ELI	6.5x65mm
MCITPSS-7x30	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 7X30	TI6AL4V-ELI	7x30mm
MCITPSS-7x35	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 7X35	TI6AL4V-ELI	7x35mm
MCITPSS-7x40	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 7X40	TI6AL4V-ELI	7x40mm
MCITPSS-7x45	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 7X45	TI6AL4V-ELI	7x45mm
MCITPSS-7x50	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 7X50	TI6AL4V-ELI	7x50mm
MCITPSS-7x55	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 7X55	TI6AL4V-ELI	7x55mm
MCITPSS-7x60	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 7X60	TI6AL4V-ELI	7x60mm
MCITPSS-7x65	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 7X65	TI6AL4V-ELI	7x65mm
MCITPSS-7.5x30	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 7.5X30	TI6AL4V-ELI	7.5x30mm
MCITPSS-7.5x35	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 7.5X35	TI6AL4V-ELI	7.5x35mm
MCITPSS-7.5x40	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 7.5X40	TI6AL4V-ELI	7.5x40mm
MCITPSS-7.5x45	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 7.5X45	TI6AL4V-ELI	7.5x45mm

# Thoracolumbar Pedicular Screws



MCITPSS-7.5x50	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 7.5X50	TI6AL4V-ELI	7.5x50mm
MCITPSS-7.5x55	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 7.5X55	TI6AL4V-ELI	7.5x55mm
MCITPSS-7.5x60	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 7.5X60	TI6AL4V-ELI	7.5x60mm
MCITPSS-7.5x65	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 7.5X65	TI6AL4V-ELI	7.5x65mm
MCITPSS-8x30	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 8X30	TI6AL4V-ELI	8x30mm
MCITPSS-8x35	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 8X35	TI6AL4V-ELI	8x35mm
MCITPSS-8x40	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 8X40	TI6AL4V-ELI	8x40mm
MCITPSS-8x45	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 8X45	TI6AL4V-ELI	8x45mm
MCITPSS-8x50	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 8X50	TI6AL4V-ELI	8x50mm
MCITPSS-8x55	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 8X55	TI6AL4V-ELI	8x55mm
MCITPSS-8x60	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 8X60	TI6AL4V-ELI	8x60mm
MCITPSS-8x65	MC I TYPE POLYAXIAL SPONDYLOLISTHESIS SCREW 8X65	TI6AL4V-ELI	8x65mm

# Thoracolumbar

## Pedicular Screws



MCITMS-45x25	MC I TYPE MONOAXIAL SCREW 4.5X25	Ti6AL4V-ELI	4.5x25mm
MCITMS-45x30	MC I TYPE MONOAXIAL SCREW 4.5X30	Ti6AL4V-ELI	4.5x30mm
MCITMS-45x35	MC I TYPE MONOAXIAL SCREW 4.5X35	Ti6AL4V-ELI	4.5x35mm
MCITMS-45x40	MC I TYPE MONOAXIAL SCREW 4.5x40	Ti6AL4V-ELI	4.5x40mm
MCITMS-45x45	MC I TYPE MONOAXIAL SCREW 4.5x45	Ti6AL4V-ELI	4.5x45mm
MCITMS-45x50	MC I TYPE MONOAXIAL SCREW 4.5x50	Ti6AL4V-ELI	4.5x50mm
MCITMS-5x30	MC I TYPE MONOAXIAL SCREW 5X30	Ti6AL4V-ELI	5x30mm
MCITMS-5x35	MC I TYPE MONOAXIAL SCREW 5X35	Ti6AL4V-ELI	5x35mm
MCITMS-5x40	MC I TYPE MONOAXIAL SCREW 5X40	Ti6AL4V-ELI	5x40mm
MCITMS-5x45	MC I TYPE MONOAXIAL SCREW 5X45	Ti6AL4V-ELI	5x45mm
MCITMS-5x50	MC I TYPE MONOAXIAL SCREW 5X50	Ti6AL4V-ELI	5x50mm
MCITMS-5x55	MC I TYPE MONOAXIAL SCREW 5X55	Ti6AL4V-ELI	5x55mm
MCITMS-5x60	MC I TYPE MONOAXIAL SCREW 5X60	Ti6AL4V-ELI	5x60mm
MCITMS-5x65	MC I TYPE MONOAXIAL SCREW 5X65	Ti6AL4V-ELI	5x65mm
MCITMS-5.5x30	MC I TYPE MONOAXIAL SCREW 5.5X30	Ti6AL4V-ELI	5.5x30mm
MCITMS-5.5x35	MC I TYPE MONOAXIAL SCREW 5.5X35	Ti6AL4V-ELI	5.5x35mm
MCITMS-5.5x40	MC I TYPE MONOAXIAL SCREW 5.5x40	Ti6AL4V-ELI	5.5x40mm
MCITMS-5.5x45	MC I TYPE MONOAXIAL SCREW 5.5x45	Ti6AL4V-ELI	5.5x45mm
MCITMS-5.5x50	MC I TYPE MONOAXIAL SCREW 5.5x50	Ti6AL4V-ELI	5.5x50mm
MCITMS-5.5x55	MC I TYPE MONOAXIAL SCREW 5.5x55	Ti6AL4V-ELI	5.5x55mm
MCITMS-5.5x60	MC I TYPE MONOAXIAL SCREW 5.5x60	Ti6AL4V-ELI	5.5x60mm
MCITMS-5.5x65	MC I TYPE MONOAXIAL SCREW 5.5x65	Ti6AL4V-ELI	5.5x65mm
MCITMS-6x30	MC I TYPE MONOAXIAL SCREW 6X30	Ti6AL4V-ELI	6x30mm
MCITMS-6x35	MC I TYPE MONOAXIAL SCREW 6X35	Ti6AL4V-ELI	6x35mm
MCITMS-6x40	MC I TYPE MONOAXIAL SCREW 6X40	Ti6AL4V-ELI	6x40mm

# Thoracolumbar

## Pedicular Screws



MCITMS-6x45	MC I TYPE MONOAXIAL SCREW 6X45	TI6AL4V-ELI	6x45mm
MCITMS-6x50	MC I TYPE MONOAXIAL SCREW 6X50	TI6AL4V-ELI	6x50mm
MCITMS-6x55	MC I TYPE MONOAXIAL SCREW 6X55	TI6AL4V-ELI	6x55mm
MCITMS-6x60	MC I TYPE MONOAXIAL SCREW 6X60	TI6AL4V-ELI	6x60mm
MCITMS-6x65	MC I TYPE MONOAXIAL SCREW 6X65	TI6AL4V-ELI	6x65mm
MCITMS-6.5x30	MC I TYPE MONOAXIAL SCREW 6.5X30	TI6AL4V-ELI	6.5x30mm
MCITMS-6.5x35	MC I TYPE MONOAXIAL SCREW 6.5X35	TI6AL4V-ELI	6.5x35mm
MCITMS-6.5x40	MC I TYPE MONOAXIAL SCREW 6.5x40	TI6AL4V-ELI	6.5x40mm
MCITMS-6.5x45	MC I TYPE MONOAXIAL SCREW 6.5x45	TI6AL4V-ELI	6.5x45mm
MCITMS-6.5x50	MC I TYPE MONOAXIAL SCREW 6.5x50	TI6AL4V-ELI	6.5x50mm
MCITMS-6.5x55	MC I TYPE MONOAXIAL SCREW 6.5x55	TI6AL4V-ELI	6.5x55mm
MCITMS-6.5x60	MC I TYPE MONOAXIAL SCREW 6.5x60	TI6AL4V-ELI	6.5x60mm
MCITMS-6.5x65	MC I TYPE MONOAXIAL SCREW 6.5x65	TI6AL4V-ELI	6.5x65mm
MCITMS-7x30	MC I TYPE MONOAXIAL SCREW 7X30	TI6AL4V-ELI	7x30mm
MCITMS-7x35	MC I TYPE MONOAXIAL SCREW 7X35	TI6AL4V-ELI	7x35mm
MCITMS-7x40	MC I TYPE MONOAXIAL SCREW 7X40	TI6AL4V-ELI	7x40mm
MCITMS-7x45	MC I TYPE MONOAXIAL SCREW 7X45	TI6AL4V-ELI	7x45mm
MCITMS-7x50	MC I TYPE MONOAXIAL SCREW 7X50	TI6AL4V-ELI	7x50mm
MCITMS-7x55	MC I TYPE MONOAXIAL SCREW 7X55	TI6AL4V-ELI	7x55mm
MCITMS-7x60	MC I TYPE MONOAXIAL SCREW 7X60	TI6AL4V-ELI	7x60mm
MCITMS-7x65	MC I TYPE MONOAXIAL SCREW 7X65	TI6AL4V-ELI	7x65mm
MCITMS-7.5x30	MC I TYPE MONOAXIAL SCREW 7.5X30	TI6AL4V-ELI	7.5x30mm
MCITMS-7.5x35	MC I TYPE MONOAXIAL SCREW 7.5X35	TI6AL4V-ELI	7.5x35mm
MCITMS-7.5x40	MC I TYPE MONOAXIAL SCREW 7.5x40	TI6AL4V-ELI	7.5x40mm
MCITMS-7.5x45	MC I TYPE MONOAXIAL SCREW 7.5x45	TI6AL4V-ELI	7.5x45mm

# Thoracolumbar

## Pedicular Screws



MCITMS-7.5x50	MC I TYPE MONOAXIAL SCREW 7.5x50	Ti6AL4V-ELI	7.5x50mm
MCITMS-7.5x55	MC I TYPE MONOAXIAL SCREW 7.5x55	Ti6AL4V-ELI	7.5x55mm
MCITMS-7.5x60	MC I TYPE MONOAXIAL SCREW 7.5x60	Ti6AL4V-ELI	7.5x60mm
MCITMS-7.5x65	MC I TYPE MONOAXIAL SCREW 7.5x65	Ti6AL4V-ELI	7.5x65mm
MCITMS-8x30	MC I TYPE MONOAXIAL SCREW 8X30	Ti6AL4V-ELI	8x30mm
MCITMS-8x35	MC I TYPE MONOAXIAL SCREW 8X35	Ti6AL4V-ELI	8x35mm
MCITMS-8x40	MC I TYPE MONOAXIAL SCREW 8X40	Ti6AL4V-ELI	8x40mm
MCITMS-8x45	MC I TYPE MONOAXIAL SCREW 8X45	Ti6AL4V-ELI	8x45mm
MCITMS-8x50	MC I TYPE MONOAXIAL SCREW 8X50	Ti6AL4V-ELI	8x50mm
MCITMS-8x55	MC I TYPE MONOAXIAL SCREW 8X55	Ti6AL4V-ELI	8x55mm
MCITMS-8x60	MC I TYPE MONOAXIAL SCREW 8X60	Ti6AL4V-ELI	8x60mm
MCITMS-8x65	MC I TYPE MONOAXIAL SCREW 8X65	Ti6AL4V-ELI	8x65mm

# Thoracolumbar

## Pedicular Screws



MCITPS-45x25	MC I TYPE POLYAXIAL SCREW 4.5X25	Ti6AL4V-ELI	4.5x25mm
MCITPS-45x30	MC I TYPE POLYAXIAL SCREW 4.5X30	Ti6AL4V-ELI	4.5x30mm
MCITPS-45x35	MC I TYPE POLYAXIAL SCREW 4.5X35	Ti6AL4V-ELI	4.5x35mm
MCITPS-45x40	MC I TYPE POLYAXIAL SCREW 4.5x40	Ti6AL4V-ELI	4.5x40mm
MCITPS-45x45	MC I TYPE POLYAXIAL SCREW 4.5x45	Ti6AL4V-ELI	4.5x45mm
MCITPS-45x50	MC I TYPE POLYAXIAL SCREW 4.5x50	Ti6AL4V-ELI	4.5x50mm
MCITPS-5x30	MC I TYPE POLYAXIAL SCREW 5X30	Ti6AL4V-ELI	5x30mm
MCITPS-5x35	MC I TYPE POLYAXIAL SCREW 5X35	Ti6AL4V-ELI	5x35mm
MCITPS-5x40	MC I TYPE POLYAXIAL SCREW 5X40	Ti6AL4V-ELI	5x40mm
MCITPS-5x45	MC I TYPE POLYAXIAL SCREW 5X45	Ti6AL4V-ELI	5x45mm
MCITPS-5x50	MC I TYPE POLYAXIAL SCREW 5X50	Ti6AL4V-ELI	5x50mm
MCITPS-5x55	MC I TYPE POLYAXIAL SCREW 5X55	Ti6AL4V-ELI	5x55mm
MCITPS-5x60	MC I TYPE POLYAXIAL SCREW 5X60	Ti6AL4V-ELI	5x60mm
MCITPS-5x65	MC I TYPE POLYAXIAL SCREW 5X65	Ti6AL4V-ELI	5x65mm
MCITPS-5.5x30	MC I TYPE POLYAXIAL SCREW 5.5X30	Ti6AL4V-ELI	5.5x30mm
MCITPS-5.5x35	MC I TYPE POLYAXIAL SCREW 5.5X35	Ti6AL4V-ELI	5.5x35mm
MCITPS-5.5x40	MC I TYPE POLYAXIAL SCREW 5.5x40	Ti6AL4V-ELI	5.5x40mm
MCITPS-5.5x45	MC I TYPE POLYAXIAL SCREW 5.5x45	Ti6AL4V-ELI	5.5x45mm
MCITPS-5.5x50	MC I TYPE POLYAXIAL SCREW 5.5x50	Ti6AL4V-ELI	5.5x50mm
MCITPS-5.5x55	MC I TYPE POLYAXIAL SCREW 5.5X55	Ti6AL4V-ELI	5.5x55mm
MCITPS-5.5x60	MC I TYPE POLYAXIAL SCREW 5.5X60	Ti6AL4V-ELI	5.5x60mm
MCITPS-5.5x65	MC I TYPE POLYAXIAL SCREW 5.5X65	Ti6AL4V-ELI	5.5x65mm
MCITPS-6x30	MC I TYPE POLYAXIAL SCREW 6X30	Ti6AL4V-ELI	6x30mm
MCITPS-6x35	MC I TYPE POLYAXIAL SCREW 6X35	Ti6AL4V-ELI	6x35mm
MCITPS-6x40	MC I TYPE POLYAXIAL SCREW 6X40	Ti6AL4V-ELI	6x40mm

# Thoracolumbar

## Pedicular Screws



MCMVS-4x25	MC MULTIFONKSIYONEL VERTEBRAL SCREW 4X25 MM	TI6AL4V-ELI	4x25mm
MCMVS-4x30	MC MULTIFONKSIYONEL VERTEBRAL SCREW 4X30 MM	TI6AL4V-ELI	4x30mm
MCMVS-4x35	MC MULTIFONKSIYONEL VERTEBRAL SCREW 4X35 MM	TI6AL4V-ELI	4x35mm
MCMVS-4x40	MC MULTIFONKSIYONEL VERTEBRAL SCREW 4X40 MM	TI6AL4V-ELI	4x40
MCMVS-4x45	MC MULTIFONKSIYONEL VERTEBRAL SCREW 4X45 MM	TI6AL4V-ELI	4x45mm
MCMVS-4x50	MC MULTIFONKSIYONEL VERTEBRAL SCREW 4X50 MM	TI6AL4V-ELI	4x50mm
MCMVS-4.5x25	MC MULTIFONKSIYONEL VERTEBRAL SCREW 4.5X25 MM	TI6AL4V-ELI	4.5x25mm
MCMVS-4.5x30	MC MULTIFONKSIYONEL VERTEBRAL SCREW 4.5X30 MM	TI6AL4V-ELI	4.5x30mm
MCMVS-4.5x35	MC MULTIFONKSIYONEL VERTEBRAL SCREW 4.5X35 MM	TI6AL4V-ELI	4.5x35mm
MCMVS-4.5x40	MC MULTIFONKSIYONEL VERTEBRAL SCREW 4.5X40 MM	TI6AL4V-ELI	4.5x40mm
MCMVS-4.5x45	MC MULTIFONKSIYONEL VERTEBRAL SCREW 4.5X45 MM	TI6AL4V-ELI	4.5x45mm
MCMVS-4.5x50	MC MULTIFONKSIYONEL VERTEBRAL SCREW 4.5X50 MM	TI6AL4V-ELI	4.5x50mm
MCMVS-5x30	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5X30 MM	TI6AL4V-ELI	5x30mm
MCMVS-5x35	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5X35 MM	TI6AL4V-ELI	5x35mm
MCMVS-5x40	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5X40 MM	TI6AL4V-ELI	5x40mm
MCMVS-5x45	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5X45 MM	TI6AL4V-ELI	5x45mm
MCMVS-5x50	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5X50 MM	TI6AL4V-ELI	5x50mm

# Thoracolumbar

## Pedicular Screws



MCMVS-5x55	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5X55 MM	Ti6AL4V-ELI	5x55mm
MCMVS-5x60	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5X60 MM	Ti6AL4V-ELI	5x60mm
MCMVS-5x65	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5X65 MM	Ti6AL4V-ELI	5x65mm
MCMVS-5.5x30	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5.5X30 MM	Ti6AL4V-ELI	5.5x30mm
MCMVS-5.5x35	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5.5X35 MM	Ti6AL4V-ELI	5.5x35mm
MCMVS-5.5x40	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5.5X40 MM	Ti6AL4V-ELI	5.5x40mm
MCMVS-5.5x45	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5.5X45 MM	Ti6AL4V-ELI	5.5x45mm
MCMVS-5.5x50	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5.5X50 MM	Ti6AL4V-ELI	5.5x50mm
MCMVS-5.5x55	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5.5X55 MM	Ti6AL4V-ELI	5.5x55mm
MCMVS-5.5x60	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5.5X60 MM	Ti6AL4V-ELI	5.5x60mm
MCMVS-5.5x65	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5.5X65 MM	Ti6AL4V-ELI	5.5x65mm
MCMVS-6x30	MC MULTIFONKSIYONEL VERTEBRAL SCREW 6X30 MM	Ti6AL4V-ELI	6x30mm
MCMVS-6x35	MC MULTIFONKSIYONEL VERTEBRAL SCREW 6X35 MM	Ti6AL4V-ELI	6x35mm
MCMVS-6x40	MC MULTIFONKSIYONEL VERTEBRAL SCREW 6X40 MM	Ti6AL4V-ELI	6x40mm
MCMVS-6x45	MC MULTIFONKSIYONEL VERTEBRAL SCREW 6X45 MM	Ti6AL4V-ELI	6x45mm
MCMVS-6x50	MC MULTIFONKSIYONEL VERTEBRAL SCREW 6X50 MM	Ti6AL4V-ELI	6x50mm
MCMVS-6x55	MC MULTIFONKSIYONEL VERTEBRAL SCREW 6X55 MM	Ti6AL4V-ELI	6x55mm

# Thoracolumbar

## Pedicular Screws



MCMVS-4x25	MC MULTIFONKSIYONEL VERTEBRAL SCREW 4X25 MM	TI6AL4V-ELI	4x25mm
MCMVS-4x30	MC MULTIFONKSIYONEL VERTEBRAL SCREW 4X30 MM	TI6AL4V-ELI	4x30mm
MCMVS-4x35	MC MULTIFONKSIYONEL VERTEBRAL SCREW 4X35 MM	TI6AL4V-ELI	4x35mm
MCMVS-4x40	MC MULTIFONKSIYONEL VERTEBRAL SCREW 4X40 MM	TI6AL4V-ELI	4x40
MCMVS-4x45	MC MULTIFONKSIYONEL VERTEBRAL SCREW 4X45 MM	TI6AL4V-ELI	4x45mm
MCMVS-4x50	MC MULTIFONKSIYONEL VERTEBRAL SCREW 4X50 MM	TI6AL4V-ELI	4x50mm
MCMVS-4.5x25	MC MULTIFONKSIYONEL VERTEBRAL SCREW 4.5X25 MM	TI6AL4V-ELI	4.5x25mm
MCMVS-4.5x30	MC MULTIFONKSIYONEL VERTEBRAL SCREW 4.5X30 MM	TI6AL4V-ELI	4.5x30mm
MCMVS-4.5x35	MC MULTIFONKSIYONEL VERTEBRAL SCREW 4.5X35 MM	TI6AL4V-ELI	4.5x35mm
MCMVS-4.5x40	MC MULTIFONKSIYONEL VERTEBRAL SCREW 4.5X40 MM	TI6AL4V-ELI	4.5x40mm
MCMVS-4.5x45	MC MULTIFONKSIYONEL VERTEBRAL SCREW 4.5X45 MM	TI6AL4V-ELI	4.5x45mm
MCMVS-4.5x50	MC MULTIFONKSIYONEL VERTEBRAL SCREW 4.5X50 MM	TI6AL4V-ELI	4.5x50mm
MCMVS-5x30	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5X30 MM	TI6AL4V-ELI	5x30mm
MCMVS-5x35	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5X35 MM	TI6AL4V-ELI	5x35mm
MCMVS-5x40	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5X40 MM	TI6AL4V-ELI	5x40mm
MCMVS-5x45	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5X45 MM	TI6AL4V-ELI	5x45mm
MCMVS-5x50	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5X50 MM	TI6AL4V-ELI	5x50mm

# Thoracolumbar

## Pedicular Screws



MCMVS-5x55	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5X55 MM	TI6AL4V-ELI	5x55mm
MCMVS-5x60	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5X60 MM	TI6AL4V-ELI	5x60mm
MCMVS-5x65	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5X65 MM	TI6AL4V-ELI	5x65mm
MCMVS-5.5x30	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5.5X30 MM	TI6AL4V-ELI	5.5x30mm
MCMVS-5.5x35	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5.5X35 MM	TI6AL4V-ELI	5.5x35mm
MCMVS-5.5x40	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5.5X40 MM	TI6AL4V-ELI	5.5x40mm
MCMVS-5.5x45	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5.5X45 MM	TI6AL4V-ELI	5.5x45mm
MCMVS-5.5x50	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5.5X50 MM	TI6AL4V-ELI	5.5x50mm
MCMVS-5.5x55	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5.5X55 MM	TI6AL4V-ELI	5.5x55mm
MCMVS-5.5x60	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5.5X60 MM	TI6AL4V-ELI	5.5x60mm
MCMVS-5.5x65	MC MULTIFONKSIYONEL VERTEBRAL SCREW 5.5X65 MM	TI6AL4V-ELI	5.5x65mm
MCMVS-6x30	MC MULTIFONKSIYONEL VERTEBRAL SCREW 6X30 MM	TI6AL4V-ELI	6x30mm
MCMVS-6x35	MC MULTIFONKSIYONEL VERTEBRAL SCREW 6X35 MM	TI6AL4V-ELI	6x35mm
MCMVS-6x40	MC MULTIFONKSIYONEL VERTEBRAL SCREW 6X40 MM	TI6AL4V-ELI	6x40mm
MCMVS-6x45	MC MULTIFONKSIYONEL VERTEBRAL SCREW 6X45 MM	TI6AL4V-ELI	6x45mm
MCMVS-6x50	MC MULTIFONKSIYONEL VERTEBRAL SCREW 6X50 MM	TI6AL4V-ELI	6x50mm
MCMVS-6x55	MC MULTIFONKSIYONEL VERTEBRAL SCREW 6X55 MM	TI6AL4V-ELI	6x55mm

# Thoracolumbar

## Pedicular Screws



MCMVS-6x60	MC MULTIFONKSIYONEL VERTEBRAL SCREW 6X60 MM	TI6AL4V-ELI	6x60mm
MCMVS-6x65	MC MULTIFONKSIYONEL VERTEBRAL SCREW 6X65 MM	TI6AL4V-ELI	6x65mm
MCMVS-6.5x30	MC MULTIFONKSIYONEL VERTEBRAL SCREW 6.5X30 MM	TI6AL4V-ELI	6.5x30mm
MCMVS-6.5x35	MC MULTIFONKSIYONEL VERTEBRAL SCREW 6.5X35 MM	TI6AL4V-ELI	6.5x35mm
MCMVS-6.5x40	MC MULTIFONKSIYONEL VERTEBRAL SCREW 6.5X40 MM	TI6AL4V-ELI	6.5x40mm
MCMVS-6.5x45	MC MULTIFONKSIYONEL VERTEBRAL SCREW 6.5X45 MM	TI6AL4V-ELI	6.5x45mm
MCMVS-6.5x50	MC MULTIFONKSIYONEL VERTEBRAL SCREW 6.5X50 MM	TI6AL4V-ELI	6.5x50mm
MCMVS-6.5x55	MC MULTIFONKSIYONEL VERTEBRAL SCREW 6.5X55 MM	TI6AL4V-ELI	6.5x55mm
MCMVS-6.5x60	MC MULTIFONKSIYONEL VERTEBRAL SCREW 6.5X60 MM	TI6AL4V-ELI	6.5x60mm
MCMVS-6.5x65	MC MULTIFONKSIYONEL VERTEBRAL SCREW 6.5X65 MM	TI6AL4V-ELI	6.5x65mm
MCMVS-7x30	MC MULTIFONKSIYONEL VERTEBRAL SCREW 7X30 MM	TI6AL4V-ELI	7x30mm
MCMVS-7x35	MC MULTIFONKSIYONEL VERTEBRAL SCREW 7X35 MM	TI6AL4V-ELI	7x35mm
MCMVS-7x40	MC MULTIFONKSIYONEL VERTEBRAL SCREW 7X40 MM	TI6AL4V-ELI	7x40mm
MCMVS-7x45	MC MULTIFONKSIYONEL VERTEBRAL SCREW 7X45 MM	TI6AL4V-ELI	7x45mm
MCMVS-7x50	MC MULTIFONKSIYONEL VERTEBRAL SCREW 7X50 MM	TI6AL4V-ELI	7x50mm
MCMVS-7x55	MC MULTIFONKSIYONEL VERTEBRAL SCREW 7X55 MM	TI6AL4V-ELI	7x55mm
MCMVS-7x60	MC MULTIFONKSIYONEL VERTEBRAL SCREW 7X60 MM	TI6AL4V-ELI	7x60mm

# Thoracolumbar

## Pedicular Screws



MCMVS-7x65	MC MULTIFONKSIYONEL VERTEBRAL SCREW 7X65 MM	Ti6AL4V-ELI	7x65mm
MCMVS-7.5x30	MC MULTIFONKSIYONEL VERTEBRAL SCREW 7.5X30 MM	Ti6AL4V-ELI	7.5x30mm
MCMVS-7.5x35	MC MULTIFONKSIYONEL VERTEBRAL SCREW 7.5X35 MM	Ti6AL4V-ELI	7.5x35mm
MCMVS-7.5x40	MC MULTIFONKSIYONEL VERTEBRAL SCREW 7.5X40 MM	Ti6AL4V-ELI	7.5x40mm
MCMVS-7.5x45	MC MULTIFONKSIYONEL VERTEBRAL SCREW 7.5X45 MM	Ti6AL4V-ELI	7.5x45mm
MCMVS-7.5x50	MC MULTIFONKSIYONEL VERTEBRAL SCREW 7.5X50 MM	Ti6AL4V-ELI	7.5x50mm
MCMVS-7.5x55	MC MULTIFONKSIYONEL VERTEBRAL SCREW 7.5X55 MM	Ti6AL4V-ELI	7.5x55mm
MCMVS-7.5x60	MC MULTIFONKSIYONEL VERTEBRAL SCREW 7.5X60 MM	Ti6AL4V-ELI	7.5x60mm
MCMVS-7.5x65	MC MULTIFONKSIYONEL VERTEBRAL SCREW 7.5X65 MM	Ti6AL4V-ELI	7.5x65mm
MCMVS-8x30	MC MULTIFONKSIYONEL VERTEBRAL SCREW 8X30 MM	Ti6AL4V-ELI	8x30mm
MCMVS-8x35	MC MULTIFONKSIYONEL VERTEBRAL SCREW 8X35 MM	Ti6AL4V-ELI	8x35mm
MCMVS-8x40	MC MULTIFONKSIYONEL VERTEBRAL SCREW 8X40 MM	Ti6AL4V-ELI	8x40mm
MCMVS-8x45	MC MULTIFONKSIYONEL VERTEBRAL SCREW 8X45 MM	Ti6AL4V-ELI	8x45mm
MCMVS-8x50	MC MULTIFONKSIYONEL VERTEBRAL SCREW 8X50 MM	Ti6AL4V-ELI	8x50mm
MCMVS-8x55	MC MULTIFONKSIYONEL VERTEBRAL SCREW 8X55 MM	Ti6AL4V-ELI	8x55mm
MCMVS-8x60	MC MULTIFONKSIYONEL VERTEBRAL SCREW 8X60 MM	Ti6AL4V-ELI	8x60mm
MCMVS-8x65	MC MULTIFONKSIYONEL VERTEBRAL SCREW 8X65 MM	Ti6AL4V-ELI	8x65mm

# Thoracolumbar

## Pedicular Screws



MCRP-06	MC RIGID PLAQUE 6 HOLES	TI6AL4V-ELI
MCRP-07	MC RIGID PLAQUE 7 HOLES	TI6AL4V-ELI
MCRP-08	MC RIGID PLAQUE 8 HOLES	TI6AL4V-ELI

# Tlif - Plif

## Lumbar Peek Cage



CURETTE



CURETTE LEFT



CURETTE STRAIGHT

# Tlif - Plif

Lumbar Peek Cage



**MALLET**



**RASP**



**T HANDLE**



**RIMMER**




**SCRAPER**



**STRAIGHT HANDLE**

# Tlif - Plif

Lumbar Peek Cage



**TLIF INSERTER**



**TLIF TRIALS**



**SCRAPER 2**

# Tlif - Plif

Lumbar Peek Cage



**PLIF INSERTER**



**PLIF TRIALS**

# Lumbar Peek Cage



MCLEPPC-06X25MM	MC LOMBER EXPANDABLE PLIF PEEK CAGE 06X25MM	PEEK POLYMER
MCLEPPC-6.5X25MM	MC LOMBER EXPANDABLE PLIF PEEK CAGE 6.5X25MM	PEEK POLYMER
MCLEPPC-07X25MM	MC LOMBER EXPANDABLE PLIF PEEK CAGE 07X25MM	PEEK POLYMER
MCLEPPC-7.5X25MM	MC LOMBER EXPANDABLE PLIF PEEK CAGE 7.5X25MM	PEEK POLYMER
MCLEPPC-08X25MM	MC LOMBER EXPANDABLE PLIF PEEK CAGE 08X25MM	PEEK POLYMER
MCLEPPC-8.5X25MM	MC LOMBER EXPANDABLE PLIF PEEK CAGE 8.5X25MM	PEEK POLYMER
MCLEPPC-09X25MM	MC LOMBER EXPANDABLE PLIF PEEK CAGE 09X25MM	PEEK POLYMER
MCLEPPC-9.5X25MM	MC LOMBER EXPANDABLE PLIF PEEK CAGE 9.5X25MM	PEEK POLYMER
MCLEPPC-10X25MM	MC LOMBER EXPANDABLE PLIF PEEK CAGE 10X25MM	PEEK POLYMER
MCLEPPC-10.5X25MM	MC LOMBER EXPANDABLE PLIF PEEK CAGE 10.5X25MM	PEEK POLYMER
MCLEPPC-11X25MM	MC LOMBER EXPANDABLE PLIF PEEK CAGE 11X25MM	PEEK POLYMER
MCLEPPC-11.5X25MM	MC LOMBER EXPANDABLE PLIF PEEK CAGE 11.5X25MM	PEEK POLYMER
MCLEPPC-12X25MM	MC LOMBER EXPANDABLE PLIF PEEK CAGE 12X25MM	PEEK POLYMER
MCLEPPC-12.5X25MM	MC LOMBER EXPANDABLE PLIF PEEK CAGE 12.5X25MM	PEEK POLYMER
MCLEPPC-13X25MM	MC LOMBER EXPANDABLE PLIF PEEK CAGE 13X25MM	PEEK POLYMER
MCLEPPC-14X25MM	MC LOMBER EXPANDABLE PLIF PEEK CAGE 14X25MM	PEEK POLYMER

# Lumbar Peek Cage




MCLPPC-06X25MM	MC LOMBER PLIF PEEK CAGE 06X25MM	PEEK POLYMER
MCLPPC-6.5X25MM	MC LOMBER PLIF PEEK CAGE 6.5X25MM	PEEK POLYMER
MCLPPC-07X25MM	MC LOMBER PLIF PEEK CAGE 07X25MM	PEEK POLYMER
MCLPPC-7.5X25MM	MC LOMBER PLIF PEEK CAGE 7.5X25MM	PEEK POLYMER
MCLPPC-08X25MM	MC LOMBER PLIF PEEK CAGE 08X25MM	PEEK POLYMER
MCLPPC-8.5X25MM	MC LOMBER PLIF PEEK CAGE 8.5X25MM	PEEK POLYMER
MCLPPC-09X25MM	MC LOMBER PLIF PEEK CAGE 09X25MM	PEEK POLYMER
MCLPPC-9.5X25MM	MC LOMBER PLIF PEEK CAGE 9.5X25MM	PEEK POLYMER
MCLPPC-10X25MM	MC LOMBER PLIF PEEK CAGE 10X25MM	PEEK POLYMER
MCLPPC-10.5X25MM	MC LOMBER PLIF PEEK CAGE 10.5X25MM	PEEK POLYMER
MCLPPC-11X25MM	MC LOMBER PLIF PEEK CAGE 11X25MM	PEEK POLYMER
MCLPPC-11.5X25MM	MC LOMBER PLIF PEEK CAGE 11.5X25MM	PEEK POLYMER
MCLPPC-12X25MM	MC LOMBER PLIF PEEK CAGE 12X25MM	PEEK POLYMER
MCLPPC-12.5X25MM	MC LOMBER PLIF PEEK CAGE 12.5X25MM	PEEK POLYMER
MCLPPC-13X25MM	MC LOMBER PLIF PEEK CAGE 13X25MM	PEEK POLYMER
MCLPPC-14X25MM	MC LOMBER PLIF PEEK CAGE 14X25MM	PEEK POLYMER

# Lumbar Peek Cage



MCTPC-7X32 MM	MC TLIF PEEK CAGE 7X32 MM	PEEK POLYMER
MCTPC-7.5X32 MM	MC TLIF PEEK CAGE 7.5X32 MM	PEEK POLYMER
MCTPC-8X32 MM	MC TLIF PEEK CAGE 8X32 MM	PEEK POLYMER
MCTPC-8.5X32 MM	MC TLIF PEEK CAGE 8.5X32 MM	PEEK POLYMER
MCTPC-9X32 MM	MC TLIF PEEK CAGE 9X32 MM	PEEK POLYMER
MCTPC-9.5X32 MM	MC TLIF PEEK CAGE 9.5X32 MM	PEEK POLYMER
MCTPC-10X32 MM	MC TLIF PEEK CAGE 10X32 MM	PEEK POLYMER
MCTPC-10.5X32 MM	MC TLIF PEEK CAGE 10.5X32 MM	PEEK POLYMER
MCTPC-11X32 MM	MC TLIF PEEK CAGE 11X32 MM	PEEK POLYMER
MCTPC-11.5X32 MM	MC TLIF PEEK CAGE 11.5X32 MM	PEEK POLYMER
MCTPC-12X32 MM	MC TLIF PEEK CAGE 12X32 MM	PEEK POLYMER
MCTPC-12.5X32 MM	MC TLIF PEEK CAGE 12.5X32 MM	PEEK POLYMER
MCTPC-13X32 MM	MC TLIF PEEK CAGE 13X32 MM	PEEK POLYMER
MCTPC-14X32 MM	MC TLIF PEEK CAGE 14X32 MM	PEEK POLYMER
MCTPC-15X32 MM	MC TLIF PEEK CAGE 15X32 MM	PEEK POLYMER

# Posterior Cervical Pedicular Screws



3NM TORK HANDLE



OCCIPITAL PLATE  
SCREW DRIVER



ANTI TORQUE 2

# Posterior Cervical Pedicular Screws



**BONE AWL**



**PROBE CURVED**



**TAB 3MM**



**PROBE FEELER CURVED**



**ANTI TORQUE**



**PROBE FEELER STRAIGHT**

# Posterior Cervical Pedicular Screws



**COMPRESSOR**



**DISTRACTOR**



**DRILL GUIDE**

# Posterior Cervical Pedicular Screws



**INSUTU BENDER LEFT-RIGHT**



**TRANSVER DRIVER**



**T HANDLE**



**STRAIGHT HANDLE**



**PROBE STRAIGHT**



**NUT INSERTER**

# Posterior Cervical Pedicular Screws



**POLYAXIAL SCREW DRIVER**



**ROCKER**



**MALLET**



**ROD BENDER**



**ROD PUSHER**

# Posterior Cervical Pedicular Screws



MCTTB-50	MC CROSSLINK 50 MM	TI6AL4V-ELI	50mm
MCTTB-60	MC CROSSLINK 60 MM	TI6AL4V-ELI	60mm

# Posterior Cervical Pedicular Screws



MCFFCS-4.5X25MM	MC FACET FIXASYON CANNULATED SCREW 4.5X25	T16AL4V-ELI	4.5x25mm
MCFFCS-4.5X30MM	MC FACET FIXASYON CANNULATED SCREW 4.5X30	T16AL4V-ELI	4.5x30mm
MCFFCS-4.5X35MM	MC FACET FIXASYON CANNULATED SCREW 4.5X35	T16AL4V-ELI	4.5x35mm
MCFFCS-4.5X40MM	MC FACET FIXASYON CANNULATED SCREW 4.5X40	T16AL4V-ELI	4.5x40mm
MCFFCS-4.5X45MM	MC FACET FIXASYON CANNULATED SCREW 4.5X45	T16AL4V-ELI	4.5x45mm
MCFFCS-4.5X50MM	MC FACET FIXASYON CANNULATED SCREW 4.5X50	T16AL4V-ELI	4.5x50mm
MCFFCS-5X25MM	MC FACET FIXASYON CANNULATED SCREW 5X25	T16AL4V-ELI	5x25mm
MCFFCS-5X30MM	MC FACET FIXASYON CANNULATED SCREW 5X30	T16AL4V-ELI	5x30mm
MCFFCS-5X35MM	MC FACET FIXASYON CANNULATED SCREW 5X35	T16AL4V-ELI	5x35mm
MCFFCS-5X40MM	MC FACET FIXASYON CANNULATED SCREW 5X40	T16AL4V-ELI	5x40mm
MCFFCS-5X45MM	MC FACET FIXASYON CANNULATED SCREW 5X45	T16AL4V-ELI	5x45mm
MCFFCS-5X50MM	MC FACET FIXASYON CANNULATED SCREW 5X50	T16AL4V-ELI	4x50mm
MCFFCS-5.5X25MM	MC FACET FIXASYON CANNULATED SCREW 5.5X25	T16AL4V-ELI	5.5x25mm
MCFFCS-5.5X30MM	MC FACET FIXASYON CANNULATED SCREW 5.5X30	T16AL4V-ELI	5.5x30mm
MCFFCS-5.5X35MM	MC FACET FIXASYON CANNULATED SCREW 5.5X35	T16AL4V-ELI	5.5x35mm
MCFFCS-5.5X40MM	MC FACET FIXASYON CANNULATED SCREW 5.5X40	T16AL4V-ELI	5.5x40mm
MCFFCS-5.5X45MM	MC FACET FIXASYON CANNULATED SCREW 5.5X45	T16AL4V-ELI	5.5x45mm
MCFFCS-5.5X50MM	MC FACET FIXASYON CANNULATED SCREW 5.5X50	T16AL4V-ELI	5.5x50mm
MCFFS-3.5X25MM	MC FACET FIXASYON SCREW 3.5X25	T16AL4V-ELI	3.5x25mm
MCFFS-3.5X30MM	MC FACET FIXASYON SCREW 3.5X30	T16AL4V-ELI	3.5x30mm

# Posterior Cervical Pedicular Screws



MCFFS-3.5X35MM	MC FACET FIXASYON SCREW 3.5X35	Ti6AL4V-ELI	3.5x35mm
MCFFS-4.5X25MM	MC FACET FIXASYON SCREW 4.5X25	Ti6AL4V-ELI	4.5x25mm
MCFFS-4.5X30MM	MC FACET FIXASYON SCREW 4.5X30	Ti6AL4V-ELI	4.5x30mm
MCFFS-4.5X35MM	MC FACET FIXASYON SCREW 4.5X35	Ti6AL4V-ELI	4.5x35mm
MCFFS-4.5X40MM	MC FACET FIXASYON SCREW 4.5X40	Ti6AL4V-ELI	4.5x40mm
MCFFS-4.5X45MM	MC FACET FIXASYON SCREW 4.5X45	Ti6AL4V-ELI	4.5x45mm
MCFFS-4.5X50MM	MC FACET FIXASYON SCREW 4.5X50	Ti6AL4V-ELI	4.5x50mm
MCFFS-5X25MM	MC FACET FIXASYON SCREW 5X25	Ti6AL4V-ELI	5x25mm
MCFFS-5X30MM	MC FACET FIXASYON SCREW 5X30	Ti6AL4V-ELI	5x30mm
MCFFS-5X35MM	MC FACET FIXASYON SCREW 5X35	Ti6AL4V-ELI	5x35mm
MCFFS-5X40MM	MC FACET FIXASYON SCREW 5X40	Ti6AL4V-ELI	5x40mm
MCFFS-5X45MM	MC FACET FIXASYON SCREW 5X45	Ti6AL4V-ELI	5x45mm
MCFFS-5X50MM	MC FACET FIXASYON SCREW 5X50	Ti6AL4V-ELI	5x50mm
MCFFS-5.5X25MM	MC FACET FIXASYON SCREW 5.5X25	Ti6AL4V-ELI	5.5x25mm
MCFFS-5.5X30MM	MC FACET FIXASYON SCREW 5.5X30	Ti6AL4V-ELI	5.5x30mm
MCFFS-5.5X35MM	MC FACET FIXASYON SCREW 5.5X35	Ti6AL4V-ELI	5.5x35mm
MCFFS-5.5X40MM	MC FACET FIXASYON SCREW 5.5X40	Ti6AL4V-ELI	5.5x40mm
MCFFS-5.5X45MM	MC FACET FIXASYON SCREW 5.5X45	Ti6AL4V-ELI	5.5x45mm
MCFFS-5.5X50MM	MC FACET FIXASYON SCREW 5.5X50	Ti6AL4V-ELI	5.5x50mm



# Posterior Cervical Pedicular Screws



# Posterior Cervical Pedicular Screws

MCLHL-5X5	MC LAMINAR HOOK OFFSET (L) 5X5 MM	5x5mm
MCLHL-5X7	MC LAMINAR HOOK OFFSET (L) 5X7 MM	5x7mm
MCLHL-5X9	MC LAMINAR HOOK OFFSET (L) 5X9 MM	5x9mm
MCLHL-5X11	MC LAMINAR HOOK OFFSET (L) 5X11 MM	5x11mm
MCLHR-5X5	MC LAMINAR HOOK OFFSET (R) 5X5 MM	5x5mm
MCLHR-5X7	MC LAMINAR HOOK OFFSET (R) 5X7 MM	5x7mm
MCLHR-5X9	MC LAMINAR HOOK OFFSET (R) 5X9 MM	5x9mm
MCLHR-5X11	MC LAMINAR HOOK OFFSET (R) 5X11 MM	5x11mm
MCLH-5X5	MC LAMINAR HOOK 5X5 MM	5x5mm
MCLH-5X7	MC LAMINAR HOOK 5X7 MM	5x7mm
MCLH-5X9	MC LAMINAR HOOK 5X9 MM	5x9mm
MCPH-5X5	MC PEDICLE HOOK 5X5 MM	5x5mm
MCPH-5X7	MC PEDICLE HOOK 5X7 MM	5x7mm
MCPH-5X9	MC PEDICLE HOOK 5X9 MM	5x9mm
MCPLH-5X5	MC POLYAXIAL LAMINAR HOOK 5X5MM	5x5mm
MCPLH-5X6	MC POLYAXIAL LAMINAR HOOK 5X6MM	5x6mm
MCPLH-5X7	MC POLYAXIAL LAMINAR HOOK 5X7MM	5x7mm
MCPPH-5X5	MC POLYAXIAL PEDICLE HOOK 5X5MM	5x5mm
MCPPH-5X6	MC POLYAXIAL PEDICLE HOOK 5X6MM	5x6mm
MCPPH-5X7	MC POLYAXIAL PEDICLE HOOK 5X7MM	5x7mm
MCTH-5X5MM	MC TRANSVERS HOOK 5X5MM	5x5mm
MCTH-5X6MM	MC TRANSVERS HOOK 5X6MM	5x6mm
MCTH-5X7MM	MC TRANSVERS HOOK 5X7MM	5x7mm
MCTH-5X9MM	MC TRANSVERS HOOK 5X9MM	5x9mm
MCTH-5X11MM	MC TRANSVERS HOOK 5X11MM	5x11mm
MCTH-5X13MM	MC TRANSVERS HOOK 5X13MM	5x13mm
MCSH-5X5MM	MC SACRAL HOOK 5X5MM	5x5mm
MCSH-5X6MM	MC SACRAL HOOK 5X6MM	5x6mm
MCSH-5X7MM	MC SACRAL HOOK 5X7MM	5x7mm
MCSH-5X9MM	MC SACRAL HOOK 5X9MM	5x9mm
MCSH-5X11MM	MC SACRAL HOOK 5X11MM	5x11mm
MCSH-5X13MM	MC SACRAL HOOK 5X13MM	5x13mm
MCOC	MC OFFSETT CONNECTOR	

# Posterior Cervical Pedicular Screws



---

MCCLC    MC Closed Lateral Connector

---



---

MCRP-03    MC RIGID PLAQUE 3 HOLES

---

---

MCRP-04    MC RIGID PLAQUE 4 HOLES

---

---

MCRP-05    MC RIGID PLAQUE 5 HOLES

---

# Posterior Cervical Pedicular Screws




MCPFR-3.5X40	MC POSTERIOR FIXATION ROD 3.5X40	Ti6AL4V-ELI	3.5x40mm
MCPFR-3.5X50	MC POSTERIOR FIXATION ROD 3.5X50	Ti6AL4V-ELI	3.5x50mm
MCPFR-3.5X60	MC POSTERIOR FIXATION ROD 3.5X60	Ti6AL4V-ELI	3.5x60mm
MCPFR-3.5X70	MC POSTERIOR FIXATION ROD 3.5X70	Ti6AL4V-ELI	3.5x70mm
MCPFR-3.5X80	MC POSTERIOR FIXATION ROD 3.5X80	Ti6AL4V-ELI	3.5x80mm
MCPFR-3.5X90	MC POSTERIOR FIXATION ROD 3.5X90	Ti6AL4V-ELI	3.5x90mm
MCPFR-3.5X100	MC POSTERIOR FIXATION ROD 3.5X100	Ti6AL4V-ELI	3.5x100mm
MCPFR-4X40	MC POSTERIOR FIXATION ROD 4X40	Ti6AL4V-ELI	4x40mm
MCPFR-4X50	MC POSTERIOR FIXATION ROD 4X50	Ti6AL4V-ELI	4x50mm
MCPFR-4X60	MC POSTERIOR FIXATION ROD 4X60	Ti6AL4V-ELI	4x60mm
MCPFR-4X70	MC POSTERIOR FIXATION ROD 4X70	Ti6AL4V-ELI	4x70mm
MCPFR-4X80	MC POSTERIOR FIXATION ROD 4X80	Ti6AL4V-ELI	4x80mm
MCPFR-4X90	MC POSTERIOR FIXATION ROD 4X90	Ti6AL4V-ELI	4x90mm
MCPFR-4X100	MC POSTERIOR FIXATION ROD 4X100	Ti6AL4V-ELI	4x100mm


# Posterior Cervical Pedicular Screws



MCPS-3X10	MC POLYAXIAL SCREW 3X10MM	TI6AL4V-ELI	3x10mm
MCPS-3X12	MC POLYAXIAL SCREW 3X12MM	TI6AL4V-ELI	3x12mm
MCPS-3X14	MC POLYAXIAL SCREW 3X14MM	TI6AL4V-ELI	3x14mm
MCPS-3X16	MC POLYAXIAL SCREW 3X16MM	TI6AL4V-ELI	3x16mm
MCPS-3X18	MC POLYAXIAL SCREW 3X18MM	TI6AL4V-ELI	3x18mm
MCPS-3X20	MC POLYAXIAL SCREW 3X20MM	TI6AL4V-ELI	3x20mm
MCPS-3X22	MC POLYAXIAL SCREW 3X22MM	TI6AL4V-ELI	3x22mm
MCPS-3X24	MC POLYAXIAL SCREW 3X24MM	TI6AL4V-ELI	3x24mm
MCPS-3X26	MC POLYAXIAL SCREW 3X26MM	TI6AL4V-ELI	3x26mm
MCPS-3X28	MC POLYAXIAL SCREW 3X28MM	TI6AL4V-ELI	3x28mm
MCPS-3X30	MC POLYAXIAL SCREW 3X30MM	TI6AL4V-ELI	3x30mm
MCPS-3.5X10	MC POLYAXIAL SCREW 3.5X10MM	TI6AL4V-ELI	3.5x10mm
MCPS-3.5X12	MC POLYAXIAL SCREW 3.5X12MM	TI6AL4V-ELI	3.5x12mm
MCPS-3.5X14	MC POLYAXIAL SCREW 3.5X14MM	TI6AL4V-ELI	3.5x14mm
MCPS-3.5X16	MC POLYAXIAL SCREW 3.5X16MM	TI6AL4V-ELI	3.5x16mm
MCPS-3.5X18	MC POLYAXIAL SCREW 3.5X18MM	TI6AL4V-ELI	3.5x18mm
MCPS-3.5X20	MC POLYAXIAL SCREW 3.5X20MM	TI6AL4V-ELI	3.5x20mm
MCPS-3.5X22	MC POLYAXIAL SCREW 3.5X22MM	TI6AL4V-ELI	3.5x22mm
MCPS-3.5X24	MC POLYAXIAL SCREW 3.5X24MM	TI6AL4V-ELI	3.5x24mm
MCPS-3.5X26	MC POLYAXIAL SCREW 3.5X26MM	TI6AL4V-ELI	3.5x26mm



# Posterior Cervical Pedicular Screws




MCPS-3.5X28	MC POLYAXIAL SCREW 3.5X28MM	Ti6AL4V-ELI	3.5x28mm
MCPS-3.5X30	MC POLYAXIAL SCREW 3.5X30MM	Ti6AL4V-ELI	3.5x30mm
MCPS-4X10	MC POLYAXIAL SCREW 4X10MM	Ti6AL4V-ELI	4x10mm
MCPS-4X12	MC POLYAXIAL SCREW 4X12MM	Ti6AL4V-ELI	4x12mm
MCPS-4X14	MC POLYAXIAL SCREW 4X14MM	Ti6AL4V-ELI	4x14mm
MCPS-4X16	MC POLYAXIAL SCREW 4X16MM	Ti6AL4V-ELI	4x16mm
MCPS-4X18	MC POLYAXIAL SCREW 4X18MM	Ti6AL4V-ELI	4x18mm
MCPS-4X20	MC POLYAXIAL SCREW 4X20MM	Ti6AL4V-ELI	4x20mm
MCPS-4X22	MC POLYAXIAL SCREW 4X22MM	Ti6AL4V-ELI	4x22mm
MCPS-4X24	MC POLYAXIAL SCREW 4X24MM	Ti6AL4V-ELI	4x24mm
MCPS-4X26	MC POLYAXIAL SCREW 4X26MM	Ti6AL4V-ELI	4x26mm
MCPS-4X28	MC POLYAXIAL SCREW 4X28MM	Ti6AL4V-ELI	4x28mm
MCPS-4X30	MC POLYAXIAL SCREW 4X30MM	Ti6AL4V-ELI	4x30mm
MCPS-4.5X25	MC POLYAXIAL SCREW 4.5X25MM	Ti6AL4V-ELI	4.5x25mm
MCPS-4.5X30	MC POLYAXIAL SCREW 4.5X30MM	Ti6AL4V-ELI	4.5x30mm
MCPS-4.5X35	MC POLYAXIAL SCREW 4.5X35MM	Ti6AL4V-ELI	4.5x35mm
MCPS-4.5X40	MC POLYAXIAL SCREW 4.5X40MM	Ti6AL4V-ELI	4.5x40mm
MCPS-4.5X45	MC POLYAXIAL SCREW 4.5X45MM	Ti6AL4V-ELI	4.5x45mm
MCPS-4.5X50	MC POLYAXIAL SCREW 4.5X50MM	Ti6AL4V-ELI	4.5x50mm



# Anterior

Cervical Prosthesis-Cervical Peek Cage Combine Set



CASPAR EKARTOR



PIN DRIVER



BONE AWL

# Anterior

## Cervical Prosthesis-Cervical Peek Cage Combine Set



**CERVICAL CAGE HOLDER**



**CERVICAL BLADE CAGE HOLDER**



**EXPANDABLE CERVICAL  
BLADE CAGE HOLDER**



**CERVICAL PROTHESIS HOLDER**



**MALLET**



**PIN**

# Anterior

## Cervical Prosthesis-Cervical Peek Cage Combine Set



**CERVICAL PROTHESIS HOLDER**



**CERVICAL CAGE TRIALS**

# Anterior

## Cervical Prosthesis-Cervical Peek Cage Combine Set



MCKPC-12X04MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 12X04 MM	PEEK POLYMER
MCKPC-12X4.5MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 12X4.5 MM	PEEK POLYMER
MCKPC-12X05 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 12X05 MM	PEEK POLYMER
MCKPC-12X5.5MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 12X5.5 MM	PEEK POLYMER
MCKPC-12X06 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 12X06 MM	PEEK POLYMER
MCKPC-12X6.5 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 12X6.5 MM	PEEK POLYMER
MCKPC-12X07 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 12X07 MM	PEEK POLYMER
MCKPC-12X7.5 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 12X7.5MM	PEEK POLYMER
MCKPC-12X08 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 12X08 MM	PEEK POLYMER
MCKPC-12X8.5 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 12X8.5 MM	PEEK POLYMER
MCKPC-13X04MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 13X04 MM	PEEK POLYMER
MCKPC-13X4.5MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 13X4.5 MM	PEEK POLYMER
MCKPC-13X05 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 13X05 MM	PEEK POLYMER
MCKPC-13X5.5MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 13X5.5 MM	PEEK POLYMER
MCKPC-13X06 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 13X06 MM	PEEK POLYMER
MCKPC-13X6.5 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 13X6.5 MM	PEEK POLYMER
MCKPC-13X07 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 13X07 MM	PEEK POLYMER
MCKPC-13X7.5 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 13X7.5MM	PEEK POLYMER
MCKPC-13X08 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 13X08 MM	PEEK POLYMER
MCKPC-13X8.5 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 13X8.5 MM	PEEK POLYMER
MCKPC-14X04MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 14X04 MM	PEEK POLYMER
MCKPC-14X4.5MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 14X4.5 MM	PEEK POLYMER
MCKPC-14X05 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 14X05 MM	PEEK POLYMER
MCKPC-14X5.5MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 14X5.5 MM	PEEK POLYMER
MCKPC-14X06 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 14X06 MM	PEEK POLYMER

# Anterior

## Cervical Prosthesis-Cervical Peek Cage Combine Set



MCCCKPC-14X6.5 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 14X6.5 MM	PEEK POLYMER
MCCCKPC-14X07 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 14X07 MM	PEEK POLYMER
MCCCKPC-14X7.5 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 14X7.5MM	PEEK POLYMER
MCCCKPC-14X08 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 14X08 MM	PEEK POLYMER
MCCCKPC-14X8.5 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 14X8.5 MM	PEEK POLYMER
MCCCKPC-15X04MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 15X04 MM	PEEK POLYMER
MCCCKPC-15X4.5MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 15X4.5 MM	PEEK POLYMER
MCCCKPC-15X05 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 15X05 MM	PEEK POLYMER
MCCCKPC-15X5.5MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 15X5.5 MM	PEEK POLYMER
MCCCKPC-15X06 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 15X06 MM	PEEK POLYMER
MCCCKPC-15X6.5 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 15X6.5 MM	PEEK POLYMER
MCCCKPC-15X07 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 15X07 MM	PEEK POLYMER
MCCCKPC-15X7.5 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 15X7.5MM	PEEK POLYMER
MCCCKPC-15X08 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 15X08 MM	PEEK POLYMER
MCCCKPC-15X8.5 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 15X8.5 MM	PEEK POLYMER
MCCCKPC-16X04MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 16X04 MM	PEEK POLYMER
MCCCKPC-16X4.5MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 16X4.5 MM	PEEK POLYMER
MCCCKPC-16X05 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 16X05 MM	PEEK POLYMER
MCCCKPC-16X5.5MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 16X5.5 MM	PEEK POLYMER
MCCCKPC-16X06 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 16X06 MM	PEEK POLYMER
MCCCKPC-16X6.5 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 16X6.5 MM	PEEK POLYMER
MCCCKPC-16X07 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 16X07 MM	PEEK POLYMER
MCCCKPC-16X7.5 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 16X7.5MM	PEEK POLYMER
MCCCKPC-16X08 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 16X08 MM	PEEK POLYMER
MCCCKPC-16X8.5 MM	MC ANTERIOR CERVICAL KNIFE PEEK CAGE 16X8.5 MM	PEEK POLYMER

# Anterior

## Cervical Prosthesis-Cervical Peek Cage Combine Set



MCCPC-12X04MM	MC ANTERIOR CERVICAL PEEK CAGE 12X04 MM	PEEK POLYMER
MCCPC-12X4.5MM	MC ANTERIOR CERVICAL PEEK CAGE 12X4.5 MM	PEEK POLYMER
MCCPC-12X05 MM	MC ANTERIOR CERVICAL PEEK CAGE 12X05 MM	PEEK POLYMER
MCCPC-12X5.5MM	MC ANTERIOR CERVICAL PEEK CAGE 12X5.5 MM	PEEK POLYMER
MCCPC-12X06 MM	MC ANTERIOR CERVICAL PEEK CAGE 12X06 MM	PEEK POLYMER
MCCPC-12X6.5 MM	MC ANTERIOR CERVICAL PEEK CAGE 12X6.5 MM	PEEK POLYMER
MCCPC-12X07 MM	MC ANTERIOR CERVICAL PEEK CAGE 12X07 MM	PEEK POLYMER
MCCPC-12X7.5 MM	MC ANTERIOR CERVICAL PEEK CAGE 12X7.5MM	PEEK POLYMER
MCCPC-12X08 MM	MC ANTERIOR CERVICAL PEEK CAGE 12X08 MM	PEEK POLYMER
MCCPC-12X8.5 MM	MC ANTERIOR CERVICAL PEEK CAGE 12X8.5 MM	PEEK POLYMER
MCCPC-13X04MM	MC ANTERIOR CERVICAL PEEK CAGE 13X04 MM	PEEK POLYMER
MCCPC-13X4.5MM	MC ANTERIOR CERVICAL PEEK CAGE 13X4.5 MM	PEEK POLYMER
MCCPC-13X05 MM	MC ANTERIOR CERVICAL PEEK CAGE 13X05 MM	PEEK POLYMER
MCCPC-13X5.5MM	MC ANTERIOR CERVICAL PEEK CAGE 13X5.5 MM	PEEK POLYMER
MCCPC-13X06 MM	MC ANTERIOR CERVICAL PEEK CAGE 13X06 MM	PEEK POLYMER
MCCPC-13X6.5 MM	MC ANTERIOR CERVICAL PEEK CAGE 13X6.5 MM	PEEK POLYMER
MCCPC-13X07 MM	MC ANTERIOR CERVICAL PEEK CAGE 13X07 MM	PEEK POLYMER
MCCPC-13X7.5 MM	MC ANTERIOR CERVICAL PEEK CAGE 13X7.5MM	PEEK POLYMER
MCCPC-13X08 MM	MC ANTERIOR CERVICAL PEEK CAGE 13X08 MM	PEEK POLYMER
MCCPC-13X8.5 MM	MC ANTERIOR CERVICAL PEEK CAGE 13X8.5 MM	PEEK POLYMER
MCCPC-14X04MM	MC ANTERIOR CERVICAL PEEK CAGE 14X04 MM	PEEK POLYMER
MCCPC-14X4.5MM	MC ANTERIOR CERVICAL PEEK CAGE 14X4.5 MM	PEEK POLYMER
MCCPC-14X05 MM	MC ANTERIOR CERVICAL PEEK CAGE 14X05 MM	PEEK POLYMER
MCCPC-14X5.5MM	MC ANTERIOR CERVICAL PEEK CAGE 14X5.5 MM	PEEK POLYMER
MCCPC-14X06 MM	MC ANTERIOR CERVICAL PEEK CAGE 14X06 MM	PEEK POLYMER

# Anterior

## Cervical Prosthesis-Cervical Peek Cage Combine Set



MCCPC-14X6.5 MM	MC ANTERIOR CERVICAL PEEK CAGE 14X6.5 MM	PEEK POLYMER
MCCPC-14X07 MM	MC ANTERIOR CERVICAL PEEK CAGE 14X07 MM	PEEK POLYMER
MCCPC-14X7.5 MM	MC ANTERIOR CERVICAL PEEK CAGE 14X7.5MM	PEEK POLYMER
MCCPC-14X08 MM	MC ANTERIOR CERVICAL PEEK CAGE 14X08 MM	PEEK POLYMER
MCCPC-14X8.5 MM	MC ANTERIOR CERVICAL PEEK CAGE 14X8.5 MM	PEEK POLYMER
MCCPC-15X04MM	MC ANTERIOR CERVICAL PEEK CAGE 15X04 MM	PEEK POLYMER
MCCPC-15X4.5MM	MC ANTERIOR CERVICAL PEEK CAGE 15X4.5 MM	PEEK POLYMER
MCCPC-15X05 MM	MC ANTERIOR CERVICAL PEEK CAGE 15X05 MM	PEEK POLYMER
MCCPC-15X5.5MM	MC ANTERIOR CERVICAL PEEK CAGE 15X5.5 MM	PEEK POLYMER
MCCPC-15X06 MM	MC ANTERIOR CERVICAL PEEK CAGE 15X06 MM	PEEK POLYMER
MCCPC-15X6.5 MM	MC ANTERIOR CERVICAL PEEK CAGE 15X6.5 MM	PEEK POLYMER
MCCPC-15X07 MM	MC ANTERIOR CERVICAL PEEK CAGE 15X07 MM	PEEK POLYMER
MCCPC-15X7.5 MM	MC ANTERIOR CERVICAL PEEK CAGE 15X7.5MM	PEEK POLYMER
MCCPC-15X08 MM	MC ANTERIOR CERVICAL PEEK CAGE 15X08 MM	PEEK POLYMER
MCCPC-15X8.5 MM	MC ANTERIOR CERVICAL PEEK CAGE 15X8.5 MM	PEEK POLYMER
MCCPC-16X04MM	MC ANTERIOR CERVICAL PEEK CAGE 16X04 MM	PEEK POLYMER
MCCPC-16X4.5MM	MC ANTERIOR CERVICAL PEEK CAGE 16X4.5 MM	PEEK POLYMER
MCCPC-16X05 MM	MC ANTERIOR CERVICAL PEEK CAGE 16X05 MM	PEEK POLYMER
MCCPC-16X5.5MM	MC ANTERIOR CERVICAL PEEK CAGE 16X5.5 MM	PEEK POLYMER
MCCPC-16X06 MM	MC ANTERIOR CERVICAL PEEK CAGE 16X06 MM	PEEK POLYMER
MCCPC-16X6.5 MM	MC ANTERIOR CERVICAL PEEK CAGE 16X6.5 MM	PEEK POLYMER
MCCPC-16X07 MM	MC ANTERIOR CERVICAL PEEK CAGE 16X07 MM	PEEK POLYMER
MCCPC-16X7.5 MM	MC ANTERIOR CERVICAL PEEK CAGE 16X7.5MM	PEEK POLYMER
MCCPC-16X08 MM	MC ANTERIOR CERVICAL PEEK CAGE 16X08 MM	PEEK POLYMER
MCCPC-16X8.5 MM	MC ANTERIOR CERVICAL PEEK CAGE 16X8.5 MM	PEEK POLYMER

# Anterior

## Cervical Prosthesis-Cervical Peek Cage Combine Set



MCECKPC-12X04 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 12X04MM	PEEK POLYMER
MCECKPC-12X4.5 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 12X4.5MM	PEEK POLYMER
MCECKPC-12X05 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 12X05MM	PEEK POLYMER
MCECKPC-12X5.5 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 12X5.5MM	PEEK POLYMER
MCECKPC-12X06 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 12X06 MM	PEEK POLYMER
MCECKPC-12 6.5 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 12X6.5 MM	PEEK POLYMER
MCECKPC-12X07 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 12X07 MM	PEEK POLYMER
MCECKPC-12X7.5 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 12X7.5 MM	PEEK POLYMER
MCECKPC-12X08 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 12X08 MM	PEEK POLYMER
MCECKPC-12X8.5 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 12X8.5 MM	PEEK POLYMER
MCECKPC-13X04 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 13X04MM	PEEK POLYMER
MCECKPC-13X4.5 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 13X4.5MM	PEEK POLYMER
MCECKPC-13X05 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 13X05MM	PEEK POLYMER
MCECKPC-13X5.5 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 13X5.5MM	PEEK POLYMER
MCECKPC-13X06 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 13X06 MM	PEEK POLYMER
MCECKPC-13X6.5 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 13X6.5 MM	PEEK POLYMER
MCECKPC-13X07 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 13X07 MM	PEEK POLYMER
MCECKPC-13X7.5 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 13X7.5 MM	PEEK POLYMER
MCECKPC-13X08 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 13X08 MM	PEEK POLYMER
MCECKPC-13X8.5 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 13X8.5 MM	PEEK POLYMER
MCECKPC-14X04 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 14X04MM	PEEK POLYMER
MCECKPC-14X4.5 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 14X4.5MM	PEEK POLYMER
MCECKPC-14X05 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 14X05MM	PEEK POLYMER
MCECKPC-14X5.5 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 14X5.5MM	PEEK POLYMER
MCECKPC-14X06 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 14X06 MM	PEEK POLYMER
MCECKPC-14X6.5 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 14X6.5 MM	PEEK POLYMER

# Anterior

## Cervical Prosthesis-Cervical Peek Cage Combine Set



MCECKPC-14X07 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 14X07 MM	PEEK POLYMER
MCECKPC-14X7.5 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 14X7.5 MM	PEEK POLYMER
MCECKPC-14X08 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 14X08 MM	PEEK POLYMER
MCECKPC-14X8.5 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 14X8.5 MM	PEEK POLYMER
MCECKPC-15X04 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 15X04MM	PEEK POLYMER
MCECKPC-15X4.5 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 15X4.5MM	PEEK POLYMER
MCECKPC-15x05 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 15X05MM	PEEK POLYMER
MCECKPC-15X5.5 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 15X5.5MM	PEEK POLYMER
MCECKPC-15X06 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 15X06 MM	PEEK POLYMER
MCECKPC-15X6.5 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 15X6.5 MM	PEEK POLYMER
MCECKPC-15X07 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 15X07 MM	PEEK POLYMER
MCECKPC-15X7.5 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 15X7.5 MM	PEEK POLYMER
MCECKPC-15X08 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 15X08 MM	PEEK POLYMER
MCECKPC-15X8.5 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 15X8.5 MM	PEEK POLYMER
MCECKPC-16X04 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 16X04MM	PEEK POLYMER
MCECKPC-16X4.5 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 16X4.5MM	PEEK POLYMER
MCECKPC-16x05 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 16X05MM	PEEK POLYMER
MCECKPC-16X5.5 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 16X5.5MM	PEEK POLYMER
MCECKPC-16X06 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 16X06 MM	PEEK POLYMER
MCECKPC-16X6.5 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 16X6.5 MM	PEEK POLYMER
MCECKPC-16X07 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 16X07 MM	PEEK POLYMER
MCECKPC-16X7.5 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 16X7.5 MM	PEEK POLYMER
MCECKPC-16X08 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 16X08 MM	PEEK POLYMER
MCECKPC-16X8.5 MM	MC EXPANDABLE CERVICAL KNIFE PEEK CAGE 16X8.5 MM	PEEK POLYMER

# Anterior

## Cervical Prosthesis-Cervical Peek Cage Combine Set



MCCDP-12X04MM	MC CERVICAL DISC PROTHESIS 12X04MM	TI6AL4V-ELI	12x4mm
MCCDP-12X4.5MM	MC CERVICAL DISC PROTHESIS 12X4.5MM	TI6AL4V-ELI	12x4.5mm
MCCDP-12X05MM	MC CERVICAL DISC PROTHESIS 12X05MM	TI6AL4V-ELI	12x5mm
MCCDP-12X5.5MM	MC CERVICAL DISC PROTHESIS 12X5.5MM	TI6AL4V-ELI	12x5.5mm
MCCDP-12X06MM	MC CERVICAL DISC PROTHESIS 12X06MM	TI6AL4V-ELI	12x6mm
MCCDP-12X6.5MM	MC CERVICAL DISC PROTHESIS 12X6.5MM	TI6AL4V-ELI	12x6.5mm
MCCDP-12X07MM	MC CERVICAL DISC PROTHESIS 12X07MM	TI6AL4V-ELI	12x7mm
MCCDP-12X7.5MM	MC CERVICAL DISC PROTHESIS 12X7.5MM	TI6AL4V-ELI	12x7.5mm
MCCDP-12X08MM	MC CERVICAL DISC PROTHESIS 12X08MM	TI6AL4V-ELI	12x8mm
MCCDP-12X8.5MM	MC CERVICAL DISC PROTHESIS 12X8.5MM	TI6AL4V-ELI	12x8.5mm
MCCDP-13X04MM	MC CERVICAL DISC PROTHESIS 13X04MM	TI6AL4V-ELI	13x4mm
MCCDP-13X4.5MM	MC CERVICAL DISC PROTHESIS 13X4.5MM	TI6AL4V-ELI	13x4.5mm
MCCDP-13X05MM	MC CERVICAL DISC PROTHESIS 13X05MM	TI6AL4V-ELI	13x5mm
MCCDP-13X5.5MM	MC CERVICAL DISC PROTHESIS 13X5.5MM	TI6AL4V-ELI	13x5.5mm
MCCDP-13X06MM	MC CERVICAL DISC PROTHESIS 13X06MM	TI6AL4V-ELI	13x6mm
MCCDP-13X6.5MM	MC CERVICAL DISC PROTHESIS 13X6.5MM	TI6AL4V-ELI	13x6.5mm
MCCDP-13X07MM	MC CERVICAL DISC PROTHESIS 13X07MM	TI6AL4V-ELI	13x7mm
MCCDP-13X7.5MM	MC CERVICAL DISC PROTHESIS 13X7.5MM	TI6AL4V-ELI	13x7.5mm
MCCDP-13X08MM	MC CERVICAL DISC PROTHESIS 13X08MM	TI6AL4V-ELI	13x8mm
MCCDP-13X8.5MM	MC CERVICAL DISC PROTHESIS 13X8.5MM	TI6AL4V-ELI	13x8.5mm
MCCDP-14X04MM	MC CERVICAL DISC PROTHESIS 14X04MM	TI6AL4V-ELI	14x4mm
MCCDP-14X4.5MM	MC CERVICAL DISC PROTHESIS 14X4.5MM	TI6AL4V-ELI	14x4.5mm
MCCDP-14X05MM	MC CERVICAL DISC PROTHESIS 14X05MM	TI6AL4V-ELI	14x5mm
MCCDP-14X5.5MM	MC CERVICAL DISC PROTHESIS 14X5.5MM	TI6AL4V-ELI	14x5.5mm
MCCDP-14X06MM	MC CERVICAL DISC PROTHESIS 14X06MM	TI6AL4V-ELI	14x6mm

# Anterior


## Cervical Prosthesis-Cervical Peek Cage Combine Set



MCCDP-14X6.5MM	MC CERVICAL DISC PROTHESIS 14X6.5MM	TI6AL4V-ELI	14x6.5mm
MCCDP-14X07MM	MC CERVICAL DISC PROTHESIS 14X07MM	TI6AL4V-ELI	14x7mm
MCCDP-14X7.5MM	MC CERVICAL DISC PROTHESIS 14X7.5MM	TI6AL4V-ELI	14x7.5mm
MCCDP-14X08MM	MC CERVICAL DISC PROTHESIS 14X08MM	TI6AL4V-ELI	14x8mm
MCCDP-14X8.5MM	MC CERVICAL DISC PROTHESIS 14X8.5MM	TI6AL4V-ELI	14x8.5mm
MCCDP-15X04MM	MC CERVICAL DISC PROTHESIS 15X04MM	TI6AL4V-ELI	15x4mm
MCCDP-15X4.5MM	MC CERVICAL DISC PROTHESIS 15X4.5MM	TI6AL4V-ELI	15x4.5mm
MCCDP-15X05MM	MC CERVICAL DISC PROTHESIS 15X05MM	TI6AL4V-ELI	15x5mm
MCCDP-15X5.5MM	MC CERVICAL DISC PROTHESIS 15X5.5MM	TI6AL4V-ELI	15x5.5mm
MCCDP-15X06MM	MC CERVICAL DISC PROTHESIS 15X06MM	TI6AL4V-ELI	15x6mm
MCCDP-15X6.5MM	MC CERVICAL DISC PROTHESIS 15X6.5MM	TI6AL4V-ELI	15x6.5mm
MCCDP-15X07MM	MC CERVICAL DISC PROTHESIS 15X07MM	TI6AL4V-ELI	15x7mm
MCCDP-15X7.5M	MC CERVICAL DISC PROTHESIS 15X7.5MM	TI6AL4V-ELI	15x7.5mm
MCCDP-15X08MM	MC CERVICAL DISC PROTHESIS 15X08MM	TI6AL4V-ELI	15x8mm
MCCDP-15X8.5MM	MC CERVICAL DISC PROTHESIS 15X8.5MM	TI6AL4V-ELI	15x8.5mm
MCCDP-16X04MM	MC CERVICAL DISC PROTHESIS 16X04MM	TI6AL4V-ELI	16x4mm
MCCDP-16X4.5MM	MC CERVICAL DISC PROTHESIS 16X4.5MM	TI6AL4V-ELI	16x4.5mm
MCCDP-16X05MM	MC CERVICAL DISC PROTHESIS 16X05MM	TI6AL4V-ELI	16x5mm
MCCDP-16X5.5MM	MC CERVICAL DISC PROTHESIS 16X5.5MM	TI6AL4V-ELI	16x5.5mm
MCCDP-16X06MM	MC CERVICAL DISC PROTHESIS 16X06MM	TI6AL4V-ELI	16x6mm
MCCDP-16X6.5MM	MC CERVICAL DISC PROTHESIS 16X6.5MM	TI6AL4V-ELI	16x6.5mm
MCCDP-16X07MM	MC CERVICAL DISC PROTHESIS 16X07MM	TI6AL4V-ELI	16x7mm
MCCDP-16X7.5M	MC CERVICAL DISC PROTHESIS 16X7.5MM	TI6AL4V-ELI	16x7.5mm
MCCDP-16X08MM	MC CERVICAL DISC PROTHESIS 16X08MM	TI6AL4V-ELI	16x8mm
MCCDP-16X8.5MM	MC CERVICAL DISC PROTHESIS 16X8.5MM	TI6AL4V-ELI	16x8.5mm

# Corpectomy

Cervical - Lumbar System



CORPECTOMY  
HOLDER



CORPECTOMY  
MEASURE



CORPECTOMY  
NUT DRIVER

# Corpectomy

Cervical - Lumbar System



**CORPECTOMY  
PUSHER CERVICAL**



**STRAIGHT HANDLE**



**CORPECTOMY  
PUSHER LUMBAR**



**T HANDLE**



**CORPECTOMY  
TURN CERVICAL**

# Corpectomy System



MCLCC-18S	MC LUMBAR CORPECTOMY CAGE 18MM(SHORT)	TI6AL4V-ELI	Ø18mm
MCLCC-18M	MC LUMBAR CORPECTOMY CAGE 18MM(MEDIUM)	TI6AL4V-ELI	Ø18mm
MCLCC-18L	MC LUMBAR CORPECTOMY CAGE 18MM(LONG)	TI6AL4V-ELI	Ø18mm
MCLCC-22S	MC LUMBAR CORPECTOMY CAGE 22MM(SHORT)	TI6AL4V-ELI	Ø22mm
MCLCC-22M	MC LUMBAR CORPECTOMY CAGE 22MM(MEDIUM)	TI6AL4V-ELI	Ø22mm
MCLCC-22L	MC LUMBAR CORPECTOMY CAGE 22MM(LONG)	TI6AL4V-ELI	Ø22mm
MCLCC-25S	MC LUMBAR CORPECTOMY CAGE 25MM(SHORT)	TI6AL4V-ELI	Ø25mm
MCLCC-25M	MC LUMBAR CORPECTOMY CAGE 25MM(MEDIUM)	TI6AL4V-ELI	Ø25mm
MCLCC-25L	MC LUMBAR CORPECTOMY CAGE 25MM(LONG)	TI6AL4V-ELI	Ø25mm



# Corpectomy

Cervical - Lumbar System



MCCCM 08-80	MC CORPECTOMY MESH Ø8 x 80mm	Ti6AL4V-ELI	Ø8 x 80mm
MCCCM 10-80	MC CORPECTOMY MESH Ø10 x 80mm	Ti6AL4V-ELI	Ø10 x 80mm
MCCCM 11-80	MC CORPECTOMY MESH Ø11 x 80mm	Ti6AL4V-ELI	Ø11 x 80mm
MCCCM 12-80	MC CORPECTOMY MESH Ø12 x 80mm	Ti6AL4V-ELI	Ø12 x 80mm
MCCCM 13-80	MC CORPECTOMY MESH Ø13 x 80mm	Ti6AL4V-ELI	Ø13 x 80mm
MCCCM 14-80	MC CORPECTOMY MESH Ø14 x 80mm	Ti6AL4V-ELI	Ø14 x 80mm
MCCCM 15-80	MC CORPECTOMY MESH Ø15 x 80mm	Ti6AL4V-ELI	Ø15 x 80mm
MCCCM 16-80	MC CORPECTOMY MESH Ø16 x 80mm	Ti6AL4V-ELI	Ø16 x 80mm
MCCCM 17-80	MC CORPECTOMY MESH Ø17 x 80mm	Ti6AL4V-ELI	Ø17 x 80mm
MCCCM 18-80	MC CORPECTOMY MESH Ø18 x 80mm	Ti6AL4V-ELI	Ø18 x 80mm
MCCCM 19-80	MC CORPECTOMY MESH Ø19 x 80mm	Ti6AL4V-ELI	Ø19 x 80mm
MCCCM 20-80	MC CORPECTOMY MESH Ø20 x 80mm	Ti6AL4V-ELI	Ø20 x 80mm
MCCCM 21-80	MC CORPECTOMY MESH Ø21 x 80mm	Ti6AL4V-ELI	Ø21 x 80mm
MCCCM 22-80	MC CORPECTOMY MESH Ø22 x 80mm	Ti6AL4V-ELI	Ø22 x 80mm
MCCCM 23-80	MC CORPECTOMY MESH Ø23 x 80mm	Ti6AL4V-ELI	Ø23 x 80mm
MCCCM 24-80	MC CORPECTOMY MESH Ø24 x 80mm	Ti6AL4V-ELI	Ø24 x 80mm
MCCCM 25-80	MC CORPECTOMY MESH Ø25 x 80mm	Ti6AL4V-ELI	Ø25 x 80mm

# Corpectomy

Cervical - Lumbar System



MCCCC-12MM(SHORT)	MC CERVICAL CORPECTOMY IMPLANT	Ti6AL4V-ELI	Ø 12 mm
MCCCC-12MM(MEDIUM)	MC CERVICAL CORPECTOMY IMPLANT	Ti6AL4V-ELI	Ø 12 mm
MCCCC-12MM(LONG)	MC CERVICAL CORPECTOMY IMPLANT	Ti6AL4V-ELI	Ø 12 mm
MCCCC-14MM(SHORT)	MC CERVICAL CORPECTOMY IMPLANT	Ti6AL4V-ELI	Ø 14 mm
MCCCC-14MM(MEDIUM)	MC CERVICAL CORPECTOMY IMPLANT	Ti6AL4V-ELI	Ø 14 mm
MCCCC-14MM(LONG)	MC CERVICAL CORPECTOMY IMPLANT	Ti6AL4V-ELI	Ø 14 mm
MCCCC-16MM(SHORT)	MC CERVICAL CORPECTOMY IMPLANT	Ti6AL4V-ELI	Ø 16 mm
MCCCC-16MM(MEDIUM)	MC CERVICAL CORPECTOMY IMPLANT	Ti6AL4V-EL	Ø 16 mm
MCCCC-16MM(LONG)	MC CERVICAL CORPECTOMY IMPLANT	Ti6AL4V-EL	Ø 16 mm

# Anterior Cervical Plate System



CERVICAL PLATE SCREW DRIVER



CERVICAL PLATE HOLDER



CERVICAL PLATE BENDER




CERVICAL BONE AWL




CERVICAL PLATE LOCKING DRIVER

# Anterior Cervical Plate System



MCCPS-3.5X12	MC CERVICAL PLATE SCREW 3.5X12 MM	TI6AL4V-ELI	3.5x12mm
MCCPS-3.5X14	MC CERVICAL PLATE SCREW 3.5X14 MM	TI6AL4V-ELI	3.5x14mm
MCCPS-3.5X16	MC CERVICAL PLATE SCREW 3.5X16 MM	TI6AL4V-ELI	3.5x16mm
MCCPS-3.5X18	MC CERVICAL PLATE SCREW 3.5X18 MM	TI6AL4V-ELI	3.5x18mm
MCCPS-3.5X20	MC CERVICAL PLATE SCREW 3.5X20 MM	TI6AL4V-ELI	3.5x20mm
MCCPS-4X12	MC CERVICAL PLATE SCREW 4X12 MM	TI6AL4V-ELI	4x12mm
MCCPS-4X14	MC CERVICAL PLATE SCREW 4X14 MM	TI6AL4V-ELI	4x14mm
MCCPS-4X16	MC CERVICAL PLATE SCREW 4X16 MM	TI6AL4V-ELI	4x16mm
MCCPS-4X18	MC CERVICAL PLATE SCREW 4X18 MM	TI6AL4V-ELI	4x18mm
MCCPS-4X20	MC CERVICAL PLATE SCREW 4X20 MM	TI6AL4V-ELI	4x20mm
MCCPS-4.5X12	MC CERVICAL PLATE SCREW 4.5X12 MM	TI6AL4V-ELI	4.5x12mm
MCCPS-4.5X14	MC CERVICAL PLATE SCREW 4.5X14 MM	TI6AL4V-ELI	4.5x14mm
MCCPS-4.5X16	MC CERVICAL PLATE SCREW 4.5X16 MM	TI6AL4V-ELI	4.5x16mm
MCCPS-4.5X18	MC CERVICAL PLATE SCREW 4.5X18 MM	TI6AL4V-ELI	4.5x18mm
MCCPS-4.5X20	MC CERVICAL PLATE SCREW 4.5X20 MM	TI6AL4V-ELI	4.5x20mm

# Craniomaxillofacial Implant System



# Craniomaxillofacial Implant System




**BENDER**



**PLATE CUTTER**




**FORCEPS**



**SCREW BOX**



**HANDLE SCREWDRIVER**



**10x10/10x12 MESH**

# Craniomaxillofacial Implant System



MC.MTAP.200-02	0,5 Micro angled T plate right hole 2	TI6AL4V-ELI
MC.MTAP.200-03	0,5 Micro angled T plate right hole 3	TI6AL4V-ELI
MC.MTAP.200-04	0,5 Micro angled T plate right hole 4	TI6AL4V-ELI
MC.MTAP.200-05	0,5 Micro angled T plate right hole 5	TI6AL4V-ELI
MC.MTAP.200-06	0,5 Micro angled T plate right hole 6	TI6AL4V-ELI
MC.MTAP.200-07	0,5 Micro angled T plate right hole 7	TI6AL4V-ELI
MC.MTAP.200-08	0,5 Micro angled T plate right hole 8	TI6AL4V-ELI
MC.MTAP.201-02	0,5 Micro angled T plate left hole 2	TI6AL4V-ELI
MC.MTAP.201-03	0,5 Micro angled T plate left hole 3	TI6AL4V-ELI
MC.MTAP.201-04	0,5 Micro angled T plate left hole 4	TI6AL4V-ELI
MC.MTAP.201-05	0,5 Micro angled T plate left hole 5	TI6AL4V-ELI
MC.MTAP.201-06	0,5 Micro angled T plate left hole 6	TI6AL4V-ELI
MC.MTAP.201-07	0,5 Micro angled T plate left hole 7	TI6AL4V-ELI
MC.MTAP.201-08	0,5 Micro angled T plate left hole 8	TI6AL4V-ELI
MC.MTP.200-02	0,5 Micro T plate hole 2	TI6AL4V-ELI
MC.MTP.200-03	0,5 Micro T plate hole 3	TI6AL4V-ELI
MC.MTP.200-04	0,5 Micro T plate hole 4	TI6AL4V-ELI
MC.MTP.200-05	0,5 Micro T plate hole 5	TI6AL4V-ELI
MC.MTP.200-06	0,5 Micro T plate hole 6	TI6AL4V-ELI
MC.MTP.200-07	0,5 Micro T plate hole 7	TI6AL4V-ELI
MC.MTP.200-08	0,5 Micro T plate hole 8	TI6AL4V-ELI
MC.MLAP.200-02	0,5 Micro Angled L plate right hole 2	TI6AL4V-ELI
MC.MLAP.200-03	0,5 Micro Angled L plate right hole 3	TI6AL4V-ELI
MC.MLAP.200-04	0,5 Micro Angled L plate right hole 4	TI6AL4V-ELI
MC.MLAP.201-02	0,5 Micro Angled L plate left hole 2	TI6AL4V-ELI
MC.MLAP.201-03	0,5 Micro Angled L plate left hole 3	TI6AL4V-ELI
MC.MLP.200-02	0,5 Micro L plate right hole 2	TI6AL4V-ELI
MC.MLP.200-03	0,5 Micro L plate right hole 3	TI6AL4V-ELI
MC.MLP.200-04	0,5 Micro L plate right hole 4	TI6AL4V-ELI
MC.MLP.201-02	0,5 Micro L plate left hole 2	TI6AL4V-ELI
MC.MLP.201-03	0,5 Micro L plate left hole 3	TI6AL4V-ELI
MC.MLP.201-04	0,5 Micro L plate left hole 4	TI6AL4V-ELI
MC.MDTP.200-02	0,5 Micro Narrow T Plate hole 2	TI6AL4V-ELI

# Craniomaxillofacial Implant System



MC.MDTP.200-03	0,5 Micro Narrow T Plate hole 3	Ti6AL4V-ELI
MC.MDTP.200-04	0,5 Micro Narrow T Plate hole 4	Ti6AL4V-ELI
MC.MDTP.200-05	0,5 Micro Narrow T Plate hole 5	Ti6AL4V-ELI
MC.MDTP.200-06	0,5 Micro Narrow T Plate hole 6	Ti6AL4V-ELI
MC.MDTP.200-07	0,5 Micro Narrow T Plate hole 7	Ti6AL4V-ELI
MC.MDTP.200-08	0,5 Micro Narrow T Plate hole 8	Ti6AL4V-ELI
MC.MYP.200-02	0.5 Micro Y Plate Hole 2	Ti6AL4V-ELI
MC.MYP.200-03	0.5 Micro Y Plate Hole 3	Ti6AL4V-ELI
MC.MYP.200-04	0.5 Micro Y Plate Hole 5	Ti6AL4V-ELI
MC.MDP.200-04	0.5 Micro Straight Plate 4	Ti6AL4V-ELI
MC.MDP.200-05	0.5 Micro Straight Plate 5	Ti6AL4V-ELI
MC.MDP.200-06	0.5 Micro Straight Plate 6	Ti6AL4V-ELI
MC.MDP.200-07	0.5 Micro Straight Plate 7	Ti6AL4V-ELI
MC.MDP.200-08	0.5 Micro Straight Plate 8	Ti6AL4V-ELI
MC.MDP.200-09	0.5 Micro Straight Plate 9	Ti6AL4V-ELI
MC.MDP.200-10	0.5 Micro Straight Plate 10	Ti6AL4V-ELI
MC.MDP.200-11	0.5 Micro Straight Plate 11	Ti6AL4V-ELI
MC.MDP.200-12	0.5 Micro Straight Plate 12	Ti6AL4V-ELI
MC.MDP.200-13	0.5 Micro Straight Plate 13	Ti6AL4V-ELI
MC.MDP.200-14	0.5 Micro Straight Plate 14	Ti6AL4V-ELI
MC.MATP.300-02	0,8 Mini Angled T plate right hole 2	Ti6AL4V-ELI
MC.MATP.300-03	0,8 Mini Angled T plate right hole 3	Ti6AL4V-ELI
MC.MATP.300-04	0,8 Mini Angled T plate right hole 4	Ti6AL4V-ELI
MC.MATP.300-05	0,8 Mini Angled T plate right hole 5	Ti6AL4V-ELI
MC.MATP.300-06	0,8 Mini Angled T plate right hole 6	Ti6AL4V-ELI
MC.MATP.300-07	0,8 Mini Angled T plate right hole 7	Ti6AL4V-ELI
MC.MATP.300-08	0,8 Mini Angled T plate right hole 8	Ti6AL4V-ELI
MC.MATP.301-02	0,8 Mini Angled T plate left hole 2	Ti6AL4V-ELI
MC.MATP.301-03	0,8 Mini Angled T plate left hole 3	Ti6AL4V-ELI
MC.MATP.301-04	0,8 Mini Angled T plate left hole 4	Ti6AL4V-ELI
MC.MATP.301-05	0,8 Mini Angled T plate left hole 5	Ti6AL4V-ELI
MC.MATP.301-06	0,8 Mini Angled T plate left hole 6	Ti6AL4V-ELI
MC.MATP.301-07	0,8 Mini Angled T plate left hole 7	Ti6AL4V-ELI




# Craniomaxillofacial Implant System



MC.MATP.301-08	0,8 Mini Angled T plate left hole 8	TI6AL4V-ELI
MC.MTP.300-02	0,8 Mini T plate Hole 2	TI6AL4V-ELI
MC.MTP.300-03	0,8 Mini T plate Hole 3	TI6AL4V-ELI
MC.MTP.300-04	0,8 Mini T plate Hole 4	TI6AL4V-ELI
MC.MTP.300-05	0,8 Mini T plate Hole 5	TI6AL4V-ELI
MC.MTP.300-06	0,8 Mini T plate Hole 6	TI6AL4V-ELI
MC.MTP.300-07	0,8 Mini T plate Hole 7	TI6AL4V-ELI
MC.MTP.300-08	0,8 Mini T plate Hole 8	TI6AL4V-ELI
MC.MALP.300-02	0,8 Mini Angled L plate Right hole 2	TI6AL4V-ELI
MC.MALP.300-03	0,8 Mini Angled L plate Right hole 3	TI6AL4V-ELI
MC.MALP.300-04	0,8 Mini Angled L plate Right hole 4	TI6AL4V-ELI
MC.MALP.301-02	0,8 Mini Angled L plate Left hole 2	TI6AL4V-ELI
MC.MALP.301-03	0,8 Mini Angled L plate Left hole 3	TI6AL4V-ELI
MC.MALP.301-04	0,8 Mini Angled L plate Left hole 4	TI6AL4V-ELI
MC.MLP.300-02	0,8 Mini L plate right hole 2	TI6AL4V-ELI
MC.MLP.300-03	0,8 Mini L plate right hole 3	TI6AL4V-ELI
MC.MLP.300-04	0,8 Mini L plate right hole 4	TI6AL4V-ELI
MC.MLP.301-02	0,8 Mini L plate left hole 2	TI6AL4V-ELI
MC.MLP.301-03	0,8 Mini L plate left hole 3	TI6AL4V-ELI
MC.MLP.301-04	0,8 Mini L plate left hole 4	TI6AL4V-ELI
MC.MDTP.300-02	0,8 Mini Narrow T plate hole 2	TI6AL4V-ELI
MC.MDTP.300-03	0,8 Mini Narrow T plate hole 3	TI6AL4V-ELI
MC.MDTP.300-04	0,8 Mini Narrow T plate hole 4	TI6AL4V-ELI
MC.MDTP.300-05	0,8 Mini Narrow T plate hole 5	TI6AL4V-ELI
MC.MDTP.300-06	0,8 Mini Narrow T plate hole 6	TI6AL4V-ELI
MC.MDTP.300-07	0,8 Mini Narrow T plate hole 7	TI6AL4V-ELI
MC.MDTP.300-08	0,8 Mini Narrow T plate hole 8	TI6AL4V-ELI
MC.MYP.300-02	0,8 Mini Y Plate Hole 2	TI6AL4V-ELI
MC.MYP.300-03	0,8 Mini Y Plate Hole 3	TI6AL4V-ELI
MC.MYP.300-04	0,8 Mini Y Plate Hole 4	TI6AL4V-ELI
MC.MDP.300-04	0,8 Mini Straight Plate Hole 4	TI6AL4V-ELI
MC.MDP.300-05	0,8 Mini Straight Plate Hole 5	TI6AL4V-ELI
MC.MDP.300-06	0,8 Mini Straight Plate Hole 6	TI6AL4V-ELI

# Craniomaxillofacial Implant System



MC.MDP.300-07	0.8 Mini Straight Plate Hole 7	Ti6AL4V-ELI
MC.MDP.300-08	0.8 Mini Straight Plate Hole 8	Ti6AL4V-ELI
MC.MDP.300-09	0.8 Mini Straight Plate Hole 9	Ti6AL4V-ELI
MC.MDP.300-10	0.8 Mini Straight Plate Hole 10	Ti6AL4V-ELI
MC.MDP.300-11	0.8 Mini Straight Plate Hole 11	Ti6AL4V-ELI
MC.MDP.300-12	0.8 Mini Straight Plate Hole 12	Ti6AL4V-ELI
MC.MDP.300-13	0.8 Mini Straight Plate Hole 13	Ti6AL4V-ELI
MC.MDP.300-14	0.8 Mini Straight Plate Hole 14	Ti6AL4V-ELI
MC.MATP.400-02.1	1,0 Mini Angled T plate right hole 2	Ti6AL4V-ELI
MC.MATP.400-03.1	1,0 Mini Angled T plate right hole 3	Ti6AL4V-ELI
MC.MATP.400-04.1	1,0 Mini Angled T plate right hole 4	Ti6AL4V-ELI
MC.MATP.400-05.1	1,0 Mini Angled T plate right hole 5	Ti6AL4V-ELI
MC.MATP.400-06.1	1,0 Mini Angled T plate right hole 6	Ti6AL4V-ELI
MC.MATP.400-07.1	1,0 Mini Angled T plate right hole 7	Ti6AL4V-ELI
MC.MATP.400-08.1	1,0 Mini Angled T plate right hole 8	Ti6AL4V-ELI
MC.MATP.401-02	1,0 Mini Angled T plate left hole 2	Ti6AL4V-ELI
MC.MATP.401-03	1,0 Mini Angled T plate left hole 3	Ti6AL4V-ELI
MC.MATP.401-04	1,0 Mini Angled T plate left hole 4	Ti6AL4V-ELI
MC.MATP.401-05	1,0 Mini Angled T plate left hole 5	Ti6AL4V-ELI
MC.MATP.401-06	1,0 Mini Angled T plate left hole 6	Ti6AL4V-ELI
MC.MATP.401-07	1,0 Mini Angled T plate left hole 7	Ti6AL4V-ELI
MC.MATP.401-08	1,0 Mini Angled T plate left hole 8	Ti6AL4V-ELI
MC.MTP.300-02.1	1,0 Mini T plate hole 2	Ti6AL4V-ELI
MC.MTP.300-03.1	1,0 Mini T plate hole 3	Ti6AL4V-ELI
MC.MTP.300-04.1	1,0 Mini T plate hole 4	Ti6AL4V-ELI
MC.MTP.300-05.1	1,0 Mini T plate hole 5	Ti6AL4V-ELI
MC.MTP.300-06.1	1,0 Mini T plate hole 6	Ti6AL4V-ELI
MC.MTP.300-07.1	1,0 Mini T plate hole 7	Ti6AL4V-ELI
MC.MTP.300-08.1	1,0 Mini T plate hole 8	Ti6AL4V-ELI
MC.MALP.300-02.1	1,0 Mini Angled L plate right hole 2	Ti6AL4V-ELI
MC.MALP.300-03.1	1,0 Mini Angled L plate right hole 3	Ti6AL4V-ELI
MC.MALP.300-04.1	1,0 Mini Angled L plate right hole 4	Ti6AL4V-ELI
MC.MALP.301-02.1	1,0 Mini Angled L plate left hole 2	Ti6AL4V-ELI

# Craniomaxillofacial Implant System






MC.MALP.301-03.1	1,0 Mini Angled L plate left hole 3	TI6AL4V-ELI
MC.MALP.301-04.1	1,0 Mini Angled L plate left hole 4	TI6AL4V-ELI
MC.MLP.300-02.1	1,0 Mini L plate right hole 2	TI6AL4V-ELI
MC.MLP.300-03.1	1,0 Mini L plate right hole 3	TI6AL4V-ELI
MC.MLP.300-04.1	1,0 Mini L plate right hole 4	TI6AL4V-ELI
MC.MLP.301-02	1,0 Mini L plate left hole 2	TI6AL4V-ELI
MC.MLP.301-03	1,0 Mini L plate left hole 3	TI6AL4V-ELI
MC.MLP.301-04	1,0 Mini L plate left hole 4	TI6AL4V-ELI
MC.MDTP.300-02.1	1,0 Mini Narrow T plate hole 2	TI6AL4V-ELI
MC.MDTP.300-03.1	1,0 Mini Narrow T plate hole 3	TI6AL4V-ELI
MC.MDTP.300-04.1	1,0 Mini Narrow T plate hole 4	TI6AL4V-ELI
MC.MDTP.300-05.1	1,0 Mini Narrow T plate hole 5	TI6AL4V-ELI
MC.MDTP.300-06.1	1,0 Mini Narrow T plate hole 6	TI6AL4V-ELI
MC.MDTP.300-07.1	1,0 Mini Narrow T plate hole 7	TI6AL4V-ELI
MC.MDTP.300-08.1	1,0 Mini Narrow T plate hole 8	TI6AL4V-ELI
MC.MYPD.300-02	1,0 Mini Y Plate hole 2	TI6AL4V-ELI
MC.MYPD.300-03	1,0 Mini Y Plate hole 3	TI6AL4V-ELI
MC.MYPD.300-04	1,0 Mini Y Plate hole 4	TI6AL4V-ELI
MC.MDP.300-04.1	1,0 Mini Straight plate 4	TI6AL4V-ELI
MC.MDP.300-05.1	1,0 Mini Straight plate 5	TI6AL4V-ELI
MC.MDP.300-06.1	1,0 Mini Straight plate 6	TI6AL4V-ELI
MC.MDP.300-07.1	1,0 Mini Straight plate 7	TI6AL4V-ELI
MC.MDP.300-08.1	1,0 Mini Straight plate 8	TI6AL4V-ELI
MC.MDP.300-09.1	1,0 Mini Straight plate 9	TI6AL4V-ELI
MC.MDP.300-10.1	1,0 Mini Straight plate 10	TI6AL4V-ELI
MC.MDP.300-11.1	1,0 Mini Straight plate 11	TI6AL4V-ELI
MC.MDP.300-12.1	1,0 Mini Straight plate 12	TI6AL4V-ELI
MC.MDP.300-13.1	1,0 Mini Straight plate 13	TI6AL4V-ELI
MC.MDP.300-14.1	1,0 Mini Straight plate 14	TI6AL4V-ELI
MC.MATP.400-02	1,5 Mini Angled T plate right hole 2	TI6AL4V-ELI
MC.MATP.400-03	1,5 Mini Angled T plate right hole 3	TI6AL4V-ELI
MC.MATP.400-04	1,5 Mini Angled T plate right hole 4	TI6AL4V-ELI
MC.MATP.400-05	1,5 Mini Angled T plate right hole 5	TI6AL4V-ELI

# Craniomaxillofacial Implant System



MC.MATP.400-06	1,5 Mini Angled T plate right hole 6	TI6AL4V-ELI
MC.MATP.400-07	1,5 Mini Angled T plate right hole 7	TI6AL4V-ELI
MC.MATP.400-08	1,5 Mini Angled T plate right hole 8	TI6AL4V-ELI
MC.MATP.401-02	1,5 Mini Angled T plate left hole 2	TI6AL4V-ELI
MC.MATP.401-03	1,5 Mini Angled T plate left hole 3	TI6AL4V-ELI
MC.MATP.401-04	1,5 Mini Angled T plate left hole 4	TI6AL4V-ELI
MC.MATP.401-05	1,5 Mini Angled T plate left hole 5	TI6AL4V-ELI
MC.MATP.401-06	1,5 Mini Angled T plate left hole 6	TI6AL4V-ELI
MC.MATP.401-07	1,5 Mini Angled T plate left hole 7	TI6AL4V-ELI
MC.MATP.401-08	1,5 Mini Angled T plate left hole 8	TI6AL4V-ELI
MC.MTP.400-02	1,5 Mini T plate hole 2	TI6AL4V-ELI
MC.MTP.400-03	1,5 Mini T plate hole 3	TI6AL4V-ELI
MC.MTP.400-04	1,5 Mini T plate hole 4	TI6AL4V-ELI
MC.MTP.400-05	1,5 Mini T plate hole 5	TI6AL4V-ELI
MC.MTP.400-06	1,5 Mini T plate hole 6	TI6AL4V-ELI
MC.MTP.400-07	1,5 Mini T plate hole 7	TI6AL4V-ELI
MC.MTP.400-08	1,5 Mini T plate hole 8	TI6AL4V-ELI
MC.MALP.400-02	1,5 Mini Angled L plate right hole 2	TI6AL4V-ELI
MC.MALP.400-03	1,5 Mini Angled L plate right hole 3	TI6AL4V-ELI
MC.MALP.400-04	1,5 Mini Angled L plate right hole 4	TI6AL4V-ELI
MC.MALP.401-02	1,5 Mini Angled L plate left hole 2	TI6AL4V-ELI
MC.MALP.401-03	1,5 Mini Angled L plate left hole 3	TI6AL4V-ELI
MC.MALP.401-04	1,5 Mini Angled L plate left hole 4	TI6AL4V-ELI
MC.MLP.400-02	1,5 Mini L plate right hole 2	TI6AL4V-ELI
MC.MLP.400-03	1,5 Mini L plate right hole 3	TI6AL4V-ELI
MC.MLP.400-04	1,5 Mini L plate right hole 4	TI6AL4V-ELI
MC.MLP.401-02	1,5 Mini L plate left hole 2	TI6AL4V-ELI
MC.MLP.401-03	1,5 Mini L plate left hole 3	TI6AL4V-ELI
MC.MLP.401-04	1,5 Mini L plate left hole 4	TI6AL4V-ELI
MC.MDTP.400-02	1,5 Mini Narrow T plate hole 2	TI6AL4V-ELI
MC.MDTP.400-03	1,5 Mini Narrow T plate hole 3	TI6AL4V-ELI
MC.MDTP.400-04	1,5 Mini Narrow T plate hole 4	TI6AL4V-ELI
MC.MDTP.400-05	1,5 Mini Narrow T plate hole 5	TI6AL4V-ELI

# Craniomaxillofacail Implant System



MC.MDTP.400-06	1,5 Mini Narrow T plate hole 6	Ti6AL4V-ELI
MC.MDTP.400-07	1,5 Mini Narrow T plate hole 7	Ti6AL4V-ELI
MC.MDTP.400-08	1,5 Mini Narrow T plate hole 8	Ti6AL4V-ELI
MC.MYP.400-02	1,5 Mini Y Plate Hole 2	Ti6AL4V-ELI
MC.MYP.400-03	1,5 Mini Y Plate Hole 3	Ti6AL4V-ELI
MC.MYP.400-04	1,5 Mini Y Plate Hole 4	Ti6AL4V-ELI
MC.MDP.400-04	1,5 Mini Straight Plate 4	Ti6AL4V-ELI
MC.MDP.400-05	1,5 Mini Straight Plate 5	Ti6AL4V-ELI
MC.MDP.400-06	1,5 Mini Straight Plate 6	Ti6AL4V-ELI
MC.MDP.400-07	1,5 Mini Straight Plate 7	Ti6AL4V-ELI
MC.MDP.400-08	1,5 Mini Straight Plate 8	Ti6AL4V-ELI
MC.MDP.400-09	1,5 Mini Straight Plate 9	Ti6AL4V-ELI
MC.MDP.400-10	1,5 Mini Straight Plate 10	Ti6AL4V-ELI
MC.MDP.400-11	1,5 Mini Straight Plate 11	Ti6AL4V-ELI
MC.MDP.400-12	1,5 Mini Straight Plate 12	Ti6AL4V-ELI
MC.MDP.400-13	1,5 Mini Straight Plate 13	Ti6AL4V-ELI
MC.MDP.400-14	1,5 Mini Straight Plate 14	Ti6AL4V-ELI
MC.KVKB.100-06	1,5 mm cortical screw square headed 6 mm	Ti6AL4V-ELI
MC.KVKB.100-07	1,5 mm cortical screw square headed 7 mm	Ti6AL4V-ELI
MC.KVKB.100-08	1,5 mm cortical screw square headed 8 mm	Ti6AL4V-ELI
MC.KVKB.100-09	1,5 mm cortical screw square headed 9 mm	Ti6AL4V-ELI
MC.KVKB.100-10	1,5 mm cortical screw square headed 10 mm	Ti6AL4V-ELI
MC.KVKB.100-11	1,5 mm cortical screw square headed 11 mm	Ti6AL4V-ELI
MC.KVKB.100-12	1,5 mm cortical screw square headed 12 mm	Ti6AL4V-ELI
MC.KVKB.100-13	1,5 mm cortical screw square headed 13 mm	Ti6AL4V-ELI
MC.KVKB.100-14	1,5 mm cortical screw square headed 14 mm	Ti6AL4V-ELI
MC.KVKB.100-15	1,5 mm cortical screw square headed 15 mm	Ti6AL4V-ELI
MC.KVKB.100-16	1,5 mm cortical screw square headed 16 mm	Ti6AL4V-ELI
MC.KVKB.100-17	1,5 mm cortical screw square headed 17 mm	Ti6AL4V-ELI
MC.KVKB.100-18	1,5 mm cortical screw square headed 18 mm	Ti6AL4V-ELI
MC.KVKB.100-19	1,5 mm cortical screw square headed 19 mm	Ti6AL4V-ELI
MC.KVKB.100-20	1,5 mm cortical screw square headed 20 mm	Ti6AL4V-ELI
MC.KVKB.100-21	1,5 mm cortical screw square headed 21 mm	Ti6AL4V-ELI

# Craniomaxillofacial Implant System



MC.KVKB.100-22	1,5 mm cortical screw square headed 22 mm	Ti6AL4V-ELI
MC.KVKB.100-23	1,5 mm cortical screw square headed 23 mm	Ti6AL4V-ELI
MC.KVKB.100-24	1,5 mm cortical screw square headed 24 mm	Ti6AL4V-ELI
MC.KVKB.100-25	1,5 mm cortical screw square headed 25 mm	Ti6AL4V-ELI
MC.KVKB.100-26	1,5 mm cortical screw square headed 26 mm	Ti6AL4V-ELI
MC.KKV.101-06	1,5 mm locked cortical screw 6 mm	Ti6AL4V-ELI
MC.KKV.101-07	1,5 mm locked cortical screw 7 mm	Ti6AL4V-ELI
MC.KKV.101-08	1,5 mm locked cortical screw 8 mm	Ti6AL4V-ELI
MC.KKV.101-09	1,5 mm locked cortical screw 9 mm	Ti6AL4V-ELI
MC.KKV.101-10	1,5 mm locked cortical screw 10 mm	Ti6AL4V-ELI
MC.KKV.101-11	1,5 mm locked cortical screw 11 mm	Ti6AL4V-ELI
MC.KKV.101-12	1,5 mm locked cortical screw 12 mm	Ti6AL4V-ELI
MC.KKV.101-13	1,5 mm locked cortical screw 13 mm	Ti6AL4V-ELI
MC.KKV.101-14	1,5 mm locked cortical screw 14 mm	Ti6AL4V-ELI
MC.KKV.101-15	1,5 mm locked cortical screw 15 mm	Ti6AL4V-ELI
MC.KKV.101-16	1,5 mm locked cortical screw 16 mm	Ti6AL4V-ELI
MC.KKV.101-17	1,5 mm locked cortical screw 17 mm	Ti6AL4V-ELI
MC.KKV.101-18	1,5 mm locked cortical screw 18 mm	Ti6AL4V-ELI
MC.KKV.101-19	1,5 mm locked cortical screw 19 mm	Ti6AL4V-ELI
MC.KKV.101-20	1,5 mm locked cortical screw 20 mm	Ti6AL4V-ELI
MC.KKV.101-21	1,5 mm locked cortical screw 21 mm	Ti6AL4V-ELI
MC.KKV.101-22	1,5 mm locked cortical screw 22 mm	Ti6AL4V-ELI
MC.KKV.101-23	1,5 mm locked cortical screw 23 mm	Ti6AL4V-ELI
MC.KKV.101-24	1,5 mm locked cortical screw 24 mm	Ti6AL4V-ELI
MC.KKV.101-25	1,5 mm locked cortical screw 25 mm	Ti6AL4V-ELI
MC.KKV.101-26	1,5 mm locked cortical screw 26 mm	Ti6AL4V-ELI
MC.KV.110-06	1,9 mm cortical screw 6 mm	Ti6AL4V-ELI
MC.KV.110-07	1,9 mm cortical screw 7 mm	Ti6AL4V-ELI
MC.KV.110-08	1,9 mm cortical screw 8 mm	Ti6AL4V-ELI
MC.KV.110-09	1,9 mm cortical screw 9 mm	Ti6AL4V-ELI
MC.KV.110-10	1,9 mm cortical screw 10 mm	Ti6AL4V-ELI
MC.KV.110-11	1,9 mm cortical screw 11 mm	Ti6AL4V-ELI
MC.KV.110-12	1,9 mm cortical screw 12 mm	Ti6AL4V-ELI

# Craniomaxillofacial Implant System



MC.KV.110-13	1,9 mm cortical screw 13 mm	Ti6AL4V-ELI
MC.KV.110-14	1,9 mm cortical screw 14 mm	Ti6AL4V-ELI
MC.KV.110-15	1,9 mm cortical screw 15 mm	Ti6AL4V-ELI
MC.KV.110-16	1,9 mm cortical screw 16 mm	Ti6AL4V-ELI
MC.KV.110-17	1,9 mm cortical screw 17 mm	Ti6AL4V-ELI
MC.KV.110-18	1,9 mm cortical screw 18 mm	Ti6AL4V-ELI
MC.KV.110-19	1,9 mm cortical screw 19 mm	Ti6AL4V-ELI
MC.KV.110-20	1,9 mm cortical screw 20 mm	Ti6AL4V-ELI
MC.KV.110-21	1,9 mm cortical screw 21 mm	Ti6AL4V-ELI
MC.KV.110-22	1,9 mm cortical screw 22 mm	Ti6AL4V-ELI
MC.KV.110-23	1,9 mm cortical screw 23 mm	Ti6AL4V-ELI
MC.KV.110-24	1,9 mm cortical screw 24 mm	Ti6AL4V-ELI
MC.KV.110-25	1,9 mm cortical screw 25 mm	Ti6AL4V-ELI
MC.KV.110-26	1,9 mm cortical screw 26 mm	Ti6AL4V-ELI
MC.KKV.111-06	1,9 mm locked cortical screw 6 mm	Ti6AL4V-ELI
MC.KKV.111-07	1,9 mm locked cortical screw 7mm	Ti6AL4V-ELI
MC.KKV.111-08	1,9 mm locked cortical screw 8 mm	Ti6AL4V-ELI
MC.KKV.111-09	1,9 mm locked cortical screw 9 mm	Ti6AL4V-ELI
MC.KKV.111-10	1,9 mm locked cortical screw 10 mm	Ti6AL4V-ELI
MC.KKV.111-11	1,9 mm locked cortical screw 11 mm	Ti6AL4V-ELI
MC.KKV.111-12	1,9 mm locked cortical screw 12 mm	Ti6AL4V-ELI
MC.KKV.111-13	1,9 mm locked cortical screw 13 mm	Ti6AL4V-ELI
MC.KKV.111-14	1,9 mm locked cortical screw 14 mm	Ti6AL4V-ELI
MC.KKV.111-15	1,9 mm locked cortical screw 15 mm	Ti6AL4V-ELI
MC.KKV.111-16	1,9 mm locked cortical screw 16 mm	Ti6AL4V-ELI
MC.KKV.111-17	1,9 mm locked cortical screw 17 mm	Ti6AL4V-ELI
MC.KKV.111-18	1,9 mm locked cortical screw 18 mm	Ti6AL4V-ELI
MC.KKV.111-19	1,9 mm locked cortical screw 19 mm	Ti6AL4V-ELI
MC.KKV.111-20	1,9 mm locked cortical screw 20 mm	Ti6AL4V-ELI
MC.KKV.111-21	1,9 mm locked cortical screw 21 mm	Ti6AL4V-ELI
MC.KKV.111-22	1,9 mm locked cortical screw 22 mm	Ti6AL4V-ELI
MC.KKV.111-23	1,9 mm locked cortical screw 23 mm	Ti6AL4V-ELI
MC.KKV.111-24	1,9 mm locked cortical screw 24 mm	Ti6AL4V-ELI

# Craniomaxillofacial Implant System

MC.KKV.111-25	1,9 mm locked cortical screw 25 mm	Ti6AL4V-ELI
MC.KKV.111-26	1,9 mm locked cortical screw 26 mm	Ti6AL4V-ELI



# Orthopedic

## Screws

MC.KV.120-08	2,0 mm cortical screw 8 mm	Ti6AL4V-ELI
MC.KV.120-10	2,0 mm cortical screw 10 mm	Ti6AL4V-ELI
MC.KV.120-12	2,0 mm cortical screw 12 mm	Ti6AL4V-ELI
MC.KV.120-14	2,0 mm cortical screw 14 mm	Ti6AL4V-ELI
MC.KV.120-16	2,0 mm cortical screw 16 mm	Ti6AL4V-ELI
MC.KV.120-18	2,0 mm cortical screw 18 mm	Ti6AL4V-ELI
MC.KV.120-20	2,0 mm cortical screw 20 mm	Ti6AL4V-ELI
MC.KV.120-22	2,0 mm cortical screw 22 mm	Ti6AL4V-ELI
MC.KV.120-24	2,0 mm cortical screw 24 mm	Ti6AL4V-ELI
MC.KV.120-26	2,0 mm cortical screw 26 mm	Ti6AL4V-ELI
MC.KV.121-08	2,0 mm locked cortical screw 8 mm	Ti6AL4V-ELI
MC.KV.121-10	2,0 mm locked cortical screw 10 mm	Ti6AL4V-ELI
MC.KV.121-12	2,0 mm locked cortical screw 12 mm	Ti6AL4V-ELI
MC.KV.121-14	2,0 mm locked cortical screw 14 mm	Ti6AL4V-ELI
MC.KV.121-16	2,0 mm locked cortical screw 16 mm	Ti6AL4V-ELI
MC.KV.121-18	2,0 mm locked cortical screw 18 mm	Ti6AL4V-ELI
MC.KV.121-20	2,0 mm locked cortical screw 20 mm	Ti6AL4V-ELI
MC.KV.121-22	2,0 mm locked cortical screw 22 mm	Ti6AL4V-ELI
MC.KV.121-24	2,0 mm locked cortical screw 24 mm	Ti6AL4V-ELI
MC.KV.121-26	2,0 mm locked cortical screw 26 mm	Ti6AL4V-ELI
MC.KV.130-08	2,5 mm cortical screw length 8 mm	Ti6AL4V-ELI
MC.KV.130-10	2,5 mm cortical screw length 10 mm	Ti6AL4V-ELI
MC.KV.130-12	2,5 mm cortical screw length 12 mm	Ti6AL4V-ELI
MC.KV.130-14	2,5 mm cortical screw length 14 mm	Ti6AL4V-ELI
MC.KV.130-16	2,5 mm cortical screw length 16 mm	Ti6AL4V-ELI
MC.KV.130-18	2,5 mm cortical screw length 18 mm	Ti6AL4V-ELI
MC.KV.130-20	2,5 mm cortical screw length 20 mm	Ti6AL4V-ELI
MC.KV.130-22	2,5 mm cortical screw length 22 mm	Ti6AL4V-ELI
MC.KV.130-24	2,5 mm cortical screw length 24 mm	Ti6AL4V-ELI
MC.KV.130-26	2,5 mm cortical screw length 26 mm	Ti6AL4V-ELI
MC.KV.140-10	cortical screw Self Tapping 2,7 mm length 10	Ti6AL4V-ELI
MC.KV.140-12	cortical screw Self Tapping 2,7 mm length 12	Ti6AL4V-ELI

# Orthopedic Screws

MC.KV.140-14	cortical screw Self Tapping 2,7 mm lenght 14	TI6AL4V-ELI
MC.KV.140-16	cortical screw Self Tapping 2,7 mm lenght 16	TI6AL4V-ELI
MC.KV.140-18	cortical screw Self Tapping 2,7 mm lenght 18	TI6AL4V-ELI
MC.KV.140-20	cortical screw Self Tapping 2,7 mm lenght 20	TI6AL4V-ELI
MC.KV.140-22	cortical screw Self Tapping 2,7 mm lenght 22	TI6AL4V-ELI
MC.KV.140-24	cortical screw Self Tapping 2,7 mm lenght 24	TI6AL4V-ELI
MC.KV.140-26	cortical screw Self Tapping 2,7 mm lenght 26	TI6AL4V-ELI
MC.KV.140-28	cortical screw Self Tapping 2,7 mm lenght 28	TI6AL4V-ELI
MC.KV.140-30	cortical screw Self Tapping 2,7 mm lenght 30	TI6AL4V-ELI
MC.KV.140-32	cortical screw Self Tapping 2,7 mm lenght 32	TI6AL4V-ELI
MC.KV.140-34	cortical screw Self Tapping 2,7 mm lenght 34	TI6AL4V-ELI
MC.KV.140-36	cortical screw Self Tapping 2,7 mm lenght 36	TI6AL4V-ELI
MC.KV.140-38	cortical screw Self Tapping 2,7 mm lenght 38	TI6AL4V-ELI
MC.KV.140-40	cortical screw Self Tapping 2,7 mm lenght 40	TI6AL4V-ELI
MC.KV.140-42	cortical screw Self Tapping 2,7 mm lenght 42	TI6AL4V-ELI
MC.KV.140-44	cortical screw Self Tapping 2,7 mm lenght 44	TI6AL4V-ELI
MC.KV.140-46	cortical screw Self Tapping 2,7 mm lenght 46	TI6AL4V-ELI
MC.KV.140-48	cortical screw Self Tapping 2,7 mm lenght 48	TI6AL4V-ELI
MC.KV.140-50	cortical screw Self Tapping 2,7 mm lenght 50	TI6AL4V-ELI

# Orthopedic Screws

MC.KV.150-10	cortical screw Self Tapping 3,5 mm length 10	Ti6AL4V-ELI
MC.KV.150-12	cortical screw Self Tapping 3,5 mm length 12	Ti6AL4V-ELI
MC.KV.150-14	cortical screw Self Tapping 3,5 mm length 14	Ti6AL4V-ELI
MC.KV.150-16	cortical screw Self Tapping 3,5 mm length 16	Ti6AL4V-ELI
MC.KV.150-18	cortical screw Self Tapping 3,5 mm length 18	Ti6AL4V-ELI
MC.KV.150-20	cortical screw Self Tapping 3,5 mm length 20	Ti6AL4V-ELI
MC.KV.150-22	cortical screw Self Tapping 3,5 mm length 22	Ti6AL4V-ELI
MC.KV.150-24	cortical screw Self Tapping 3,5 mm length 24	Ti6AL4V-ELI
MC.KV.150-26	cortical screw Self Tapping 3,5 mm length 26	Ti6AL4V-ELI
MC.KV.150-28	cortical screw Self Tapping 3,5 mm length 28	Ti6AL4V-ELI
MC.KV.150-30	cortical screw Self Tapping 3,5 mm length 30	Ti6AL4V-ELI
MC.KV.150-32	cortical screw Self Tapping 3,5 mm length 32	Ti6AL4V-ELI
MC.KV.150-34	cortical screw Self Tapping 3,5 mm length 34	Ti6AL4V-ELI
MC.KV.150-36	cortical screw Self Tapping 3,5 mm length 36	Ti6AL4V-ELI
MC.KV.150-38	cortical screw Self Tapping 3,5 mm length 38	Ti6AL4V-ELI
MC.KV.150-40	cortical screw Self Tapping 3,5 mm length 40	Ti6AL4V-ELI
MC.KV.150-45	cortical screw Self Tapping 3,5 mm length 45	Ti6AL4V-ELI
MC.KV.150-50	cortical screw Self Tapping 3,5 mm length 50	Ti6AL4V-ELI
MC.KV.160-14	cortical screw Self Tapping 4,0 mm length 14	Ti6AL4V-ELI
MC.KV.160-16	cortical screw Self Tapping 4,0 mm length 16	Ti6AL4V-ELI
MC.KV.160-18	cortical screw Self Tapping 4,0 mm length 18	Ti6AL4V-ELI
MC.KV.160-20	cortical screw Self Tapping 4,0 mm length 20	Ti6AL4V-ELI
MC.KV.160-22	cortical screw Self Tapping 4,0 mm length 22	Ti6AL4V-ELI
MC.KV.160-24	cortical screw Self Tapping 4,0 mm length 24	Ti6AL4V-ELI
MC.KV.160-26	cortical screw Self Tapping 4,0 mm length 26	Ti6AL4V-ELI
MC.KV.160-28	cortical screw Self Tapping 4,0 mm length 28	Ti6AL4V-ELI
MC.KV.160-30	cortical screw Self Tapping 4,0 mm length 30	Ti6AL4V-ELI
MC.KV.160-32	cortical screw Self Tapping 4,0 mm length 32	Ti6AL4V-ELI
MC.KV.160-34	cortical screw Self Tapping 4,0 mm length 34	Ti6AL4V-ELI
MC.KV.160-36	cortical screw Self Tapping 4,0 mm length 36	Ti6AL4V-ELI
MC.KV.160-38	cortical screw Self Tapping 4,0 mm length 38	Ti6AL4V-ELI
MC.KV.160-40	cortical screw Self Tapping 4,0 mm length 40	Ti6AL4V-ELI
MC.KV.160-42	cortical screw Self Tapping 4,0 mm length 42	Ti6AL4V-ELI

# Orthopedic Screws



MC.KV.160-44	cortical screw Self Tapping 4,0 mm length 44	Ti6AL4V-ELI
MC.KV.160-46	cortical screw Self Tapping 4,0 mm length 46	Ti6AL4V-ELI
MC.KV.160-48	cortical screw Self Tapping 4,0 mm length 48	Ti6AL4V-ELI
MC.KV.160-50	cortical screw Self Tapping 4,0 mm length 50	Ti6AL4V-ELI
MC.KV.170-14	cortical screw Self Tapping 4,5 mm length 14	Ti6AL4V-ELI
MC.KV.170-16	cortical screw Self Tapping 4,5 mm length 16	Ti6AL4V-ELI
MC.KV.170-18	cortical screw Self Tapping 4,5 mm length 18	Ti6AL4V-ELI
MC.KV.170-20	cortical screw Self Tapping 4,5 mm length 20	Ti6AL4V-ELI
MC.KV.170-22	cortical screw Self Tapping 4,5 mm length 22	Ti6AL4V-ELI
MC.KV.170-24	cortical screw Self Tapping 4,5 mm length 24	Ti6AL4V-ELI
MC.KV.170-26	cortical screw Self Tapping 4,5 mm length 26	Ti6AL4V-ELI
MC.KV.170-28	cortical screw Self Tapping 4,5 mm length 28	Ti6AL4V-ELI
MC.KV.170-30	cortical screw Self Tapping 4,5 mm length 30	Ti6AL4V-ELI
MC.KV.170-32	cortical screw Self Tapping 4,5 mm length 32	Ti6AL4V-ELI
MC.KV.170-34	cortical screw Self Tapping 4,5 mm length 34	Ti6AL4V-ELI
MC.KV.170-36	cortical screw Self Tapping 4,5 mm length 36	Ti6AL4V-ELI
MC.KV.170-38	cortical screw Self Tapping 4,5 mm length 38	Ti6AL4V-ELI
MC.KV.170-40	cortical screw Self Tapping 4,5 mm length 40	Ti6AL4V-ELI
MC.KV.170-42	cortical screw Self Tapping 4,5 mm length 42	Ti6AL4V-ELI
MC.KV.170-44	cortical screw Self Tapping 4,5 mm length 44	Ti6AL4V-ELI
MC.KV.170-46	cortical screw Self Tapping 4,5 mm length 46	Ti6AL4V-ELI
MC.KV.170-48	cortical screw Self Tapping 4,5 mm length 48	Ti6AL4V-ELI
MC.KV.170-50	cortical screw Self Tapping 4,5 mm length 50	Ti6AL4V-ELI
MC.KV.131-08	2,5 mm locked cortical screw length 8 mm	Ti6AL4V-ELI
MC.KV.131-10	2,5 mm locked cortical screw length 10 mm	Ti6AL4V-ELI
MC.KV.131-12	2,5 mm locked cortical screw length 12 mm	Ti6AL4V-ELI
MC.KV.131-14	2,5 mm locked cortical screw length 14 mm	Ti6AL4V-ELI
MC.KV.131-16	2,5 mm locked cortical screw length 16 mm	Ti6AL4V-ELI
MC.KV.131-18	2,5 mm locked cortical screw length 18 mm	Ti6AL4V-ELI
MC.KV.131-20	2,5 mm locked cortical screw length 20 mm	Ti6AL4V-ELI
MC.KV.131-22	2,5 mm locked cortical screw length 22 mm	Ti6AL4V-ELI
MC.KV.131-24	2,5 mm locked cortical screw length 24 mm	Ti6AL4V-ELI
MC.KV.131-26	2,5 mm locked cortical screw length 26 mm	Ti6AL4V-ELI

# Orthopedic Screws



MC.KKV.141-10	locked cortical screw Self Tapping 2,7 mm lenght 10	Ti6AL4V-ELI
MC.KKV.141-12	locked cortical screw Self Tapping 2,7 mm lenght 12	Ti6AL4V-ELI
MC.KKV.141-14	locked cortical screw Self Tapping 2,7 mm lenght 14	Ti6AL4V-ELI
MC.KKV.141-16	locked cortical screw Self Tapping 2,7 mm lenght 16	Ti6AL4V-ELI
MC.KKV.141-18	locked cortical screw Self Tapping 2,7 mm lenght 18	Ti6AL4V-ELI
MC.KKV.141-20	locked cortical screw Self Tapping 2,7 mm lenght 20	Ti6AL4V-ELI
MC.KKV.141-22	locked cortical screw Self Tapping 2,7 mm lenght 22	Ti6AL4V-ELI
MC.KKV.141-24	locked cortical screw Self Tapping 2,7 mm lenght 24	Ti6AL4V-ELI
MC.KKV.141-26	locked cortical screw Self Tapping 2,7 mm lenght 26	Ti6AL4V-ELI
MC.KKV.141-28	locked cortical screw Self Tapping 2,7 mm lenght 28	Ti6AL4V-ELI
MC.KKV.141-30	locked cortical screw Self Tapping 2,7 mm lenght 30	Ti6AL4V-ELI
MC.KKV.141-32	locked cortical screw Self Tapping 2,7 mm lenght 32	Ti6AL4V-ELI
MC.KKV.141-34	locked cortical screw Self Tapping 2,7 mm lenght 34	Ti6AL4V-ELI
MC.KKV.141-36	locked cortical screw Self Tapping 2,7 mm lenght 36	Ti6AL4V-ELI
MC.KKV.141-38	locked cortical screw Self Tapping 2,7 mm lenght 38	Ti6AL4V-ELI
MC.KKV.141-40	locked cortical screw Self Tapping 2,7 mm lenght 40	Ti6AL4V-ELI
MC.KKV.141-42	locked cortical screw Self Tapping 2,7 mm lenght 42	Ti6AL4V-ELI
MC.KKV.141-44	locked cortical screw Self Tapping 2,7 mm lenght 44	Ti6AL4V-ELI
MC.KKV.141-46	locked cortical screw Self Tapping 2,7 mm lenght 46	Ti6AL4V-ELI
MC.KKV.141-48	locked cortical screw Self Tapping 2,7 mm lenght 48	Ti6AL4V-ELI
MC.KKV.141-50	locked cortical screw Self Tapping 2,7 mm lenght 50	Ti6AL4V-ELI
MC.KKV.151-10	locked cortical screw Self Tapping 3,5 mm lenght 10	Ti6AL4V-ELI
MC.KKV.151-12	locked cortical screw Self Tapping 3,5 mm lenght 12	Ti6AL4V-ELI
MC.KKV.151-14	locked cortical screw Self Tapping 3,5 mm lenght 14	Ti6AL4V-ELI
MC.KKV.151-16	locked cortical screw Self Tapping 3,5 mm lenght 16	Ti6AL4V-ELI
MC.KKV.151-18	locked cortical screw Self Tapping 3,5 mm lenght 18	Ti6AL4V-ELI
MC.KKV.151-20	locked cortical screw Self Tapping 3,5 mm lenght 20	Ti6AL4V-ELI
MC.KKV.151-22	locked cortical screw Self Tapping 3,5 mm lenght 22	Ti6AL4V-ELI
MC.KKV.151-24	locked cortical screw Self Tapping 3,5 mm lenght 24	Ti6AL4V-ELI
MC.KKV.151-26	locked cortical screw Self Tapping 3,5 mm lenght 26	Ti6AL4V-ELI
MC.KKV.151-28	locked cortical screw Self Tapping 3,5 mm lenght 28	Ti6AL4V-ELI
MC.KKV.151-30	locked cortical screw Self Tapping 3,5 mm lenght 30	Ti6AL4V-ELI
MC.KKV.151-32	locked cortical screw Self Tapping 3,5 mm lenght 32	Ti6AL4V-ELI

# Orthopedic Screws

MC.KKV.151-34	locked cortical screw Self Tapping 3,5 mm lenght 34	Ti6AL4V-ELI
MC.KKV.151-36	locked cortical screw Self Tapping 3,5 mm lenght 36	Ti6AL4V-ELI
MC.KKV.151-38	locked cortical screw Self Tapping 3,5 mm lenght 38	Ti6AL4V-ELI
MC.KKV.151-40	locked cortical screw Self Tapping 3,5 mm lenght 40	Ti6AL4V-ELI
MC.KKV.151-45	locked cortical screw Self Tapping 3,5 mm lenght 45	Ti6AL4V-ELI
MC.KKV.151-50	locked cortical screw Self Tapping 3,5 mm lenght 50	Ti6AL4V-ELI
MC.KKV.161.14	locked cortical screw Self Tapping 4,0 mm lenght 14	Ti6AL4V-ELI
MC.KKV.161.16	locked cortical screw Self Tapping 4,0 mm lenght 16	Ti6AL4V-ELI
MC.KKV.161.18	locked cortical screw Self Tapping 4,0 mm lenght 18	Ti6AL4V-ELI
MC.KKV.161.20	locked cortical screw Self Tapping 4,0 mm lenght 20	Ti6AL4V-ELI
MC.KKV.161.22	locked cortical screw Self Tapping 4,0 mm lenght 22	Ti6AL4V-ELI
MC.KKV.161.24	locked cortical screw Self Tapping 4,0 mm lenght 24	Ti6AL4V-ELI
MC.KKV.161.26	locked cortical screw Self Tapping 4,0 mm lenght 26	Ti6AL4V-ELI
MC.KKV.161.28	locked cortical screw Self Tapping 4,0 mm lenght 28	Ti6AL4V-ELI
MC.KKV.161.30	locked cortical screw Self Tapping 4,0 mm lenght 30	Ti6AL4V-ELI
MC.KKV.161.32	locked cortical screw Self Tapping 4,0 mm lenght 32	Ti6AL4V-ELI
MC.KKV.161.34	locked cortical screw Self Tapping 4,0 mm lenght 34	Ti6AL4V-ELI
MC.KKV.161.36	locked cortical screw Self Tapping 4,0 mm lenght 36	Ti6AL4V-ELI
MC.KKV.161.38	locked cortical screw Self Tapping 4,0 mm lenght 38	Ti6AL4V-ELI
MC.KKV.161.40	locked cortical screw Self Tapping 4,0 mm lenght 40	Ti6AL4V-ELI
MC.KKV.161.42	locked cortical screw Self Tapping 4,0 mm lenght 42	Ti6AL4V-ELI
MC.KKV.161.44	locked cortical screw Self Tapping 4,0 mm lenght 44	Ti6AL4V-ELI
MC.KKV.161.46	locked cortical screw Self Tapping 4,0 mm lenght 46	Ti6AL4V-ELI
MC.KKV.161.48	locked cortical screw Self Tapping 4,0 mm lenght 48	Ti6AL4V-ELI
MC.KKV.161.50	locked cortical screw Self Tapping 4,0 mm lenght 50	Ti6AL4V-ELI
MC.KKV.171-14	locked cortical screw Self Tapping 4,5 mm lenght 14	Ti6AL4V-ELI
MC.KKV.171-16	locked cortical screw Self Tapping 4,5 mm lenght 16	Ti6AL4V-ELI
MC.KKV.171-18	locked cortical screw Self Tapping 4,5 mm lenght 18	Ti6AL4V-ELI
MC.KKV.171-20	locked cortical screw Self Tapping 4,5 mm lenght 20	Ti6AL4V-ELI
MC.KKV.171-22	locked cortical screw Self Tapping 4,5 mm lenght 22	Ti6AL4V-ELI
MC.KKV.171-24	locked cortical screw Self Tapping 4,5 mm lenght 24	Ti6AL4V-ELI
MC.KKV.171-26	locked cortical screw Self Tapping 4,5 mm lenght 26	Ti6AL4V-ELI

# Orthopedic Screws

MC.KKV.171-28	locked cortical screw Self Tapping 4,5 mm length 28	Ti6AL4V-ELI
MC.KKV.171-30	locked cortical screw Self Tapping 4,5 mm length 30	Ti6AL4V-ELI
MC.KKV.171-32	locked cortical screw Self Tapping 4,5 mm length 32	Ti6AL4V-ELI
MC.KKV.171-34	locked cortical screw Self Tapping 4,5 mm length 34	Ti6AL4V-ELI
MC.KKV.171-36	locked cortical screw Self Tapping 4,5 mm length 36	Ti6AL4V-ELI
MC.KKV.171-38	locked cortical screw Self Tapping 4,5 mm length 38	Ti6AL4V-ELI
MC.KKV.171-40	locked cortical screw Self Tapping 4,5 mm length 40	Ti6AL4V-ELI
MC.KKV.171-42	locked cortical screw Self Tapping 4,5 mm length 42	Ti6AL4V-ELI
MC.KKV.171-44	locked cortical screw Self Tapping 4,5 mm length 44	Ti6AL4V-ELI
MC.KKV.171-46	locked cortical screw Self Tapping 4,5 mm length 46	Ti6AL4V-ELI
MC.KKV.171-48	locked cortical screw Self Tapping 4,5 mm length 48	Ti6AL4V-ELI
MC.KKV.171-50	locked cortical screw Self Tapping 4,5 mm length 50	Ti6AL4V-ELI
MC.KKKV.140-10	locked cannulated cortical screw 2,7 mm length 10	Ti6AL4V-ELI
MC.KKKV.140-12	locked cannulated cortical screw 2,7 mm length 12	Ti6AL4V-ELI
MC.KKKV.140-14	locked cannulated cortical screw 2,7 mm length 14	Ti6AL4V-ELI
MC.KKKV.140-16	locked cannulated cortical screw 2,7 mm length 16	Ti6AL4V-ELI
MC.KKKV.140-18	locked cannulated cortical screw 2,7 mm length 18	Ti6AL4V-ELI
MC.KKKV.140-20	locked cannulated cortical screw 2,7 mm length 20	Ti6AL4V-ELI
MC.KKKV.140-22	locked cannulated cortical screw 2,7 mm length 22	Ti6AL4V-ELI
MC.KKKV.140-24	locked cannulated cortical screw 2,7 mm length 24	Ti6AL4V-ELI
MC.KKKV.140-26	locked cannulated cortical screw 2,7 mm length 26	Ti6AL4V-ELI
MC.KKKV.140-28	locked cannulated cortical screw 2,7 mm length 28	Ti6AL4V-ELI
MC.KKKV.140-30	locked cannulated cortical screw 2,7 mm length 30	Ti6AL4V-ELI
MC.KKKV.140-32	locked cannulated cortical screw 2,7 mm length 32	Ti6AL4V-ELI
MC.KKKV.140-34	locked cannulated cortical screw 2,7 mm length 34	Ti6AL4V-ELI
MC.KKKV.140-36	locked cannulated cortical screw 2,7 mm length 36	Ti6AL4V-ELI
MC.KKKV.140-38	locked cannulated cortical screw 2,7 mm length 38	Ti6AL4V-ELI
MC.KKKV.140-40	locked cannulated cortical screw 2,7 mm length 40	Ti6AL4V-ELI
MC.KKKV.140-42	locked cannulated cortical screw 2,7 mm length 42	Ti6AL4V-ELI
MC.KKKV.140-44	locked cannulated cortical screw 2,7 mm length 44	Ti6AL4V-ELI
MC.KKKV.140-46	locked cannulated cortical screw 2,7 mm length 46	Ti6AL4V-ELI
MC.KKKV.140-48	locked cannulated cortical screw 2,7 mm length 48	Ti6AL4V-ELI
MC.KKKV.140-50	locked cannulated cortical screw 2,7 mm length 50	Ti6AL4V-ELI

# Orthopedic Screws

MC.KKKV.150-10	locked cannulated cortical screw 3,5 mm length 10	Ti6AL4V-ELI
MC.KKKV.150-12	locked cannulated cortical screw 3,5 mm length 12	Ti6AL4V-ELI
MC.KKKV.150-14	locked cannulated cortical screw 3,5 mm length 14	Ti6AL4V-ELI
MC.KKKV.150-16	locked cannulated cortical screw 3,5 mm length 16	Ti6AL4V-ELI
MC.KKKV.150-18	locked cannulated cortical screw 3,5 mm length 18	Ti6AL4V-ELI
MC.KKKV.150-20	locked cannulated cortical screw 3,5 mm length 20	Ti6AL4V-ELI
MC.KKKV.150-22	locked cannulated cortical screw 3,5 mm length 22	Ti6AL4V-ELI
MC.KKKV.150-24	locked cannulated cortical screw 3,5 mm length 24	Ti6AL4V-ELI
MC.KKKV.150-26	locked cannulated cortical screw 3,5 mm length 26	Ti6AL4V-ELI
MC.KKKV.150-28	locked cannulated cortical screw 3,5 mm length 28	Ti6AL4V-ELI
MC.KKKV.150-30	locked cannulated cortical screw 3,5 mm length 30	Ti6AL4V-ELI
MC.KKKV.150-32	locked cannulated cortical screw 3,5 mm length 32	Ti6AL4V-ELI
MC.KKKV.150-34	locked cannulated cortical screw 3,5 mm length 34	Ti6AL4V-ELI
MC.KKKV.150-36	locked cannulated cortical screw 3,5 mm length 36	Ti6AL4V-ELI
MC.KKKV.150-38	locked cannulated cortical screw 3,5 mm length 38	Ti6AL4V-ELI
MC.KKKV.150-40	locked cannulated cortical screw 3,5 mm length 40	Ti6AL4V-ELI
MC.KKKV.150-42	locked cannulated cortical screw 3,5 mm length 42	Ti6AL4V-ELI
MC.KKKV.150-44	locked cannulated cortical screw 3,5 mm length 44	Ti6AL4V-ELI
MC.KKKV.150-46	locked cannulated cortical screw 3,5 mm length 46	Ti6AL4V-ELI
MC.KKKV.150-48	locked cannulated cortical screw 3,5 mm length 48	Ti6AL4V-ELI
MC.KKKV.150-50	locked cannulated cortical screw 3,5 mm length 50	Ti6AL4V-ELI
MC.KKKV.170-14	locked cannulated cortical screw 4,5 mm length 14	Ti6AL4V-ELI
MC.KKKV.170-16	locked cannulated cortical screw 4,5 mm length 16	Ti6AL4V-ELI
MC.KKKV.170-18	locked cannulated cortical screw 4,5 mm length 18	Ti6AL4V-ELI
MC.KKKV.170-20	locked cannulated cortical screw 4,5 mm length 20	Ti6AL4V-ELI
MC.KKKV.170-22	locked cannulated cortical screw 4,5 mm length 22	Ti6AL4V-ELI
MC.KKKV.170-24	locked cannulated cortical screw 4,5 mm length 24	Ti6AL4V-ELI
MC.KKKV.170-26	locked cannulated cortical screw 4,5 mm length 26	Ti6AL4V-ELI
MC.KKKV.170-28	locked cannulated cortical screw 4,5 mm length 28	Ti6AL4V-ELI
MC.KKKV.170-30	locked cannulated cortical screw 4,5 mm length 30	Ti6AL4V-ELI
MC.KKKV.170-32	locked cannulated cortical screw 4,5 mm length 32	Ti6AL4V-ELI
MC.KKKV.170-34	locked cannulated cortical screw 4,5 mm length 34	Ti6AL4V-ELI
MC.KKKV.170-36	locked cannulated cortical screw 4,5 mm length 36	Ti6AL4V-ELI

# Orthopedic Screws



MC.KKKV.170-38	locked cannulated cortical screw 4,5 mm length 38	Ti6AL4V-ELI
MC.KKKV.170-40	locked cannulated cortical screw 4,5 mm length 40	Ti6AL4V-ELI
MC.KKKV.170-42	locked cannulated cortical screw 4,5 mm length 42	Ti6AL4V-ELI
MC.KKKV.170-44	locked cannulated cortical screw 4,5 mm length 44	Ti6AL4V-ELI
MC.KKKV.170-46	locked cannulated cortical screw 4,5 mm length 46	Ti6AL4V-ELI
MC.KKKV.170-48	locked cannulated cortical screw 4,5 mm length 48	Ti6AL4V-ELI
MC.KKKV.170-50	locked cannulated cortical screw 4,5 mm length 50	Ti6AL4V-ELI
MC.KÜV.140-14	cannulated screw 2,7 mm length 14	Ti6AL4V-ELI
MC.KÜV.140-16	cannulated screw 2,7 mm length 15	Ti6AL4V-ELI
MC.KÜV.140-18	cannulated screw 2,7 mm length 16	Ti6AL4V-ELI
MC.KÜV.140-20	cannulated screw 2,7 mm length 17	Ti6AL4V-ELI
MC.KÜV.140-22	cannulated screw 2,7 mm length 18	Ti6AL4V-ELI
MC.KÜV.140-24	cannulated screw 2,7 mm length 19	Ti6AL4V-ELI
MC.KÜV.140-26	cannulated screw 2,7 mm length 20	Ti6AL4V-ELI
MC.KÜV.140-28	cannulated screw 2,7 mm length 21	Ti6AL4V-ELI
MC.KÜV.140-30	cannulated screw 2,7 mm length 22	Ti6AL4V-ELI
MC.KÜV.140-32	cannulated screw 2,7 mm length 23	Ti6AL4V-ELI
MC.KÜV.140-34	cannulated screw 2,7 mm length 24	Ti6AL4V-ELI
MC.KÜV.140-36	cannulated screw 2,7 mm length 25	Ti6AL4V-ELI
MC.KÜV.140-38	cannulated screw 2,7 mm length 26	Ti6AL4V-ELI
MC.KÜV.140-40	cannulated screw 2,7 mm length 27	Ti6AL4V-ELI
MC.KÜV.140-42	cannulated screw 2,7 mm length 28	Ti6AL4V-ELI
MC.KÜV.140-44	cannulated screw 2,7 mm length 29	Ti6AL4V-ELI
MC.KÜV.140-46	cannulated screw 2,7 mm length 30	Ti6AL4V-ELI
MC.KÜV.150-14	cannulated screw 3,5 mm length 14	Ti6AL4V-ELI
MC.KÜV.150-16	cannulated screw 3,5 mm length 16	Ti6AL4V-ELI
MC.KÜV.150-18	cannulated screw 3,5 mm length 18	Ti6AL4V-ELI
MC.KÜV.150-20	cannulated screw 3,5 mm length 20	Ti6AL4V-ELI
MC.KÜV.150-22	cannulated screw 3,5 mm length 22	Ti6AL4V-ELI
MC.KÜV.150-24	cannulated screw 3,5 mm length 24	Ti6AL4V-ELI
MC.KÜV.150-26	cannulated screw 3,5 mm length 26	Ti6AL4V-ELI
MC.KÜV.150-28	cannulated screw 3,5 mm length 28	Ti6AL4V-ELI
MC.KÜV.150-30	cannulated screw 3,5 mm length 30	Ti6AL4V-ELI


# Orthopedic Screws



MC.KÜV.150-32	cannulated screw 3,5 mm length 32	Ti6AL4V-ELI
MC.KÜV.150-34	cannulated screw 3,5 mm length 34	Ti6AL4V-ELI
MC.KÜV.150-36	cannulated screw 3,5 mm length 36	Ti6AL4V-ELI
MC.KÜV.150-38	cannulated screw 3,5 mm length 38	Ti6AL4V-ELI
MC.KÜV.150-40	cannulated screw 3,5 mm length 40	Ti6AL4V-ELI
MC.KÜV.150-42	cannulated screw 3,5 mm length 42	Ti6AL4V-ELI
MC.KÜV.150-44	cannulated screw 3,5 mm length 44	Ti6AL4V-ELI
MC.KÜV.150-46	cannulated screw 3,5 mm length 46	Ti6AL4V-ELI
MC.KÜV.150-48	cannulated screw 3,5 mm length 48	Ti6AL4V-ELI
MC.KÜV.150-50	cannulated screw 3,5 mm length 50	Ti6AL4V-ELI
MC.KÜV.170-20	cannulated screw 4,5 mm length 20	Ti6AL4V-ELI
MC.KÜV.170-22	cannulated screw 4,5 mm length 22	Ti6AL4V-ELI
MC.KÜV.170-24	cannulated screw 4,5 mm length 24	Ti6AL4V-ELI
MC.KÜV.170-26	cannulated screw 4,5 mm length 26	Ti6AL4V-ELI
MC.KÜV.170-28	cannulated screw 4,5 mm length 28	Ti6AL4V-ELI
MC.KÜV.170-30	cannulated screw 4,5 mm length 30	Ti6AL4V-ELI
MC.KÜV.170-32	cannulated screw 4,5 mm length 32	Ti6AL4V-ELI
MC.KÜV.170-34	cannulated screw 4,5 mm length 34	Ti6AL4V-ELI
MC.KÜV.170-36	cannulated screw 4,5 mm length 36	Ti6AL4V-ELI
MC.KÜV.170-38	cannulated screw 4,5 mm length 38	Ti6AL4V-ELI
MC.KÜV.170-40	cannulated screw 4,5 mm length 40	Ti6AL4V-ELI
MC.KÜV.170-42	cannulated screw 4,5 mm length 42	Ti6AL4V-ELI
MC.KÜV.170-44	cannulated screw 4,5 mm length 44	Ti6AL4V-ELI
MC.KÜV.170-46	cannulated screw 4,5 mm length 46	Ti6AL4V-ELI
MC.KÜV.170-48	cannulated screw 4,5 mm length 48	Ti6AL4V-ELI
MC.KÜV.170-50	cannulated screw 4,5 mm length 50	Ti6AL4V-ELI
MC.KÜV.180.14	cannulated screw 5,0 mm length 14	Ti6AL4V-ELI
MC.KÜV.180.16	cannulated screw 5,0 mm length 16	Ti6AL4V-ELI
MC.KÜV.180.18	cannulated screw 5,0 mm length 18	Ti6AL4V-ELI
MC.KÜV.180-20	cannulated screw 5,0 mm length 20	Ti6AL4V-ELI
MC.KÜV.180-22	cannulated screw 5,0 mm length 22	Ti6AL4V-ELI
MC.KÜV.180-24	cannulated screw 5,0 mm length 24	Ti6AL4V-ELI
MC.KÜV.180-26	cannulated screw 5,0 mm length 26	Ti6AL4V-ELI



# Orthopedic Screws



MC.KÜV.180-28	cannulated screw 5,0 mm length 28	Ti6AL4V-ELI
MC.KÜV.180-30	cannulated screw 5,0 mm length 30	Ti6AL4V-ELI
MC.KÜV.180-32	cannulated screw 5,0 mm length 32	Ti6AL4V-ELI
MC.KÜV.180-34	cannulated screw 5,0 mm length 34	Ti6AL4V-ELI
MC.KÜV.180-36	cannulated screw 5,0 mm length 36	Ti6AL4V-ELI
MC.KÜV.180-38	cannulated screw 5,0 mm length 38	Ti6AL4V-ELI
MC.KÜV.180-40	cannulated screw 5,0 mm length 40	Ti6AL4V-ELI
MC.KÜV.180-42	cannulated screw 5,0 mm length 42	Ti6AL4V-ELI
MC.KÜV.180-44	cannulated screw 5,0 mm length 44	Ti6AL4V-ELI
MC.KÜV.180-46	cannulated screw 5,0 mm length 46	Ti6AL4V-ELI
MC.KÜV.180-48	cannulated screw 5,0 mm length 48	Ti6AL4V-ELI
MC.KÜV.180-50	cannulated screw 5,0 mm length 50	Ti6AL4V-ELI
MC.KÜV.190-14	cannulated screw 6,5 mm length 14	Ti6AL4V-ELI
MC.KÜV.190-16	cannulated screw 6,5 mm length 16	Ti6AL4V-ELI
MC.KÜV.190-18	cannulated screw 6,5 mm length 18	Ti6AL4V-ELI
MC.KÜV.190-20	cannulated screw 6,5 mm length 20	Ti6AL4V-ELI
MC.KÜV.190-22	cannulated screw 6,5 mm length 22	Ti6AL4V-ELI
MC.KÜV.190-24	cannulated screw 6,5 mm length 24	Ti6AL4V-ELI
MC.KÜV.190-26	cannulated screw 6,5 mm length 26	Ti6AL4V-ELI
MC.KÜV.190-28	cannulated screw 6,5 mm length 28	Ti6AL4V-ELI
MC.KÜV.190-30	cannulated screw 6,5 mm length 30	Ti6AL4V-ELI
MC.KÜV.190-35	cannulated screw 6,5 mm length 35	Ti6AL4V-ELI
MC.KÜV.190-40	cannulated screw 6,5 mm length 40	Ti6AL4V-ELI
MC.KÜV.190-45	cannulated screw 6,5 mm length 45	Ti6AL4V-ELI
MC.KÜV.190-50	cannulated screw 6,5 mm length 50	Ti6AL4V-ELI
MC.BKV.100-08	headless cannulated compression screw 2.7X8MM	Ti6AL4V-ELI
MC.BKV.100-10	headless cannulated compression screw 2.7X10MM	Ti6AL4V-ELI
MC.BKV.100-12	headless cannulated compression screw 2.7X12MM	Ti6AL4V-ELI
MC.BKV.100-14	headless cannulated compression screw 2.7X14MM	Ti6AL4V-ELI
MC.BKV.100-16	headless cannulated compression screw 2.7X16MM	Ti6AL4V-ELI
MC.BKV.100-18	headless cannulated compression screw 2.7X18MM	Ti6AL4V-ELI
MC.BKV.100-20	headless cannulated compression screw 2.7X20MM	Ti6AL4V-ELI
MC.BKV.100-22	headless cannulated compression screw 2.7X22MM	Ti6AL4V-ELI

# Orthopedic Screws



MC.BKV.100-24	headless cannulated compression screw 2.7X24MM	Ti6AL4V-ELI
MC.BKV.100-26	headless cannulated compression screw 2.7X26MM	Ti6AL4V-ELI
MC.BKV.100-28	headless cannulated compression screw 2.7X28MM	Ti6AL4V-ELI
MC.BKV.100-30	headless cannulated compression screw 2.7X30MM	Ti6AL4V-ELI
MC.BKV.110-12	headless cannulated compression screw 3.0X12MM	Ti6AL4V-ELI
MC.BKV.110-14	headless cannulated compression screw 3.0X14MM	Ti6AL4V-ELI
MC.BKV.110-16	headless cannulated compression screw 3.0X16MM	Ti6AL4V-ELI
MC.BKV.110-18	headless cannulated compression screw 3.0X18MM	Ti6AL4V-ELI
MC.BKV.110-20	headless cannulated compression screw 3.0X20MM	Ti6AL4V-ELI
MC.BKV.110-22	headless cannulated compression screw 3.0X22MM	Ti6AL4V-ELI
MC.BKV.110-24	headless cannulated compression screw 3.0X24MM	Ti6AL4V-ELI
MC.BKV.110-26	headless cannulated compression screw 3.0X26MM	Ti6AL4V-ELI
MC.BKV.110-28	headless cannulated compression screw 3.0X28MM	Ti6AL4V-ELI
MC.BKV.110-30	headless cannulated compression screw 3.0X30MM	Ti6AL4V-ELI
MC.BKV.110-32	headless cannulated compression screw 3.0X32MM	Ti6AL4V-ELI
MC.BKV.110-34	headless cannulated compression screw 3.0X34MM	Ti6AL4V-ELI
MC.BKV.110-36	headless cannulated compression screw 3.0X36MM	Ti6AL4V-ELI
MC.BKV.110-38	headless cannulated compression screw 3.0X38MM	Ti6AL4V-ELI
MC.BKV.110-40	headless cannulated compression screw 3.0X40MM	Ti6AL4V-ELI
MC.BKV.120-14	headless cannulated compression screw 3.5X14MM	Ti6AL4V-ELI
MC.BKV.120-16	headless cannulated compression screw 3.5X16MM	Ti6AL4V-ELI
MC.BKV.120-18	headless cannulated compression screw 3.5X18MM	Ti6AL4V-ELI
MC.BKV.120-20	headless cannulated compression screw 3.5X20MM	Ti6AL4V-ELI
MC.BKV.120-22	headless cannulated compression screw 3.5X22MM	Ti6AL4V-ELI
MC.BKV.120-24	headless cannulated compression screw 3.5X24MM	Ti6AL4V-ELI
MC.BKV.120-26	headless cannulated compression screw 3.5X26MM	Ti6AL4V-ELI
MC.BKV.120-28	headless cannulated compression screw 3.5X28MM	Ti6AL4V-ELI
MC.BKV.120-30	headless cannulated compression screw 3.5X30MM	Ti6AL4V-ELI
MC.BKV.120-32	headless cannulated compression screw 3.5X32MM	Ti6AL4V-ELI
MC.BKV.120-34	headless cannulated compression screw 3.5X34MM	Ti6AL4V-ELI
MC.BKV.120-36	headless cannulated compression screw 3.5X36MM	Ti6AL4V-ELI
MC.BKV.120-38	headless cannulated compression screw 3.5X38MM	Ti6AL4V-ELI
MC.BKV.120-40	headless cannulated compression screw 3.5X40MM	Ti6AL4V-ELI

# Orthopedic Screws



MC.BKV.120-42	headless cannulated compression screw 3,5X45MM	Ti6AL4V-ELI
MC.BKV.120-44	headless cannulated compression screw 3,5X50MM	Ti6AL4V-ELI
MC.BKV.130-20	headless cannulated compression screw 4,0X20MM	Ti6AL4V-ELI
MC.BKV.130-22	headless cannulated compression screw 4,0X22MM	Ti6AL4V-ELI
MC.BKV.130-24	headless cannulated compression screw 4,0X24MM	Ti6AL4V-ELI
MC.BKV.130-26	headless cannulated compression screw 4,0X26MM	Ti6AL4V-ELI
MC.BKV.130-28	headless cannulated compression screw 4,0X28MM	Ti6AL4V-ELI
MC.BKV.130-30	headless cannulated compression screw 4,0X30MM	Ti6AL4V-ELI
MC.BKV.130-32	headless cannulated compression screw 4,0X32MM	Ti6AL4V-ELI
MC.BKV.130-34	headless cannulated compression screw 4,0X34MM	Ti6AL4V-ELI
MC.BKV.130-36	headless cannulated compression screw 4,0X36MM	Ti6AL4V-ELI
MC.BKV.130-38	headless cannulated compression screw 4,0X38MM	Ti6AL4V-ELI
MC.BKV.130-40	headless cannulated compression screw 4,0X40MM	Ti6AL4V-ELI
MC.BKV.130-42	headless cannulated compression screw 4,0X42MM	Ti6AL4V-ELI
MC.BKV.130-44	headless cannulated compression screw 4,0X44MM	Ti6AL4V-ELI
MC.BKV.130-46	headless cannulated compression screw 4,0X46MM	Ti6AL4V-ELI
MC.BKV.130-48	headless cannulated compression screw 4,0X48MM	Ti6AL4V-ELI
MC.BKV.130-50	headless cannulated compression screw 4,0X50MM	Ti6AL4V-ELI
MC.BKV.140-24	headless cannulated compression screw 5,0X24MM	Ti6AL4V-ELI
MC.BKV.140-26	headless cannulated compression screw 5,0X26MM	Ti6AL4V-ELI
MC.BKV.140-28	headless cannulated compression screw 5,0X28MM	Ti6AL4V-ELI
MC.BKV.140-30	headless cannulated compression screw 5,0X30MM	Ti6AL4V-ELI
MC.BKV.140-35	headless cannulated compression screw 5,0X35MM	Ti6AL4V-ELI
MC.BKV.140-40	headless cannulated compression screw 5,0X40MM	Ti6AL4V-ELI
MC.BKV.140-45	headless cannulated compression screw 5,0X45MM	Ti6AL4V-ELI
MC.BKV.140-50	headless cannulated compression screw 5,0X50MM	Ti6AL4V-ELI
MC.BKV.150-24	headless cannulated compression screw 6,5X24MM	Ti6AL4V-ELI
MC.BKV.150-26	headless cannulated compression screw 6,5X26MM	Ti6AL4V-ELI
MC.BKV.150-28	headless cannulated compression screw 6,5X28MM	Ti6AL4V-ELI
MC.BKV.150-30	headless cannulated compression screw 6,5X30MM	Ti6AL4V-ELI
MC.BKV.150-32	headless cannulated compression screw 6,5X32MM	Ti6AL4V-ELI
MC.BKV.150-34	headless cannulated compression screw 6,5X34MM	Ti6AL4V-ELI
MC.BKV.150-36	headless cannulated compression screw 6,5X36MM	Ti6AL4V-ELI

# Orthopedic Screws

MC.BKV.150-38	headless cannulated compression screw 6,5X38MM	Ti6AL4V-ELI
MC.BKV.150-40	headless cannulated compression screw 6,5X40MM	Ti6AL4V-ELI
MC.BKV.150-45	headless cannulated compression screw 6,5X45MM	Ti6AL4V-ELI
MC.BKV.150-50	headless cannulated compression screw 6,5X50MM	Ti6AL4V-ELI
MC.KDSV.100-16	cannulated Double compression Screw 3.0 x 16mm	Ti6AL4V-ELI
MC.KDSV.100-18	cannulated Double compression Screw 3.0 x 18mm	Ti6AL4V-ELI
MC.KDSV.100-20	cannulated Double compression Screw 3.0 x 20mm	Ti6AL4V-ELI
MC.KDSV.100-22	cannulated Double compression Screw 3.0 x 22mm	Ti6AL4V-ELI
MC.KDSV.100-24	cannulated Double compression Screw 3.0 x 24mm	Ti6AL4V-ELI
MC.KDSV.100-26	cannulated Double compression Screw 3.0 x 26mm	Ti6AL4V-ELI
MC.KDSV.100-28	cannulated Double compression Screw 3.0 x 28mm	Ti6AL4V-ELI
MC.KDSV.100-30	cannulated Double compression Screw 3.0 x 30mm	Ti6AL4V-ELI
MC.KDSV.110-16	cannulated Double compression Screw 4.5 x 16mm	Ti6AL4V-ELI
MC.KDSV.110-18	cannulated Double compression Screw 4.5 x 18mm	Ti6AL4V-ELI
MC.KDSV.110-20	cannulated Double compression Screw 4.5 x 20mm	Ti6AL4V-ELI
MC.KDSV.110-22	cannulated Double compression Screw 4.5 x 22mm	Ti6AL4V-ELI
MC.KDSV.110-24	cannulated Double compression Screw 4.5 x 24mm	Ti6AL4V-ELI
MC.KDSV.110-26	cannulated Double compression Screw 4.5 x 26mm	Ti6AL4V-ELI
MC.KDSV.110-28	cannulated Double compression Screw 4.5 x 28mm	Ti6AL4V-ELI
MC.KDSV.110-30	cannulated Double compression Screw 4.5 x 30mm	Ti6AL4V-ELI
MC.KDSV.110-32	cannulated Double compression Screw 4.5 x 32mm	Ti6AL4V-ELI
MC.KDSV.110-34	cannulated Double compression Screw 4.5 x 34mm	Ti6AL4V-ELI
MC.KDSV.110-36	cannulated Double compression Screw 4.5 x 36mm	Ti6AL4V-ELI
MC.KDSV.110-38	cannulated Double compression Screw 4.5 x 38mm	Ti6AL4V-ELI
MC.KDSV.110-40	cannulated Double compression Screw 4.5 x 40mm	Ti6AL4V-ELI
MC.KDSV.110-42	cannulated Double compression Screw 4.5 x 42mm	Ti6AL4V-ELI
MC.KDSV.110-44	cannulated Double compression Screw 4.5 x 44mm	Ti6AL4V-ELI
MC.KDSV.110-46	cannulated Double compression Screw 4.5 x 46mm	Ti6AL4V-ELI
MC.KDSV.110-48	cannulated Double compression Screw 4.5 x 48mm	Ti6AL4V-ELI
MC.KDSV.110-50	cannulated Double compression Screw 4.5 x 50mm	Ti6AL4V-ELI
MC.KDSV.120-16	cannulated Double compression Screw 5.5 x 16mm	Ti6AL4V-ELI
MC.KDSV.120-18	cannulated Double compression Screw 5.5 x 18mm	Ti6AL4V-ELI
MC.KDSV.120-20	cannulated Double compression Screw 5.5 x 20mm	Ti6AL4V-ELI

# Orthopedic Screws



MC.KDSV.120-22	cannulated Double compression Screw 5.5 x 22mm	Ti6AL4V-ELI
MC.KDSV.120-24	cannulated Double compression Screw 5.5 x 24mm	Ti6AL4V-ELI
MC.KDSV.120-26	cannulated Double compression Screw 5.5 x 26mm	Ti6AL4V-ELI
MC.KDSV.120-28	cannulated Double compression Screw 5.5 x 28mm	Ti6AL4V-ELI
MC.KDSV.120-30	cannulated Double compression Screw 5.5 x 30mm	Ti6AL4V-ELI
MC.KDSV.120-32	cannulated Double compression Screw 5.5 x 32mm	Ti6AL4V-ELI
MC.KDSV.120-34	cannulated Double compression Screw 5.5 x 34mm	Ti6AL4V-ELI
MC.KDSV.120-36	cannulated Double compression Screw 5.5 x 36mm	Ti6AL4V-ELI
MC.KDSV.120-38	cannulated Double compression Screw 5.5 x 38mm	Ti6AL4V-ELI
MC.KDSV.120-40	cannulated Double compression Screw 5.5 x 40mm	Ti6AL4V-ELI
MC.KDSV.120-42	cannulated Double compression Screw 5.5 x 42mm	Ti6AL4V-ELI
MC.KDSV.120-44	cannulated Double compression Screw 5.5 x 44mm	Ti6AL4V-ELI
MC.KDSV.120-46	cannulated Double compression Screw 5.5 x 46mm	Ti6AL4V-ELI
MC.KDSV.120-48	cannulated Double compression Screw 5.5 x 48mm	Ti6AL4V-ELI
MC.KDSV.120-50	cannulated Double compression Screw 5.5 x 50mm	Ti6AL4V-ELI
MC.KDSV.130-16	cannulated Double compression Screw 6.5 x 16mm	Ti6AL4V-ELI
MC.KDSV.130-18	cannulated Double compression Screw 6.5 x 18mm	Ti6AL4V-ELI
MC.KDSV.130-20	cannulated Double compression Screw 6.5 x 20mm	Ti6AL4V-ELI
MC.KDSV.130-22	cannulated Double compression Screw 6.5 x 22mm	Ti6AL4V-ELI
MC.KDSV.130-24	cannulated Double compression Screw 6.5 x 24mm	Ti6AL4V-ELI
MC.KDSV.130-26	cannulated Double compression Screw 6.5 x 26mm	Ti6AL4V-ELI
MC.KDSV.130-28	cannulated Double compression Screw 6.5 x 28mm	Ti6AL4V-ELI
MC.KDSV.130-30	cannulated Double compression Screw 6.5 x 30mm	Ti6AL4V-ELI
MC.KDSV.130-32	cannulated Double compression Screw 6.5 x 32mm	Ti6AL4V-ELI
MC.KDSV.130-34	cannulated Double compression Screw 6.5 x 34mm	Ti6AL4V-ELI
MC.KDSV.130-36	cannulated Double compression Screw 6.5 x 36mm	Ti6AL4V-ELI
MC.KDSV.130-38	cannulated Double compression Screw 6.5 x 38mm	Ti6AL4V-ELI
MC.KDSV.130-40	cannulated Double compression Screw 6.5 x 40mm	Ti6AL4V-ELI
MC.KDSV.130-42	cannulated Double compression Screw 6.5 x 42mm	Ti6AL4V-ELI
MC.KDSV.130-44	cannulated Double compression Screw 6.5 x 44mm	Ti6AL4V-ELI
MC.KDSV.130-46	cannulated Double compression Screw 6.5 x 46mm	Ti6AL4V-ELI
MC.KDSV.130-48	cannulated Double compression Screw 6.5 x 48mm	Ti6AL4V-ELI
MC.KDSV.130-50	cannulated Double compression Screw 6.5 x 50mm	Ti6AL4V-ELI

